

Product datasheet for **SC331856**

BRF1 (NM_001242786) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	BRF1 (NM_001242786) Human Untagged Clone
Tag:	Tag Free
Symbol:	BRF1
Synonyms:	BRF; BRF-1; CFDS; GTF3B; hBRF; HEL-S-76p; TAF3B2; TAF3C; TAFIII90; TF3B90; TFIIB90
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC331856 representing NM_001242786. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGCCGTGAGCAGGCACCTGACCCGCGGCCGGAAGATGGCCACGTGATTGCTGCCTGCCTCTACCTG
GTCTGCCGTACGGAGGGCAGCCGCACATGCTCCTGGACCTCAGCGACCTGCTCCAGGTGAATGTGTAC
GTGCTTGGAAAGACGTTTCTTCTTGGCAAGAGAGCTCTGCATCAATGCGCCGGCCATAGACCCGTGC
CTGTATATTCCACGCTTTGCGCACCTGCTGGAATTCGGGGAGAAGAACCACGAGGTGCCATGACTGCC
CTGAGGCTCTACAGAGGATGAAGCGGGACTGGATGCACACAGGCCGGCGCCCTCGGGCCTCTGCGGA
GCAGCGCTCCTGGTTGCAGCCAGAATGCATGACTTCAGGAGGACTGTGAAGGAGTTCATCAGTGTGGTC
AAAGTGTGTGAGTCCACGCTGCGGAAGAGGCTCACGGAATTTGAAGACACCCACCAGTCAGTTGACC
ATTGATGAGTTCATGAAGATCGACCTGGAGGAGGAGTGCACCCCTCGTACACAGCTGGGCAGAGG
AAGCTGCGGATGAAGCAGCTTGAACAAGTCTGTCAAAAAACTGGAGGAGGTTGAAGGTGAAATATCC
AGTTACCAGGATGCAATTGAGATTGAACTAGAAAAACAGCCGGCCAAAGGCCAAGGGGGCCCTGGCCAGC
CTGGCAAAAGATGGCTCCACCGAGGACACCCGCTCCAGCTTGTGTGGCAGGAGGACACAGAGGACGAG
GAGCTGGAAGCCGCGCCAGCCACCTGAACAAAGACTTATACCGGGAGCTCCTTGGTGGTCCCGCCGGC
AGCTCGGAAGCAGCAGGAAGCCCGAGTGGGGCGCAGACCTCCGGCCCTGGGGTCCCTGCTGACCCC
CTCCCCACTGCAGCCAGCCTGGGCATCTCAGACTCCATCCGCGAATGCATCTCCTCTCAGAGCAGCGAC
CCCAAAGATGCTTCAGGAGACGGTGGAGCTGGACCTCAGTGGCATTGATGACCTGGAGATTGACAGGCGC
GACCTTAGCATGCCCCGGTGTGCGAAGGCTAAAAGCCAGCCCACTTCCCTGTGCTCGCCAGTACATC
CTGAATGAGTCGGAAGCCCGCTGAAGGCCGAGCTGTGGATGAGGGAGAACGCCGAGTACCTGCGGGAA
CAGAGGGAAAAAGAAGCAAGAATAGCGAAAGAGAAGGAGCTCGGCATCTACAAGGAACACAAGCCCAAG
AAGTCTTGCAAGCGACGGGAGCCAATTCAGGCCAGTACCGCCAGGGAGGCCATCGAGAAGATGCTGGAG
CAGAAGAAGATCTCCAGCAAGATCAATTATAGCGTGTCCGGGGCCTCAGCAGCGCCGGCGGGGAGT
CCGCACAGGGAGGATGCACAGCCGAGCATAGCGCCAGTGCAGGAAGCTGTACGAAGGAGGACGCCG
GCCAGCAGAAGTGGGGCTGACCCTGTGACCAGTGTGGGAAAAGTTGAGCCTCTGGTGTCTACGCAG
CCAGCAAAGAAGGTGGCCACGGGAGAGGCTTTGCTCCAAAGCTCTCCACCTCGGAGCTGAGCCTGCC
AGGCCCCAGGCGGTGCTGGTGGAGAGCGGGCCCGTGTACATACCACGCCGACGAGGAGGCTGACGAGGAG
GAGCCTGACGAGGAGGACGGGGAGCCCTGCGTCACTGCCCTGCAGATGATGGGCAGCAACGACTATGGC
TGTGATGGCGATGAGGACGACGGCTACTGA

```

Restriction Sites: SgfI-MluI



[View online »](#)

ACCN:	NM_001242786
Insert Size:	1755 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001242786.1
RefSeq Size:	3452 bp
RefSeq ORF:	1755 bp
Locus ID:	2972
UniProt ID:	Q92994
Cytogenetics:	14q32.33
Protein Families:	Transcription Factors
MW:	64.3 kDa
Gene Summary:	<p>This gene encodes one of the three subunits of the RNA polymerase III transcription factor complex. This complex plays a central role in transcription initiation by RNA polymerase III on genes encoding tRNA, 5S rRNA, and other small structural RNAs. The gene product belongs to the TF2B family. Several alternatively spliced variants encoding different isoforms, that function at different promoters transcribed by RNA polymerase III, have been identified. [provided by RefSeq, Jun 2011]</p> <p>Transcript Variant: This variant (4) differs in the 5' UTR and coding sequence and uses an alternate in-frame splice junction at the 5' end of an exon compared to variant (1). The resulting isoform (4) is shorter at the N-terminus and contains an insert compared to isoform 1.</p>