

## Product datasheet for **SC331861**

### AP180 (SNAP91) (NM\_001242792) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	AP180 (SNAP91) (NM_001242792) Human Untagged Clone
Tag:	Tag Free
Symbol:	AP180
Synonyms:	AP180; CALM
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

**Fully Sequenced ORF:** >SC331861 representing NM\_001242792.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGTCGGGCCAAACGCTCACGGATCGGATCGCCGCCGCTCAGTACAGCGTTACAGGCTCTGCTGTAGCA
AGAGCGGTCTGCAAAGCCACTACTCATGAAGTAATGGGCCCAAGAAAAAGCACCTGGACTATTTGATC
CAGGCTACCAACGAGACCAATGTTAATATTCTCAGATGGCCGACACTCTCTTTGAGCGGGCAACAAAC
AGTAGCTGGGTGGTGTGTTAAGGCTTTAGTGACAACACATCATCTCATGGTGCATGGAATGAGAGA
TTTATTCAATATTTGGCTTCTAGAAATACACTATTCAATCTCAGCAATTTTTTGGACAAAAGTGGATCC
CATGGTTATGATATGTCTACCTTCATAAGGCGCTATAGTAGATATTTGAATGAAAAGGCTTTTTCTTAC
AGACAGATGGCCTTTGATTTTGCCAGGGTGAAGAAAGGGGCCGATGGTGAATGAGGACAATGGCTCCC
GAAAAGCTGCTAAAGAGTATGCCAATACTACAGGGACAAATTGATGCACTGCTTGAATTTGATGTGCAT
CCAAATGAACTAACAAATGGTGCATAAATGCAGCATTATGCTTCTTTCAAAGATCTTATCAAACCTT
TTTGCTTGCTACAATGATGGTGTATTAACCTACTCGAAAAGTTTTTGAATGAAGAAAGGACAATGT
AAAGATGCTCTAGAAATTTACAAACGATTTCTAACTAGAATGACACGAGTGTCTGAATTTCTCAAGTT
GCAGAGCAAGTTGGTATTGATAAAGGTGACATTCTGACCTCACACAGGCTCCCAGCAGTCTTATGGAG
ACGCTTGAACAGCATCTAAATACATTAGAAGGAAAGAAACCTGGAAACAATGAAGGATCTGGTGTCCC
TCTCCATTAAGTAAGTCTTCTCCAGCCACAACCTGTTACGTCTCCTAATTCTACACCAGCTAAAACCTATT
GACACATCCCCACCGTTGATTTATTTGCAACTGCATCTGCGGCTGTCCCAGTCAGCACTTCTAAACCA
TCTAGTGATCTCCTGGACCTCCAGCCAGACTTTTCTCTGGAGGGGCGAGCAGCAGCCGAGCACCAGCA
CCACCACCACCTGCTGGAGGAGCCACTGCATGGGAGACCTTTTGGGAGAGGATTCTTTGGCTGCACTT
TCCTCTGTTCCCTCTGAAGCACAGATTTTCAGATCCATTTGACCAGAACCTACCCCTCTACTACAACT
GCTGAAATTTGCAACTGCCTCAGCTTCTGCTCCACTACTACAACCTGTTACTGCTGCTACTGGAAGTG
GATCTCTTTGGAGATGCCTTTGCAGCTTCTCCTGGGAGGCCCTGCAGCATCCGAAGGGGCCCGCA
CCAGCTACCCCAACCCCTGTAGCAGCAGCACTGATGCATGTTTCAGGAAATGACCCCTTTGCCCCGCT
GAAGGTAGTGACAGGCTGCACCTGAGCTGGACCTCTTTGCAATGAAGCCACCTGAGACCAGTGTTCCT
GTAGTTACCCCTACAGCTAGCACAGCCCCTCCAGTTCCCGCAACTGCTCCTTCTCCTGCTCCTGCCGTT
GCAGCTGCTGCTGCTGCCACTACTGCTGCCACCGCGCTGCCACCACCTACCACCACCTCCGCTGCC
ACCGCCACCCTGCTCCTCCTGCTCTAGATATCTTTGGTGATTTATTTGAGTCCACTCCTGAAGTTGCT
GCAGCGCCTAAGCCAGATGCTGCTCCTAGCATAGACCTGTTTAGTACAGATGCTTTCTCCTCCACCA
CAAGGGGCTCTCCTGTGCTGAGAGTTCTCTCACTGCTGACCTTTATCTGTGGATGCATTTGCAGCA
CCATCTCCTGCAACCCTGCCTGCCAGCAAAGGTGGATTCTTCAGGTGCATAGACCTTTTGGGGAT
GCATTTGGAAGTAGTGCTTCTGAACCCCAACCTGCATCTCAGGCTGCTTCTAGTTTCATCAGCATCGCA
GACCTACTAGCTGGATTTGGGGTTCTTTCATGGCGCCTTCCCATCTCCAGTGACTCCAGCTCAGAAT
AACCTGCTACAGCCCAATTTTGGAGCAGCTTTTGGGACAACGCCTTCAACTTCCAGCAGCAGCTCCTTT
GATCCATCAGTGTGTTGATGGTCTAGGTGATCTTTGATGCCAACCATGGCACCAGCTGGGAGCCTGCA
CCTGTCTCAATGGTACCACCCAGTCTGCAATGGCAGCCAGCAAAGCCCTTGAAGTGTCTTGATTCA
TCTCTTGGCAGCTTAGTAGGCAATCTTGGAAATTTCTGGTACCACAACAAAAAGGGAGATCTTCAAGTG
AATGCTGGAGAGAAAAAGTTGACTGGTGGAGCCAACTGGCAGCCTAAAGTAGTCCAGAACCTGGTCA
GCAGGCGTTCCACCAAGTGCACCTTTGCAAGGAGCTGTACCTCCAACCAAGTTCAAGTCTCCTGTTGCC
GGGGCCCCATCGGTTGGACAACCTGGAGCAGGATTTGGAATGCCTCCTGCTGGGACAGGCATGCCATG
ATGCCCTAGCAGCCGGTATGTTTGACAGCCCATGATGAGGCCCCCTTTGGAGCTGCCCTGTACCT
GGCAGCAGCTTTCTCCAAGCCCTACACCTGCCAGTCAGAGTCCCAAGAAACCTCCAGCAAAGGCCCA
TTAGCGGATCTTAACATCAAGGATTTCTGTAA
  
```

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001242792  
**Insert Size:** 2724 bp

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_001242792.1</a></u>
<b>RefSeq Size:</b>	4462 bp
<b>RefSeq ORF:</b>	2724 bp
<b>Locus ID:</b>	9892
<b>UniProt ID:</b>	<u><a href="#">O60641</a></u>
<b>Cytogenetics:</b>	6q14.2
<b>MW:</b>	92.5 kDa
<b>Gene Summary:</b>	<p>Adaptins are components of the adapter complexes which link clathrin to receptors in coated vesicles. Clathrin-associated protein complexes are believed to interact with the cytoplasmic tails of membrane proteins, leading to their selection and concentration. Binding of AP180 to clathrin triskelia induces their assembly into 60-70 nm coats (By similarity).[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (3) has an alternate 5' UTR exon, compared to variant 1. Variants 1 and 3 encode the same isoform (a). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>