

## Product datasheet for **SC331863**

### AP180 (SNAP91) (NM\_001242794) Human Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** AP180 (SNAP91) (NM\_001242794) Human Untagged Clone  
**Tag:** Tag Free  
**Symbol:** AP180  
**Synonyms:** AP180; CALM  
**Vector:** pCMV6-Entry (PS100001)  
**Fully Sequenced ORF:** >SC331863 representing NM\_001242794.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGTCGGGCCAAACGCTCACGGATCGGATCGCCGCCGCTCAGTACAGCGTTACAGGCTCTGCTGTAGCA
AGAGCGGTCTGCAAAGCCACTACTCATGAAGTAATGGGCCCAAGAAAAAGCACCTGGACTATTTGATC
CAGGCTACCAACGAGACCAATGTTAATATTCCTCAGATGGCCGACACTCTCTTTGAGCGGCAACAAAC
AGTAGCTGGGTGGTTGTGTTAAGGCTTAGTGACAACACATCATCTCATGGTGCATGGAATGAGAGA
TTTATTCAATATTTGGCTTCTAGAAATACACTATTCATCTCAGCAATTTTTGGACAAAAGTGGATCC
CATGGTTATGATATGTCTACCTTCATAAGCGCTATAGTAGATATTTGAATGAAAAGGCTTTTTCTTAC
AGACAGATGGCCTTTGATTTGCCAGGGTGAAGAAAGGGGCCGATGGTGTAAATGAGGACAATGGCTCCC
GAAAAGCTGCTAAAGAGTATGCCAATACTACAGGGACAAATTGATGCACTGCTTGAATTTGATGTGCAT
CCAAATGAACTAACAAATGGTGTCAATAATGCAGCATTATGCTTCTTTCAAAGATCTTATCAAACCT
TTTGCTTGCTACAATGATGGTGTATTAACTTACTCGAAAAGTTTTTGAATGAAGAAAGGACAATGT
AAAGATGCTCTAGAAATTTACAAACGATTTCTAACTAGAATGACACGAGTGTCTGAATTTCTCAAGGT
GCAGAGGCTCCCAGCAGTCTTATGGAGACGCTTGAACAGCATCTAAATACATTAGAAGGAAAGAAACCT
GGAAACAAATCTGGTGTCCCTCTCCATTAAGTAAGTCTTCTCCAGCCACAACCTGTACGTCTCCTAAT
TCTACACCAGCTAAAATATTGACACATCCCACCGTTGATTTATTTGCAACTGCATCTGCGGCTGTC
CCAGTCAGCACTTCTAAACCATCTAGTGATCTCCTGGACCTCCAGCCAGACTTTTCTCTGGAGGGGCA
GCAGCAGCCGACGACCAGCACCACCACCTGCTGGAGGAGCCACTGCATGGGAGGATTTGGGGGT
TCTTTTATGGCGCTTCCCCTCTCCAGTGACTCCAGCTCAGAATAACCTGCTACAGCCCAATTTGAG
GCAGCTTTTGGACAACGCCTTCAACTCCAGCAGCAGCTCCTTTGATCCATCAGTGTGTTGATGGTCTA
GGTGTCTTTTGTGCAACCATGGCACCAGCTGGGCAGCCTGCACCTGTCTCAATGGTACCACCCAGT
CCTGCAATGGCAGCCAGCAAAGCCCTTGGAAAGTATCTTGATTCTCTTCCAGCTTAGTAGGCAAT
CTTGGAATTTCTGGTACCACAACAAAAAGGGAGATCTTCAAGTGAATGCTGGAGAGAAAAAGTTGACT
GGTGGAGCCAACGGCAGCCTAAAGTAGCTCCAGCAACCTGGTCCAGCAGCGTTCCACCAAGTGCACCT
TTGCAAGGAGCTGTACCTCCAACCAAGTTCAGTTCTCCTGTTGCCGGGGCCCATCGGTTGGACAACCT
GGAGCAGGATTTGGAATGCCTCCTGCTGGGACAGGCATGCCATGATGCCTCAGCAGCCGGTTCATGTTT
GCACAGCCATGATGAGGCCCTTTGGAGCTGCCGCTGTACCTGGCACGCAGCTTTCTCCAAGCCCT
ACACCTGCCAGTCAGAGTCCAAGAAACCTCCAGCAAAGGACCCATTAGCGGATCTTAACATCAAGGAT
TTCTTGTA
  
```



[View online »](#)

<b>Restriction Sites:</b>	Sgfl-Mlul
<b>ACCN:</b>	NM_001242794
<b>Insert Size:</b>	1803 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001242794.1</a></u>
<b>RefSeq Size:</b>	3340 bp
<b>RefSeq ORF:</b>	1803 bp
<b>Locus ID:</b>	9892
<b>UniProt ID:</b>	<u><a href="#">O60641</a></u>
<b>Cytogenetics:</b>	6q14.2
<b>MW:</b>	63.1 kDa
<b>Gene Summary:</b>	<p>Adaptins are components of the adapter complexes which link clathrin to receptors in coated vesicles. Clathrin-associated protein complexes are believed to interact with the cytoplasmic tails of membrane proteins, leading to their selection and concentration. Binding of AP180 to clathrin triskelia induces their assembly into 60-70 nm coats (By similarity).[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (5) lacks an internal segment in the 5' UTR and multiple exons in the CDS, compared to variant 1. The resulting isoform (c) lacks three internal segments, compared to isoform a.</p>