

## Product datasheet for **SC331867**

### ZNF790 (NM\_001242800) Human Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** ZNF790 (NM\_001242800) Human Untagged Clone  
**Tag:** Tag Free  
**Symbol:** ZNF790  
**Vector:** pCMV6-Entry (PS100001)  
**Fully Sequenced ORF:** >SC331867 representing NM\_001242800.  
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGGCCCATTTGATGATGTTTCAGGGATGTGGCTGTAGATTTCTCTCAGGAGGAGTGGGAGTGCCTGGAC
CTGGAACAGAGGGATTTATATAGAGATGTGATGTTGGAGAACTACAGCAACATGGTCTCACTGGGTTTT
TGCATTTATCAGCCAGAAGCGTTCTCTTTATTGGAGAAAGGAAAGAGCCCTGGAAGATTTTGAGGGAT
GAGACAAGAGGACCTTGCCAGACATGCAGTCGAGGTGTCAGACCAAGAAGTTATTACCAAAAAATGGC
ATTTTTGAGAGAGAAAATAGCCCAATTGGAATAATGAGAATTTGTAAAAACCACAGCCTTGACTGTTTA
GTTTTTAGAGGTGACTGGGAAGGCAACACTCAGTTTCAAACACTCAAGATAATCAAGAGGAATGCCTTC
AAGCAGGTGATACGCACCTGTGAAAAAAGGCCACTTTTAACCAGCATACAGTGTTAATCTACACCAG
AGACTTAATACAGGAGACAACTGAATGAATTTAAAGAACTGGGAAAAGCCTTTATTTCTGGTTCAGAT
CATACTCAACATCAGTTAATTCACACAAGTGAGAAATCTGTGGAGATAAAGAATGTGGGAATACCTTT
CTTCTGATTCAGAAGTTATTCAATATCAGACAGTTCACACTGTTAAGAAAAATATGAATGTAAAGAA
TGTGGGAAGTCTTTAGTTTACGTTTCGAGTCTTACTGGTCATAAGAGAATTCATACCGGTGAGAAACCT
TTTAAATGTAAGGATTGTGGAAAAGCCTTTAGATTTTCATTCACAACCTAGTGCCATAAGCGAATTCAT
ACTGGTGAGAAATCTTATGAATGTAAGGAATGTGGGAAGCCTTTAGTTGTGGCTCAGATCTTACTCGA
CATCAGAGAATTCATACTGGTGAAAAACCTATGAATGTAATGAATGTAGAAAGCCTTTAGTCAGCGA
TCACATCTTATTAACATCAGAGAATTCACACTGGTGAAAAACCTTATGAATGTAAGGAGTGTGGGAAA
GCTTTTACTCGTGGATCACACCTAACTCAGCATCAGAGAATTCATACTGGTGAGAAATCTCATGAGTGT
AAGGAATGTGGAAAAGCCTTTATTCGTGGTTCAAATCTTGCTCAACATCAGAATGTTTCATGTTGGTAGG
AAACCTTATAAATGTGAGAAATGTGGAAAAGCCTATATTTGGAGCTCACACCTTGCTCGACATCAGCGA
ATTCATACTGGCAGGAAACCTTATGAATGTAAGCAATGCGGGAAGACTTTTACTTGGGCTTCGTATCTT
GCTCAACATGAGAAAATTCACAATGAGAGGAAATCCTATGAATGTAAGGAATGTGGAAAAGACCTTTCTT
CATGGCTCAGAGTTTAAATCGACATCAGAAAATTCATACTGGTGAGAGAACTATGAATGTAAGGAATGT
GGAAAGACCTTTTTTCGTGGTTTCAAGAACTTAAATCGACACCAGAAAATTCATACTGGAAAAGGCCATAT
GAATGTGAAGAATGTGGAAAAGCCTTTCTCTGGGTTTCAACAACCTTACTCGACATCAGAGAATGCATACT
GGTGAGGAAACCTTACGTATGTAAGAATGTGGGAAATCTTTTATCTGGGTTTACAGCTTACACGACAT
AAGAAAATTCATACTGATGCAGAACCTTATGGATGCAAGAAAAGTAGCCACATCTTTAGTACCATTCA
TATTTTACTGAACAAAAAATTCATAATAGTGCAAACTCTGTGAATGGACAGACTATGGGAACACCTTT
AGTCATGAGTCAAACCTTGTCAACACCAGAATATTTACTCTTTGAGAAATCCTATGAATTTAAAGAT
TTTGAGAAAGCATTCTTCAAGCTCTCACTTCATTTCACTCTTGTGA
```

**Restriction Sites:** Sgfl-Mlul



<b>ACCN:</b>	NM_001242800
<b>Insert Size:</b>	1911 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_001242800.1</a></u>
<b>RefSeq Size:</b>	2183 bp
<b>RefSeq ORF:</b>	1911 bp
<b>Locus ID:</b>	388536
<b>UniProt ID:</b>	<u><a href="#">Q6PG37</a></u>
<b>Cytogenetics:</b>	19q13.12
<b>MW:</b>	74.6 kDa
<b>Gene Summary:</b>	<p>May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (2) represents a distinct first exon and 5' UTR.</p>