

## Product datasheet for **SC331875**

### **BRDT (NM\_001242810) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	BRDT (NM_001242810) Human Untagged Clone
Tag:	Tag Free
Symbol:	BRDT
Synonyms:	BRD6; CT9; SPGF21
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

**Fully Sequenced ORF:** >SC331875 representing NM\_001242810.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGATTTAAATACAATTAAGAAGCGCTTGGAGAATAAATATTATGCGAAGGCTTCAGAATGTATAGAA
GACTTCAATACAATGTTCTCAAATTGTTATTTATATAACAAGCCTGGAGATGACATTGTTCTTATGGCA
CAAGCTCTAGAGAAGCTGTTTATGCAGAAATTATCTCAGATGCCACAAGAAGAGCAAGTTGTGGTGTT
AAGGAAAGAATCAAGAAAGGCACCTCAACAGAATATAGCTGTTTCTCTGCTAAAGAAAAATCATCACCC
AGCGCAACAGAAAAAGTATTTAAGCAGCAAGAAATTCCTTCTGTATTTCTAAGACATCTATTTCTCCC
TTGAACGTGGTACAGGGAGCTTCAGTCAACTCCAGTTCACAAACTGCGGCCCAAGTTACAAAAGGTGTG
AAGAGGAAAGCAGATACAACAACCTCTGCAACTTCAGCAGTTAAAGCAAGTAGTGAATTTTCTCCAACA
TTCACAGAAAAATCAGTGGCACTGCCACCTATAAAAGAAAAATAGCCAAAGAATGTTTTGCCAGATTCT
CAGCAACAATATAATGTTGTGAAGACTGTTAAAGTAACTGAACAATTAAGGCACTGTAGTGAGATTCTT
AAAGAAATGCTTGCAAAGAAACATTTTTTCATATGCATGGCCCTTTATAATCCTGTTGACGTTAATGCT
TTGGGACTCCATAACTACTATGACGTTGTCAAAAATCCGATGGATCTTGGAACTATTAAGGAGAAAAATG
GATAACCAAGAATATAAGGATGCATACAAATTTGCGGCAGATGTTAGATTAATGTTTCATGAATTGCTAC
AAGTACAATCCTCCAGATCACGAAGTTGTGACAATGGCAAGAATGCTTCAGGATGTTTTCGAAACGCAT
TTTTCAAAGATCCCGATTGAACCTGTTGAGAGTATGCCTTTATGTTACATCAAAAACAGATATCACAGAA
ACCACTGGTAGAGAGAACAATAAGGCTCCTCTGAAGGGAACCTTCTGATGATTCTGAAGATGAG
CGAGTTAAGCGTCTTGCAAAGCTTCAGGAGCAGCTTAAAGCTGTACATCAACAGCTCCAGGTTTTGTCC
CAAGTACCTTTCCGTAAGCTAAATAAAAAGAAAGAGAAGTCTAAAAGGAAAAAGAAAAAGAAAGGTT
AATAACAGCAATGAAAAATCAAGAAAAATGTGTGAGCAATGAGGCTAAAGGAAAAAGTCCAAGAGAAAT
CAGCCAAAAGAAAGGAAACAACAGTTTATTGGTCTAAAATCTGAAGATGAAGATAATGCTAAACCTATG
AACTATGATGAGAAAAGGCAGTTAAGTCTGAATATAAACAACCTCCCTGGAGATAAACTGGGCGAGTA
GTTACATAATACAATCAAGAGAGCCTTCTGAGCAATTCCTGATGAGATAGAGATAGACTTTT
GAAACTGAAAGCATCAACTAAGAGAATTAGAAAAATATGTTTCGGCATGTCTAAGAAAGAGACCA
TTAAAACCTCCTGCTAAGAAAAATATGATGTCCAAGAAGAACCTCACTCACAGAAAAACAGGAATTG
GAAAAGCGGTTACTGGATGTTAATAATCAGTTAAATTCTAGAAAACGTCAACAAAAATCTGATAAAACG
CAACCATCCAAGCTGTTGAAAATGTTTCCGACTGAGTGAGAGCAGCAGCAGCAGCAGCTCATCA
GAGTCTGAAAGTAGCAGCAGTACTTAAGCTCTCAGACAGCAGTATTCTGAATCAGAAATGTTCCCT
AAGTTTACAGAAGTAAAACCAATGATTCTCCTTCTAAAGAGAATGTAAGAAAAATGAAGAATGAATGC
ATACTGCCTGAAGGAAGAACAGGCGTCACACAGATAGGATATTGTGTGCAAGACACAACCTCTGCCAAT
ACTACCCTTGTTTCATCAGACCACACCTTACATGTAATGCCACCAATCACCACCAATTAGCATTAAAT
TATCAAGAATTAGAACATTTACAGACTGTGAAAAACATTTACCTTTACAATTTGCTCCCTCAGGT
GATTCTGAACAGCTCTCAAAATGGCATAACTGTGATGCATCCATCTGGTGTAGTGACACAACGATGTTA
GAATCTGAATGTCAAGCTCCTGTACAGAAGGATATAAAGATTAAGAATGCAGATTGCAAGAAAGTTTA
GGCAAACCAAGTAAACCATCAGGTGTAATGAAATCCTCAGATGAGCTCTTCAACCAATTTAGAAAAGCA
GCCATAGAAAAGGAAGTAAAAGCTCGGACACAGGAACCTATACGGAAGCATTGGAAACAAAATACAAAG
GAACTAAAAGCATCTCAAGAAAAATCAGAGGGATCTTGGGAATGGATTGACTGTAGAATCTTTTTCAAT
AAAATACAAAACAAGTCTCTGGAGAAGAGCAGAAAGAACATCAGCAGTCATCAGAAGCTCAAGATAAA
TCCAAACTCTGGCTTCTCAAGACCGTGATTTAGCAAGGCAGAAAGAACAAGAGAGGAGGAGAGAGAA
GCAATGGTGGTACCATTGATATGACCCTTCAAAGTGACATTATGACAATGTTTGAACAACTTTGAT
TAA
  
```

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001242810  
**Insert Size:** 2625 bp

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001242810.1</a>
<b>RefSeq Size:</b>	3717 bp
<b>RefSeq ORF:</b>	2625 bp
<b>Locus ID:</b>	676
<b>UniProt ID:</b>	<a href="#">Q58F21</a>
<b>Cytogenetics:</b>	1p22.1
<b>Protein Families:</b>	Protein Kinase, Transcription Factors
<b>MW:</b>	99.4 kDa
<b>Gene Summary:</b>	<p>BRDT is similar to the RING3 protein family. It possesses 2 bromodomain motifs and a PEST sequence (a cluster of proline, glutamic acid, serine, and threonine residues), characteristic of proteins that undergo rapid intracellular degradation. The bromodomain is found in proteins that regulate transcription. Several transcript variants encoding multiple isoforms have been found for this gene. [provided by RefSeq, Jun 2011]</p> <p>Transcript Variant: This variant (7) differs in the 5' UTR and coding sequence compared to variant 4. The resulting isoform (d) is shorter at the N-terminus and lacks a short internal segment compared to isoform a.</p>