

# Product datasheet for SC331889

### ELOVL5 (NM\_001242830) Human Untagged Clone

### **Product data:**

#### **Product Type: Expression Plasmids Product Name:** ELOVL5 (NM 001242830) Human Untagged Clone Tag: Tag Free ELOVL5 Symbol: Synonyms: dJ483K16.1; HELO1; SCA38 Vector: pCMV6-Entry (PS100001) >SC331889 representing NM\_001242830. **Fully Sequenced ORF:** Blue=Insert sequence Red=Cloning site Green=Tag(s) ATGGAACATTTTGATGCATCACTTAGTACCTATTTCAAGGCATTGCTAGGCCCTCGAGATACTAGAGTA AAAGGATGGTTTCTTCTGGACAATTATATACCCACATTTATCTGCTCTGTCATATATTTACTAATTGTA TGGCTGGGACCAAAATACATGAGGAATAAACAGCCATTCTCTTGCCGGGGGATTTTAGTGGTGTATAAC CTTGGACTCACACTGCTGTCTCTGTATATGTTCTGTGAGTTAGTAACAGGAGTATGGGAAGGCAAATAC AACTTCTTCTGTCAGGGCACACGCACCGCAGGAGAATCAGATATGAAGATTATCCGTGTCCTCTGGTGG TACTACTTCTCCAAACTCATAGAATTTATGGACACTTTCTTCTTCATCCTGCGCAAGAACAACCACCAG ATCACGGTCCTGCACGTCTACCACCATGCCTCGATGCTGAACATCTGGTGGTTTGTGATGAACTGGGTC CCCTGCGGCCACTCTTCAGTTTGTGCTGACAATCATCCAGACCAGCTGCGGGGTCATCTGGCCGTGCAC ATTCCCTCTTGGTTGGTTGTATTTCCAGATTGGATACATGATTTCCCTGATTGCTCTCTCACAAACTT CTACATTCAGACCTACAACAAGAAAGGGGCCTCCCGAAGGAAAGACCACCTGAAGGACCACCAGAATGG GTCCATGGCTGCTGTGAATGGACACACCAACAGCTTTTCACCCCTGGAAAACAATGTGAAGCCAAGGAA GCTGCGGAAGGATTGAAGTCAAAGAATTGA **Restriction Sites:** Sgfl-Mlul ACCN: NM 001242830 Insert Size: 789 bp **OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP). The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube **Components:** containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



View online »

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## Science ELOVL5 (NM\_001242830) Human Untagged Clone – SC331889

| Reconstitution Method: | <ol> <li>Centrifuge at 5,000xg for 5min.</li> <li>Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>Close the tube and incubate for 10 minutes at room temperature.</li> <li>Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>   |
|------------------------|--|
| RefSeq:                | <u>NM 001242830.1</u>  |
| RefSeq Size:           | 2913 bp  |
| RefSeq ORF:            | 789 bp   |
| Locus ID:              | 60481  |
| UniProt ID:            | Q9NYP7   |
| Cytogenetics:          | 6p12.1   |
| Protein Families:      | Transmembrane  |
| Protein Pathways:      | Biosynthesis of unsaturated fatty acids  |
| MW:                    | 30.7 kDa   |
| Gene Summary:          | This gene belongs to the ELO family. It is highly expressed in the adrenal gland and testis,<br>and encodes a multi-pass membrane protein that is localized in the endoplasmic reticulum.<br>This protein is involved in the elongation of long-chain polyunsaturated fatty acids. Mutations<br>in this gene have been associated with spinocerebellar ataxia-38 (SCA38). Alternatively spliced<br>transcript variants have been found for this gene. [provided by RefSeq, Sep 2014]<br>Transcript Variant: This variant (3) lacks an internal coding exon, which causes a frame-shift<br>compared to variant 1. The resulting isoform (3) is shorter with a distinct C-terminus<br>compared to isoform 1. |

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US