

Product datasheet for SC331911

CARKD (NAXD) (NM_001242882) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CARKD (NAXD) (NM_001242882) Human Untagged Clone
Tag:	Tag Free
Symbol:	CARKD
Synonyms:	CARKD; LP3298; PEBEL2
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC331911 representing NM_001242882. Blue=Insert sequence Red=Cloning site Green=Tag(s)

ATGGCCCTGGGTCTCGCTGTGGGGCAATCCGGGCTTGCAGACGAGTTTTAGAAAGAGCGTTTTCGCTA
CGTAAAGCACATTCGATAAAGGATATGGAATACTTTGCAGCTGGTGAGAAATATCATACCTCCTCTG
TCTTCCACAAAGCACAAAGGGCAAGATGGAAGAATAGGCGTAGTTGGAGGCTGTCAGGAGTACACTGGA
GCCCCATATTTTGCAGCAATCTCAGCTCTCAAAGTGGGCGCAGACTTGTCACCGTGTCTGTGCCAGT
GCGGCCGCACCTGTGATTAAGGCCTACAGCCCGGAGCTGATCGTCCACCCAGTTCTTGACAGCCCCAAT
GCTGTTTCATGAGGTGAGAAGTGGCTGCCCGGCTGCATGCTCTTGTCGTAGGACCTGGCTTGGGTAGA
GATGATGCGCTTCTCAGAAATGTCCAGGGCATTGGAAGTGTCAAAGGCCAGGGACATCCCTGTTGTC
ATCGACGCGGATGGCCTGTGGCTGGTCGCTCAGCAGCCGGCCCTCATCCATGGCTACCGGAAGGCTGTG
CTCACTCCCAACCACGTGGAGTTCAGCAGACTGTATGACGCTGTGCTCAGAGGCCCTATGGACAGCGAT
GACAGCCATGGATCTGTGCTAAGACTCAGCCAAGCCCTGGGCAACGTGACGGTGGTCCAGAAAGGAGAG
CGCGACATCCTCTCCAACGGCCAGCAGGTGCTTGTGTGCAGCCAGGAAGGCAGCAGCCGAGGTGTGGA
GGGCAAGGGGACCTCCTGTCGGGCTCCCTGGGCGTCTGGTACACTGGGCGCTCCTTGCTGGACCACAG
AAAACAAATGGGTCCAGCCCTCTCCTGGTGGCCGCGTTTGGCGCCTGCTCTCTACCAGGCAGTGCAAC
CACCAAGCCTTCCAGAAGCACGGTCGCTCCACCACCACCTCCGACATGATCGCCGAGGTGGGGGCCGCG
TTCAGCAAGCTCTTTGAAACCTGA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001242882
Insert Size:	990 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001242882.1</u>
RefSeq Size:	2617 bp
RefSeq ORF:	990 bp
Locus ID:	55739
UniProt ID:	<u>Q8IW45</u>
Cytogenetics:	13q34
Protein Families:	Transmembrane
MW:	35.1 kDa
Gene Summary:	Catalyzes the dehydration of the S-form of NAD(P)HX at the expense of ATP, which is converted to ADP. Together with NAD(P)HX epimerase, which catalyzes the epimerization of the S- and R-forms, the enzyme allows the repair of both epimers of NAD(P)HX, a damaged form of NAD(P)H that is a result of enzymatic or heat-dependent hydration.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (3) differs in the 5' UTR and 5' coding region, uses an alternate start codon, and uses an alternate splice site that causes a frameshift in the 3' coding region, compared to variant 1. The encoded isoform (c) has distinct N- and C-termini and is shorter than isoform a.