

Product datasheet for SC331917

ZKSCAN3 (NM_001242894) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: ZKSCAN3 (NM_001242894) Human Untagged Clone
Tag: Tag Free
Symbol: ZKSCAN3
Synonyms: dj874C20.1; dj874C20.1.; ZF47; zfp-47; Zfp47; ZFP306; ZNF306; ZNF309; ZSCAN13; ZSCAN35
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC331917 representing NM_001242894.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGCTAGAGAATTAAGTAAAGCACAGCCCTGGATGCCAGTCTACAGAAGACCAGATGGAGCTTCTG
GTCATAAAGGTGGAGGAAGAAGAAGCCGGTTTTCCAGTAGCCAGATCTGGGTTCTGAGGGCTCCCGC
GAGCGCTCCGAGGCTCCGCTACCCGGAGGCTGCAGGCCCGCGAGGCGCTGAGTCGGCTCCGAGAG
CTCTGCCGACAGTGGCTGCAGCCTGAGATGCACAGCAAGGAGCAGATCCTGGAGCTGCTGGTCTGGAG
CAGTTCCTGACCATCCTGCCGGGAATCTGCAGAGCTGGGTGCCGGAGCAGCATCCAGAGAGCGGGGAG
GAGGTGGTGGTGTATTGGAGTATTTGGAGAGGACAGTGGATGAGCCGGCGCCGAGTTTCAGGTGTT
GACCAGGGGCAAGAACTGCTCTGTTGCAAGATGGCACTATTGACACCAGCCCAGGGTCACAAAGTAGC
CAATTTACAGTAATGAAGGCTCTGCTCAAGCATGAATCTGTGGGATCCCAGCCTTACAAGATAGAGTT
CTCCAGGTCCCGTGCTTGCCCATGGAGGATGCTGCAGAGAAGATAAAGTGGTAGCTTCTAGGCTTACT
CCAGAGTCCCAGGGGTTGTTGAAAGTGGAAAGATGTGGCCCTGACCCTCACCCCTGAATGGACACAGCAG
GATTCATCTCAGGGGAATCTCTGTAGAGATGAAAAGCAGGAGAACCATGGCAGCCTGGTCTCCCTGGGT
GATGAAAAACAGACTAAGAGCAGGGACTTGCCCTCAGCTGAGGAGCTCCAGAAAAGGAGCATGGGAAG
ATATCGTGCCACCTGAGAGAAGACATTGCCAGATTCCACATGTGCAGAAGCTGGTGAACAGGAGGGC
AGGCTACAAAGAAAGCAGAAAAATGCCACAGGAGGGAGGCGGCACATCTGCCATGAATGTGGAAAGAGT
TTTGCTCAAAGCTCAGGCCTGAGTAAACACAGGAGAATCCACACTGGTGAGAAAACCCTACGAATGTGAA
GAGTGTGGCAAAGCCTTCATTGGGAGCTCTGCCCTTGTCATTCAGAGAGTCCACACTGGTGAGAAG
CCATATGAGTGTGAAGAATGTGGTAAGGCCTTCAGTCATAGCTCAGACCTTATCAAGCATCAGAGAACC
CACACTGGGGAGAAGCCCTATGAGTGTGATGACTGTGGGAAGACCTTCAGCCAGAGCTGCAGCCTCCTT
GAACATCACAGAATCCACACTGGGGAGAAGCCGTATCAGTGCAGTATGTGTGGCAAAGCCTTTAGGCGA
AGTTCACATCTCCTGAGACATCAGAGGATCCATACTGGGGATAAAAAATGTTCAGGAACCTGAGCAGGGA
GAGGCCTGAAAAGTAGGATGAAAAGCCAGTTGAAAAATGTTGAAACTCCCATGTCTTATAAATGTAAT
GAGTGTGAAAAGAGTTTCACTCAGAATACAGGCCTCATTGAACATCAAAAAATCCACACTGGTGAGAAA
CCCTATCAGTGAATGCGTGTGGAAAAGGCTTCACCCGAATTTACATACCTTGTTCAACATCAGAGAAGC
CATGTAGGAAAAAACATCCTATCACAGTGA
  
```

Restriction Sites: SgfI-MluI
ACCN: NM_001242894



[View online »](#)

Insert Size:	1617 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001242894.1</u>
RefSeq Size:	4842 bp
RefSeq ORF:	1617 bp
Locus ID:	80317
UniProt ID:	<u>Q9BRR0</u>
Cytogenetics:	6p22.1
Protein Families:	Transcription Factors
MW:	60.6 kDa
Gene Summary:	<p>Transcriptional factor that binds to the consensus sequence 5'-[GT][AG][AGT]GGGG-3' and acts as a repressor of autophagy. Specifically represses expression of genes involved in autophagy and lysosome biogenesis/function such as MAP1LC3B, ULK1 or WIPI2. Associates with chromatin at the ITGB4 and VEGF promoters. Also acts as a transcription activator and promotes cancer cell progression and/or migration in various tumors and myelomas. [UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) represents the longest transcript and encodes the longer isoform (1). Variants 1 and 2 encode the same protein. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>