

Product datasheet for **SC331940**

ABCB9 (NM_001243014) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ABCB9 (NM_001243014) Human Untagged Clone
Tag:	Tag Free
Symbol:	ABCB9
Synonyms:	EST122234; TAPL
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

Fully Sequenced ORF: >SC331940 representing NM_001243014.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGCGGCTGTGGAAGGCGGTGGTGGTACTTTGGCCTTCATGAGTGTGGACATCTGCGTGACCACGGCC
ATCTATGTCTTCAGCCACCTGGACCGCAGCCTCCTGGAGGACATCCGCCACTTCAACATCTTTGACTCG
GTGCTGGATCTCTGGCAGCCTGCCTGTACCGCAGCTGCCTGCTGCTGGGAGCCACCATTGGTGTGGCC
AAGAACAGTGCCTGGGCCCCGGCGGCTGCGGGCCTCGTGGTGCATCACCCTCGTGCCTCTTC
GTGGCATCTATGCCATGGTGAAGCTGCTGCTTCTCAGAGGTGCGCAGGCCCATCCGGGACCCCTGG
TTTTGGCCCTGTTGCTGTGGACGTACATTTCACTCGGCCATCCTTCTGCTCTGGTGGCTGCTGTCC
ACCGTGCAGCCAGGCACCCAGGCCCTGGAGCCAGGGGCGGCCACCGAGGCTGAGGGCTTCCCTGGGAGC
GGCCGGCCACCGCCGAGCAGGCGTCTGGGGCCACGCTGCAGAAGCTGCTCTCTACACCAAGCCCGAC
GTGGCCTTCTCGTGGCCGCTCCTTCTCCTCATCGTGGCAGCTCTGGGAGAGACCTTCTGCCCTAC
TACACGGCCCGCCATTGATGGCATCGTCATCCAGAAAAGCATGGATCAGTTCAGCACGGCTGTCGTC
ATCGTGTGCTGCTGGCCATTGGCAGCTCATTTGCCGAGGATTCGGGGCGGCATTTTACCCTCATA
TTTGCCAGACTGAACATTCGCCTTCGAACTGTCTCTCCGCTCACTGGTGTCCCAGGAGACAAGCTTC
TTTGATGAGAACCACAGGGGACCTCATCTCCGCTGACCTCGGACACCACCATGGTACAGCAGCTG
GTCTCCAGAACATCAATGTCTTCTGCGGAACACAGTCAAGGTACGGGCGTGGTGGTCTTCATGTTC
AGCCTCTCATGGCAGCTCTCCTTGGTACCTTCATGGGCTTCCCACATCATCATGATGGTGTCCAACATC
TACGGCAAGTACTACAAGAGGCTCTCCAAGAGGTCCAGAATGCCCTGGCCAGAGCGAGCAACACGGCG
GAGGAGACCATCAGTGCCATGAAGACTGTCCGGAGCTTCGCAATGAGGAGGAGGAGGCAGAGGTGTAC
CTGCGGAAGTGCAGCAGGTGTACAAGCTGAACAGGAAGGAGGAGCAGCTGCCTACATGTACTACGTCTGG
GGCAGCGGGCTCACACTGCTGGTGGTCCAGGTGAGTCCCTCTACTACGGGGCCACCTTGTGCATCTCA
GGCAGATGACCAGCGGCAACCTCATCGCTTCATCATCTACGAGTTTGTCTGGGAGATTGATGGAG
TCCGTGGGCTCCGTCTACAGTGGCTGATGCAGGAGTGGGGGCTGCTGAGAAGGTGTTGAGTTTATC
GACCGGACGCGACCATGGTGCACGATGGCAGCTTGGCCCCGACCACCTGGAGGGCCGGGTGGACTTT
GAGAATGTGACCTTACCTACCGCACTCGGCCCCACACCAGGTCTGCAGAATGTCTCCTTCAGCCTG
TCCCCCGCAAGGTGACGGCCCTGGTGGGGCCCTCGGGCAGTGGGAAGAGCTCCTGTGTCAACATCTG
GAGAACTTCTACCCCTGGAGGGGGCCGGGTGCTGCTGGACGGCAAGCCATCAGCGCTACGACCAC
AAGTACTGCACCGTGTGATCTCCTGGTGGAGCCAGGAGCCGCTGCTGTTGCCCCCTCCATCACGGAT
AACATCTCTACGGCTGCCACTGTGCCTTTCGAGATGGTGGTGGAGGCCGACAGAAGGCCAATGCC
CACGGCTTCATCATGGAATCCAGGACGGCTACAGCACAGAGACAGGGGAGAAGGGCGCCAGCTGTCA
GGTGGCCAGAAGCAGCGGTGGCCATGGCCCGGCTCTGGTGGGAACCCCCAGTCTCATCTGGAT
GAAGCCACCAGCGCTTTGGATGCCGAGAGCGAGTATCTGATCTGA
  
```

Restriction Sites: Sgfl-Mlul

ACCN: NM_001243014

Insert Size: 2046 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001243014.1](#)

RefSeq Size: 2443 bp

RefSeq ORF: 2046 bp

Locus ID: 23457

UniProt ID: [Q9NP78](#)

Cytogenetics: 12q24.31

Protein Families: Druggable Genome, Transmembrane

Protein Pathways: ABC transporters, Lysosome

MW: 75.3 kDa

Gene Summary: The membrane-associated protein encoded by this gene is a member of the superfamily of ATP-binding cassette (ABC) transporters. ABC proteins transport various molecules across extra- and intra-cellular membranes. ABC genes are divided into seven distinct subfamilies (ABC1, MDR/TAP, MRP, ALD, OABP, GCN20, White). This protein is a member of the MDR/TAP subfamily. Members of the MDR/TAP subfamily are involved in multidrug resistance as well as antigen presentation. This family member functions in the translocation of peptides from the cytosol into the lysosomal lumen. Alternative splicing of this gene results in distinct isoforms which are likely to have different substrate specificities. [provided by RefSeq, Jul 2011]
Transcript Variant: This variant (5, also known as 12C) has an alternate 3' exon and thus differs in the 3' coding region and 3' UTR, compared to variant 1. The encoded isoform (5) has a distinct and shorter C-terminus, compared to isoform 1.