

Product datasheet for **SC331946**

FYB1 (NM_001243093) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FYB1 (NM_001243093) Human Untagged Clone
Tag:	Tag Free
Symbol:	FYB1
Synonyms:	ADAP; FYB; PRO0823; SLAP-130; SLAP130; THC3
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

Fully Sequenced ORF: >SC331946 representing NM_001243093.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGATGGAAAGGCAGATGTAAGTCCCTCATGGCGAAATATAACACGGGGGGCAACCCGACAGAGGAT
GTCTCAGTCAATAGCCGACCCTTCAGAGTCACAGGGCCAACTCATCTTCAGGAATACAAGCAAGAAAG
AACTTATTCAACAACCAAGGAAATGCCAGCCCTCCTGCAGGACCCAGCAATGTACCTAAGTTTGGGTCC
CCAAAGCCACCTGTGGCAGTCAAACCTTCTTCTGAGGAAAAGCCTGACAAGGAACCAAGCCCCCGTTT
CTAAAGCCCACTGGAGCAGGCCAAAGATTCGGAACACCAGCCAGCTTGACCACCAGAGACCCCGAGGCG
AAAGTGGGATTTCTGAAACCTGTAGGCCCAAGCCATCAACTTGCCCAAAGAAGATTCCAAACCTACA
TTTCCCTGGCCTCCTGGAAACAAGCCATCTCTTCACAGTGTAACCAAGACCATGACTTAAAGCCACTA
GGCCCGAAATCTGGGCTACTCTCCAACCTCAGAAAAAAGAACAAGCAAGCGTTTCCCAAATTGACT
GGGGTTAAAGGGAAATTTATGTCAGCATACAAGATCTTGAACCAAGCCCTCTTCCCAAACCCGCC
TTTGGCCAGAAGCCGCCCTAAGTACCGAGAATCCCATGAAGACGAAAGCCCCATGAAGAATGTGTCT
TCATCAAAGGGTCCCAGCTCCCCTGGGAGTCAGGTCCAAAAGCGGCCCTTTAAAACAGCAAGGGAA
GACTCAGAAAATAAAGACCATGCAGGGGAGATTTCAAGTTTGCCTTTCTGGAGTGGTTTTGAAACCT
GCTGCGAGCAGGGGAGGCCAGGTCTCCAAAAATGGTGAAGAAAAAAGGAAGATAGGAAGATAGAT
GCTGCTAAGAACACCTTCCAGAGCAAAAATAATCAGGAAGAGTTGGCCTCAGGGACTCCTCTGCCAGG
TTCCCTAAGGCCCTTCTAAGCTGACAGTGGGGGGCCATGGGGCCAAAGTCAGGAAAAGGAAAAGGGA
GACAAGAATTCAGCCACCCCGAAACAGAAGCCATTGCCTCCCTTGTTTACCTTGGGTCCACCTCCACCA
AAACCCACAGACCACCAAAATGTTGACCTGACGAAATTCACAAAACCTTCTGAAAACAGTACTAGC
AAAGGCCAGAGCTTACTCAACAATTCCTGCCACCACCTCCACCATCCCATCCGCCAGCCAACCA
CCATTGCCAGCATCTCACCCATCACAACCACAGTCCCAAGCCTACCTCCAGAAAACATTAACCTCCG
CTAGATGAGGAACAAGACAGTGAAGGAGAAACATATGAAGACATAGAAGCATCCAAAGAAAGAGAGAAG
AAAAGGGAAAAGGAAGAAAAGAGAGGTTAGAGCTGGAGAAAAGGAACGAAAGAGAAAAGAAAAGAAA
GAACAAGAAAATAAAGAAGAAATTTAAACTAACAGGCCCTATTCAAGTCATCCATCTTGCAAAAGCTTGT
TGTGATGTCAAAGGAGGAAAGAATGAACTGAGCTTCAAGCAAGGAGAGCAAAATGAAATCATCCGCATC
ACAGACAACCCAGAAGGAAAATGGTTGGGCAGAACAGCAAGGGGTTCATATGGCTATATTAACAACCT
GCTGTAGAGATTGACTATGATCTTTGAAACTGAAAAAAGACTCTTGGTGCCCTTCAAGACCTATT
GAAGATGACCAAGAAGTATATGATGATGTTGCAGAGCAGGATGATATTAGCAGCCACAGTCAGAGTGA
AGTGGAGGGATATCCCTCCACCACCAGATGATGACATTTATGATGGGATTGAAGAGGAAGATGCTGAT
GATGGCTCCACACTACAGGTTCAAGAGAAGAGTAAACGTTGGTCCCTGGGGGATTTTGAAGATGTTAAG
GGAAAAGATGACAGAAAGAAAAGTATACGAGAGAAACCTAAAGTCTCTGACTCAGACAATAATGAAGGT
TCATCTTCCCTGCTCCTCCTAAACAATTGGACATGGGAGATGAAGTTTACGATGATGTGGATACCTCT
GATTTCCCTGTTTCATCAGCAGAGATGAGTCAAGGAACTAATGTTGGAAAAGCTAAGACAGAAAGAAAAG
GACCTTAAAGAGCTAAAAAGCAGGAAAAAGAAAAGAAAGAAAGACTTCAGGAAAAAATTTAAATATGATGGT
GAAATTAGAGTCTATATTCAACTAAAGTTACAACCTCCATAACTTCTAAAAAGTGGGGAACCAGAGAT
CTACAGGTAACCTGGTGAATCTCTAGAAGTTATACAAACCACAGATGACACAAAAGTTCTCTGCAGA
AATGAAGAAGGGAAATATGGTTATGTCCTTCGGAGTTACCTAGCCGACAATGATGGAGAGATCTATGAT
GATATTGCTGATGGCTGCATCTATGACAATGACTAG
  
```

Restriction Sites: Sgfl-Mlul
ACCN: NM_001243093
Insert Size: 2520 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001243093.1
RefSeq Size:	4811 bp
RefSeq ORF:	2520 bp
Locus ID:	2533
UniProt ID:	O15117
Cytogenetics:	5p13.1
Protein Families:	Druggable Genome
MW:	91.7 kDa
Gene Summary:	<p>The protein encoded by this gene is an adapter for the FYN protein and LCP2 signaling cascades in T-cells. The encoded protein is involved in platelet activation and controls the expression of interleukin-2. Three transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2011]</p> <p>Transcript Variant: This variant (3) encodes the longest isoform (3). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>