

Product datasheet for **SC331987**

ACSF3 (NM_001243279) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: ACSF3 (NM_001243279) Human Untagged Clone
Tag: Tag Free
Symbol: ACSF3
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC331987 representing NM_001243279.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGCTGCCCCATGTGGTCTCACCTTCCGGCGCCTGGGCTGCGCCTTGGCGTCTGCCGGCTGGCGCCT
GCGAGACACAGAGGAAGTGGTCTTCTGCACACAGCCCCAGTGGCCCGCTCGGACAGGAGCGCCCCGGTG
TTCACCCGTGCCCTGGCCTTTGGGGACAGAATCGCCCTGGTTGACCAGCACGGCCGCCACACGTACAGG
GAGCTTTATCCCGCAGCCTTCGCTGTCCAGGAGATCTGCAGGCTCTGCGGGTGTGTGCGGGGGAC
TCCGGGAGGAGAGGGTCTCCTTCTATGCGCTAACGATGCCTCCTACGTCGTGGCCACGTGGCGTCA
TGGATGAGTGGCGGTGTGGCAGTCCCCCTACAGGAAGCATCCCGCGGCCAGTGGAGTATGTCATC
TGGCACTCCAGAGCTCTGTGGTCTTGGCAGCCAGGAGTACCTGGAGCTCTGAGCCCGTGGTCAAGG
AAGCTGGGGTCCCGCTGCTGCCCTCACACCAGCCATCTACACTGGAGCAGTAGAGGAACCGGCAGAG
GTCCCGTCCCAGAGCAGGGATGGAGGAACAAGGGCGCCATGATCATCTACACCAGTGGGACCACGGGG
AGGCCAAAGGGCGTGTGAGCACGCACCAAAACATCAGGGCTGTGGTACCAGGGTGGTCCACAAGTGG
GCATGGACCAAGACGACGTGATCCTCCACGTGCTCCCGTGCACCACGTCCATGGTGTGGTCAACGCG
CTGCTCTGCTCTCTGGGTGGGAGCCACCTGTGTGATGATGCCTGAGTTCAGCCCTCAGCAGGTTTGG
GAAAAGTCTTAAGTTCTGAAACGCCGGATCAATGTCTTTATGGCAGTGCCTACAATATACCCAAG
CTGATGGAGTACTACGACAGGCATTTTACCAGCCGCACGCCAGGATTTCTTGGTGCAGTTTGTGAA
GAAAAAATTAGGCTGATGGTCTCAGGCTCAGCTGCCCTGCCCTCCCAGTGTGGAGAAGTGAAGAAC
ATCACGGGCCACACCCTGCTGGAGCGGTATGGCATGACCGAGATCGGCATGGCTCTGTCCGGGCCCTG
ACCACTGCCGTGCGCCTGCCAGTTCCTGTTGGGGACCCACTGCCTGGAGTACAGGTGCGCATTGTCTCA
GAAAACCCACAGAGGGAAGCCTGCTCTACACCATCCACGCAGAGGGAGACGAGAGGGGGACCAAGGTG
ACCCAGGGTTTGAAGAAAAGGAGGGGGAGCTGCTGGTGGAGGGGACCCCTCCGTGTTTCGAGAATACTGG
AATAAACCAAGAACTAAGAGTGCATTCACCCTGATGGTGGTAAAGACAGGGGACACCGTGGTG
TTAAGGATGGCCAGTACTGGATCCGAGCCGGACCTCAGTGGACATCATCAAGACTGGAGGCTACAAG
GTCAGCGCCCTGGAGGTGGAGTGGACCTGCTGGCCACCCAGCATCACAGATGGCTGTGATTGGA
GTTCCGGATATGACATGGGGCCAGCGGTCACTGCTGTGGTACCCTCCGAGAAGGACACTCACTGTCC
CACAGGGAGCTCAAAGAGTGGGCCAGAAATGCTCCTGGCCCGTACGCGGTGCCCTCGGAGCTGGTCTG
GTGGAGGAGATCCCGGGAACAGATGGCAAGATTGACAAGAAGGCGCTCATCAGGCACTTCCACCCC
TCATGA
  
```

Restriction Sites: SgfI-MluI
ACCN: NM_001243279



[View online >](#)

Insert Size:	1731 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001243279.2
RefSeq Size:	3751 bp
RefSeq ORF:	1731 bp
Locus ID:	197322
UniProt ID:	Q4G176
Cytogenetics:	16q24.3
MW:	64.1 kDa
Gene Summary:	<p>This gene encodes a member of the acyl-CoA synthetase family of enzymes that activate fatty acids by catalyzing the formation of a thioester linkage between fatty acids and coenzyme A. The encoded protein is localized to mitochondria, has high specificity for malonate and methylmalonate and possesses malonyl-CoA synthetase activity. Mutations in this gene are a cause of combined malonic and methylmalonic aciduria. Alternatively spliced transcript variants have been observed for this gene. [provided by RefSeq, Sep 2013]</p> <p>Transcript Variant: This variant (4) differs in the 5' UTR compared to variant 1. Variants 1, 2 and 4 encode the same isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>