

Product datasheet for **SC332004**

STARD13 (NM_001243466) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	STARD13 (NM_001243466) Human Untagged Clone
Tag:	Tag Free
Symbol:	STARD13
Synonyms:	ARHGAP37; DLC2; GT650; LINC00464
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

Fully Sequenced ORF: >SC332004 representing NM_001243466.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGCTGGAACCCCTCCTCAGTTCTCCATGCAAACGTTAACCAGGCCCTTTGTGGTGCTTGGTGCTGCGC
TGGTGCAGAGAATGCAAAGACACTGTCTGTGGTGGGAAACAGAAAAGCAGAGTGAACCACACATCCAG
CGCCGGGAAATTGAGGCAAAAAGAAGCATGTGACTGGCTCCGTGCTGCCGGGTTCCCGCAATACGCTCAG
TTATATGAGGATTCACAATTTCCCATCAACATTGTGGCTGTCAAGAATGATCATGATTTTCTTGAAAAG
GACCTTGTAGAACCTCTTTCGAGACGACTAAATACGTTGAACAAGTGTGCCTCAATGAAAACCTTGATGTG
AACTTCCAAGGAAAAAGGGTGACGACTCCGATGAGGAAGATCTTTGTATCAGCAACAAATGGACTTTC
CAAAGAACCAGTCGCAGGTGGTCTCGTGTGGACGACCTCTACACGCTGCTCCCTCGAGGAGACAGAAAT
GGGTACCAGGGAGGCACGGGATGAGGAACACGACCAGCAGTGTGAGAGCGTCTCACAGACCTGAGCGAG
CCTGAGGTCTGCTCCATTACAGCGAAAGCAGTGGAGGCAGCGACAGTGCAGCCAGCCGGGCCAGTGC
TGTACAGACAACCCGGTCTGCTGGATGCCCACTCGTACGAGCAGCCTCCACAGCCCCCAGAGAT
GTCCTCAACCACCCCTTCCACCCAAGAATGAGAAGCCACGAGGGCTAGGGCCAAATCATTTTTGAAA
CGCATGGAACACTCCGAGGGAAGGGAGCCACGGGAGGCATAAGGGGTCTGGGCGGACAGTGGCCTG
GTGATCAGTGGGCCATGTTGCAGCAGGAGCCAGAGTCCTTAAGGCTATGCAGTGCATCCAATACCA
AATGGAGATCTCCAGAATTCGCCGCCACCTGCCTGCAGAAAAGGGCTCCCATGCTCTGGCAAGTCGAGT
GGCAGAGCAGCCCGTCGGAGCACAGCAGCAGCGGGGTGAGCACGCCCTGCCTGAAGGAACGCAAGTGC
CACGAGGCCAACAAAGCGCGGGGGCATGTACTTGGAGGACCTAGATGTGCTGGCGGGGACAGCACTGCCG
GATGCAGGGGACCAAGCCGATGCATGAATTTCACTCCCAAGAGAATTTGGTGGTGCATATCCCAAG
GATCACAACCCAGGAACATTCACCAAGGCACTTCTATTGAAAGCCTCTCTCCACAGATAGTAGCAAT
GGGGTTAATTGGAGGACCGGTAGCATCTCCCTGGCAGAGAGCAGTCCCTGGTGCAGGGAGCCCGG
CTCATGGCGTCTGCCACAGAGCCAGCCGAGTATCTATGACAATGTCCTTGGCTCCCATCTGTAT
GCCAGCAGGAGATCTTTTGGACTTGGAGAAAGATGACCTTTTCCCTCACTTGGATGACATTCTGCAG
CATGTCAATGGGCTCCAAGAGGTAGTCGATGACTGGTCCAAGATGTCTTGCCTGAACTGAAAACCTCAT
GATACATTGGTTGGGAAACCTGGCTTATCCACCTTCCATCTCCTAATCAGATCACCTTAGATTTTGAA
GGTAACCTGTCTCAGAAGGTCGGACGACACCCAGTGTGTGAAAGAGATGTAACATCTCTAATGAA
TCTGAGCCTCTGGGTCAGAGACAGGAGGATTCTGGTGTAGGGGCTCTCTGACCAGGCCAAACAGG
CGACTCCGATGGAACAGTTTCCAGCTGTCGCACCAGCCCCGGCCGCCCCAGCATCGCCCCACATCAGC
AGCCAGACGGCCAGCCAGCTGAGCCTGCTCCAGCGTTTCTACTGCTCCGCTCACGCCATCATGGAG
AAGCACTCCATGTCCAACAAGCACGGCTGGACATGGTCAGTTCCAAAGTTCATGAAGAGGATGAAAGTT
CCCGACTACAAAGACAAGGCTGTCTTTGGCGTTCTCTCATAGTCCACGTCCAAAGAACGGGACAGCCC
CTGCCTCAAAGTATTCAGCAAGCACTGAGATATCTACGCAGCAACTGCCTCGATCAGGAGTGA
  
```

Restriction Sites: Sgfl-MluI

ACCN: NM_001243466

Insert Size: 2064 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001243466.1](#)

RefSeq Size: 2825 bp

RefSeq ORF: 2064 bp

Locus ID: 90627

UniProt ID: [Q9Y3M8](#)

Cytogenetics: 13q13.1-q13.2

MW: 76.4 kDa

Gene Summary: This gene encodes a protein which contains an N-terminal sterile alpha motif (SAM) for protein-protein interactions, followed by an ATP/GTP-binding motif, a GTPase-activating protein (GAP) domain, and a C-terminal STAR-related lipid transfer (START) domain. It may be involved in regulation of cytoskeletal reorganization, cell proliferation, and cell motility, and acts as a tumor suppressor in hepatoma cells. The gene is located in a region of chromosome 13 that is associated with loss of heterozygosity in hepatocellular carcinomas. Alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Aug 2011]

Transcript Variant: This variant (5) contains an alternate 5' terminal exon, and differs at the 3' end compared to variant 1. This results in a shorter isoform (5) with distinct N- and C- termini compared to isoform 1.