

Product datasheet for **SC332118**

CTAGE5 (MIA2) (NM_001247989) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CTAGE5 (MIA2) (NM_001247989) Human Untagged Clone
Tag:	Tag Free
Symbol:	CTAGE5
Synonyms:	CTAGE5; MEA6; MGEA; MGEA6; MGEA11; TALI
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

Fully Sequenced ORF: >SC332118 representing NM_001247989.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGAGGAGCCCGGGTTACCCCTCAACCGTATTTGGGGCTGCTCCTGGAGGAGCTACGCAGGGTTGTG
GCAGCACTGCCTGAAGGTATGAGACCAGATTCTAATCTTTATGGTTTTCCATGGGAATTGGTGATATGT
GCAGCTGTTGTTGGATTTTTGCTGTTCTCTTTTTTTTGTGGAGAAGTTTAGATCGGTTAGGAGTCGG
CTTTATGTGGGACGAGAGAAAAAGCTTGCTCTAATGCTTTCTGGACTAATTGAAGAAAAAGTAAACTA
CTTGAAAAATTTAGCCTTGTTCAAAAAGAGTATGAAGGCTATGAAGTAGAGTCATCTTTAAAGGATGCC
AGCTTTGAGAAGGAGGCAACAGAAGCACAAAGTTTGGAGGTAGAAAAATCAATGGCAACCTGTGAAAAG
CTGAACAGGTCCAATTCTGAACCTGAGGATGAAATACTCTGTCTAGAAAAAGAGTTAAAAGAAGAGAAA
TCCAAACATTCTGAACAAGATGAATTGATGGCGGATATTTCAAAAAGGATACAGTCTCTAGAAGATGAG
TCAAAAATCCCTCAATCACAAGTAGCTGAAGCCAAAATGACCTTCAAGATATTTCAATGAATGAAGAA
CGACTGAAGATAGCAATAAAAGATGCTTTGAATGAAAATCTCAACTTCAGGAAAGCCAGAAACAGCTT
TTGCAAGAAGCTGAAGTATGGAAGAACAAGTGAAGTAAATAAACAGAAAGTAACATTTGAAGAC
TCCAAAGTACATGCAGAACAAGTTCTAATGATAAAGAAAAGTACATCAAGACTCTGACTGAACGCTTG
TTAAAGATGAAAGATTGGGCTGCTATGCTTGGAGAAGACATAACGGATGATGATAACTTGAATAGAA
ATGAACAGTGAATCGGAAAATGGTGCTTACTTAGATAATCCTCCAAAAGGAGCTTTGAAGAACTGATT
CATGCTGCTAAGTTAAATGCTTCTTTAAAAACCTTAGAAGGAGAAAAGAAACCAATTTATATTCAGTTG
TCTGAAGTTGATAAAAACAAGGAAGAGCTTACAGAGCATATTAATAATCTTCAGACTGAACAAGCATCT
TTGCAGTCAGAAAAACACACATTTTGAATGAGAATCAGAAGCTTCAACAGAACTTAAAGTAATGACT
GAATTAATCAAGAAAATGAAATGAAACTCCACAGGAAATTAACAGTAGAGGAAAATATCGGTTAGAG
AAAGAAGAGAACTTTCTAAAGTAGATGAAAAGATCAGCCATGCCACTGAAGAGCTGGAGACCTATAGA
AAGCGAGCCAAAGATCTTGAAGAAGAATTGGAGAGAATTCATTCTTATCAAGGCAGATTATTTCC
CATGAGAAAAAAGCACATGATAATTGGTTGGCAGCTCGGAATGCTGAAAAGAAACCTCAATGATTTAAGG
AAAGAAAATGCTCACAACAGACAAAAATTAACGAAACAGAGCTTAAATTTGAACTTTTAGAAAAAGAT
CCTTATGCACTCGATGTTCCAAATACAGCATTTGGCAGAGAGCATTCCCCATATGGTCCCTCACCATTG
GGTTGGCCTTCTGAAACAAGAGCTTTTCTCTCCTCCAACCTTTGTTGGAGGGTCCACTCAGACTC
TCACCTTTGCTTCCAGGGGAGGAGGAAGAGGCTCAGAGGCCAGGGAATCCTCTGGACCATCAGATT
ACCAATGAAAGAGGAGAATCAAGCTGTGATAGGTTAACCGATCCTCATAGGGCTCCCTCTGACACTGGG
TCTCTGTACCTCCATGGGACCAGGACCGTAGGATGATGTTTCTCCGCCAGGACAATCATATCCTGAT
TCAGCCCTTCTCCAAAAGGCAAGACAGATTTTGTCTAATTCTGGTAGACTGTCTGGACCAGCAGAA
CTCAGAAGTTTTAATATGCCTTCTTTGGATAAAAATGGATGGGTCAATGCCTTCAAGAAATGGAATCCAGT
AGAAATGATACCAAAGATGATCTTGTAATTTAAATGTGCCTGATTCATCTCCTGCTGAAAAATGAA
GCCACTGGCCCTGGCTTTGTTCCCTCCACCTCTTGTCTCAATCAGAGGTCCATTGTTTCCAGTGGATGCA
AGAGGCCCATCTTGAGAAGAGGACCTCTTTCCCCCACCTCCTCCAGGAGCCATGTTTGGAGCTTCT
CGAGATTATTTTCCACCAGGGGATTTCCAGGTCCACCACCTGCTCCATTTGCAATGAGAATGTCTAT
CCACCGAGGGGTTTTCTCCTTACCTTCCCCAAGACCTGGATTTTTCCCCCACCCACATTCTGAA
GGTAGAAGTGAAGTCCCCTCAGGTTTGATTCCACCTTCAATGAGCCTGCTACTGAACATCCAGAACCA
CAGCAAGAAACCTGA
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001247989

Insert Size: 2430 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001247989.1
RefSeq Size:	3695 bp
RefSeq ORF:	2430 bp
Locus ID:	4253
UniProt ID:	Q96PC5
Cytogenetics:	14q21.1
Protein Families:	Transmembrane
MW:	91.6 kDa
Gene Summary:	<p>This gene encodes s receptor in the endoplasmic reticulum, which plays a role in the export of large pre-chylomicrons and pre-very low density lipoproteins (pre-VLDLs). Three major classes of transcripts are generated from this gene- melanoma inhibitory activity 2-specific transcripts, cTAGE family member 5-specific transcripts and transcripts that include exons from both these transcript species (TANGO1-like or TALI). Additionally, alternative splicing in these transcripts results in multiple transcript variants encoding multiple isoforms. [provided by RefSeq, Sep 2016]</p> <p>Transcript Variant: This variant (6) represents the longest transcript and encodes the longest isoform (6, also known as TALI). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>