

Product datasheet for **SC332152**

SPOCK3 (NM_001251967) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: SPOCK3 (NM_001251967) Human Untagged Clone
Tag: Tag Free
Symbol: SPOCK3
Synonyms: HSAJ1454; TES-3; TICN3
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC332152 representing NM_001251967.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGCTCAAGGTGTCAGCCGACTGTGTGTGTGTGCAGCCGCTTGGTGCAGTCAGTCTCTCGCAGCTGCC
GCGGCGGTGGCTGCAGCCGGGGGCGGTTCGGACGGCGGTAATTTTCTGGATGATAAAACAATGGCTCACC
ACAATCTCTCAGTATGACAAGGAAGTCGGACAGTGAACAATCCGAGACTGCAAAGTGAATATCAG
GCATGTGTCTTAGGAAAACAGATCTCAGTCAAATGTGAAGGACATTGCCCATGTCCTTCAGATAAGCCC
ACCAAGTACAAGCAGAAAATGTTAAGAGAGCATGCAGTGACCTGGAGTTCAGGGAAGTGCCAAACAGATTG
CGGGACTGGTTCAAGGCCCTTCATGAAAGTGAAGTCAAAAACAAGAAGACAAAAACATTGCTGAGGCCT
GAGAGAAGCAGATTCGATACCAGCATCTTGCCAAATTTGCAAGGACTCACTTGGCTGGATGTTAACAAGA
CTTGATACAAACTATGACCTGCTATTGGACCAGTCAGAGCTCAGAAGCATTTACCTTGATAAGAATGAA
CAGTGTACCAAGGCATTCTCAATTCTTGTGACACATACAAGGACAGTTAATATCTAATAATGAGTGG
TGCTACTGCTTCCAGAGACAGCAAGACCCACCTTGCCAGACTGAGCTCAGCAATATTCAGAAGCGGCAA
GGGGTAAAGAAGCTCCTAGGACAGTATATCCCCCTGTGTGATGAAGATGGTTACTACAAGCCAACACAA
TGTCATGGCAGTGTGGACAGTGTGGTGTGTTGACAGATATGGAATGAAGTATGGGATCCAGAATA
AATGGTGTGCAGATTGTGCTATAGATTTGAGATCTCCGGAGATTTGCTAGTGGCGATTTTCATGAA
TGGACTGATGATGAGGATGATGAAGACGATATTATGAATGATGAAGATGAAATTTGAAGATGATGAA
GATGAAGGGGATGATGATGATGGTGGTATGACCATGATGTATACATTGA
```

Restriction Sites: Sgfl-Mlul
ACCN: NM_001251967
Insert Size: 1017 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001251967.1](#)

RefSeq Size: 2695 bp

RefSeq ORF: 1017 bp

Locus ID: 50859

UniProt ID: [Q9BQ16](#)

Cytogenetics: 4q32.3

Protein Families: Secreted Protein, Transmembrane

MW: 38.3 kDa

Gene Summary: This gene encodes a member of a novel family of calcium-binding proteoglycan proteins that contain thyroglobulin type-1 and Kazal-like domains. The encoded protein and may play a role in adult T-cell leukemia by inhibiting the activity of membrane-type matrix metalloproteinases. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Oct 2011]

Transcript Variant: This variant (11) lacks four consecutive coding exons but maintains the reading frame, compared to variant 2. The encoded isoform (11) is shorter than isoform 2.