

## Product datasheet for SC332195

### TNIP1 (NM\_001252393) Human Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** TNIP1 (NM\_001252393) Human Untagged Clone  
**Tag:** Tag Free  
**Symbol:** TNIP1  
**Synonyms:** ABIN-1; NAF1; nip40-1; VAN  
**Vector:** pCMV6-Entry (PS100001)  
**Fully Sequenced ORF:** >SC332195 representing NM\_001252393.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGAAGGGAGAGGACCGTACCGGATCTACGACCTGGGGGACGCGTGCCTCAGGAGAGGCATCCGCA
GCTTTTGAGCGCTAGTGAAGGAGAATCCCGGCTGAAGGAAAAATGCAAGGGATAAAGATGTTAGGG
GAGCTTTTGAAGAGTCCAGATGGAAGCGACAGGCTCCGGCAGAAGGCAGAGGAGCTAGTGAAGGAC
AACGAGCTGCTCCACCACCTTCTCCCTCCTTGGGCTCCTTCGACCCCTGGCTGAGCTCACAGGAAAG
GACTCAAATGTCACAGCATCTCCACAGCCCTGCATGCCCCAGTGACAAGCCAGCACCAGTCCAGAAG
CCTCCATCCAGTGGCACCTCCTCTGAATTTGAAGTGGTCACTCCTGAGGAGCAGAATTCACCAGAGAGC
AGCAGCCATGCCAATGCGATGGCGCTGGGCCCTGCCCCGTGAGGACGGCAACCTGATGCTGCACCTG
CAGCGCCTGGAGACCACGCTGAGTGTGTGTGCCGAGGAGCCGGACCACGGCCAGCTTTCACCCACCTG
GGCCGCATGGCCCTGGAGTTCAACCGACTGGCATCCAAGGTGCACAAGAATGAGCAGCGCACCTCCATT
CTGCAGACCCTGTGTGAGCAGCTTCGGAAGGAGAACGAGGCTCTGAAGGCCAAGTTGGATAAGGGCCTG
GAACAGCGGGATCAGGCTGCCGAGAGGCTGCCGGAGGAAAAATTTGGAGCTCAAGAAGTTGTTGATGAGC
AATGGCAACAAAGAGGGTGCCTGCGCGGCCAGGCTCACCGAAGATGGAAGGGACAGGCAAGAAGGCA
GTGGCTGGACAGCAGCAGGCTAGTGTGACGGCAGGTAAGTCCCAGAGGTGGTGGCTTGGGCGCAGCC
GAGAAGAAGGTGAAGATGCTGGAGCAGCAGCGCAGTGAAGTGAACAAGCAGTGGGACCAG
CATTTCCGGTCCATGAAGCAGCAGTATGAGCAGAAGATCACTGAGCTGCGTCAGAAGCTGGCTGATTTG
CAGAAGCAGGTGACTGACCTGGAGGCCGAGCGGGAGCAGAAGCAGCGTGACTTTGACCGAAGCTCCTC
CTGGCCAAGTCCAAGATTGAAATGGAGGAGACCGACAAGGAGCAGCTGACAGCAGAGGCCAAGGAGCTG
CGCCAAAAGGTCAAGTACCTGCAGGATCAGCTGAGCCCACTACCCGACAGCGTGAGTACCAGGAAAAG
GAGATCCAGCGGCTCAACAAGGCCCTGGAGGAAGCACTGAGCATCAAACCCCGCCATCATCTCCACCA
ACAGCATTTGGGAGCCCAGAAGGAGCAGGGCCCTCCTAAGGAAACAGGAGCTGGTCACGCAGAATGAG
TTGCTGAAACAGCAGGTGAAGATCTTCGAGGAGGACTTCCAGAGGGAGCGCAGTATCGTGAAGCAGT
AATGAGGAGAAGGAAGAGCTGAAGAAGCAAGTGGAGAAGCTGCAGGCCAGGTACCCTGTCAAATGCC
CAGCTAAAAGCATTCAAAGATGAGGAGAAGGCAAGAGAAGCCCTCAGACAGCAGAAGAGGAAAGCAAAG
GCCTCAGGAGAGCGTTACCATGTGGAGCCCCACCAGAACATCTCTGCGGGCCCTACCCCTACGCTAC
CCGCCATGCCAGCCATGGTGCCACACCATGGCTTCGAGGACTGGTCCCAGATCCGCTACCCCTCCC
CCCATGGCCATGGAGCACCCGCCCACTCCCAACTCGCGCCTTCCATCTGCCGGAATACACCTGG
CGTCTACCCTGTGGAGGGTTTCAAATCCAAATCAGAGCTCCCAAGTGATGGACCCTCCACAGCCAGG
CCTACAGAACCAGAGCCAGCTGATCTCAGATTGCCAAGAACTAG
  
```



[View online »](#)

|                               |  |
|-------------------------------|--|
| <b>Restriction Sites:</b>     | Sgfl-Mlul  |
| <b>ACCN:</b>                  | NM_001252393   |
| <b>Insert Size:</b>           | 1908 bp  |
| <b>OTI Disclaimer:</b>        | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).   |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).   |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>  |
| <b>RefSeq:</b>                | <a href="#">NM_001252393.1</a>   |
| <b>RefSeq Size:</b>           | 2859 bp  |
| <b>RefSeq ORF:</b>            | 1908 bp  |
| <b>Locus ID:</b>              | 10318  |
| <b>UniProt ID:</b>            | <a href="#">Q15025</a>   |
| <b>Cytogenetics:</b>          | 5q33.1   |
| <b>MW:</b>                    | 71.8 kDa   |
| <b>Gene Summary:</b>          | <p>This gene encodes an A20-binding protein which plays a role in autoimmunity and tissue homeostasis through the regulation of nuclear factor kappa-B activation. Mutations in this gene have been associated with psoriatic arthritis, rheumatoid arthritis, and systemic lupus erythematosus. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2011]</p> <p>Transcript Variant: This variant (7) uses two alternate splice sites in the 3' coding region, compared to variant 1. The resulting isoform (4) is longer and contains a distinct C-terminus, compared to isoform 1. Variants 6 and 7 encode the same protein (isoform 4).</p> |