

Product datasheet for **SC332210**

ACSS1 (NM_001252675) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ACSS1 (NM_001252675) Human Untagged Clone
Tag:	Tag Free
Symbol:	ACSS1
Synonyms:	ACAS2L; ACECS1; AceCS2L
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

Fully Sequenced ORF: >SC332210 representing NM_001252675.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGCGGCGCGCACCCCTGGGCCGCGGCGTGGGAGGCTGCTGGGCAGCCTGCGAGGGCTCTCGGGGCGAG
CCC GCGCGGCCGCGTGC GGGGTGAGCGCGCCGCGCAGGGCGGCTCGGGACCCCTCGGGCAGCGCTCCC
GCAGTTGCAGCAGCAGCAGCACAGCCAGGCTCGTATCCCGCGCTGAGTGCACAGGCAGCCCGGGAGCCG
GCCGCTTCTGGGGCCCTCTGGCGCGGGACACTCTCGTGTGGGACACCCCTACCACACCGTCTGGGAC
TGGCACTTCAGCACTGGCAAGATCGGCTGGTTCCTGGGAGGCCAGTTAAATGTCTCTGTCAACTGCTTG
GACCAGCATGTTTCGGAAGTCCCCGAGAGCGTTGCTTTGATCTGGGAGCGCGATGAGCCTGGAACGGAA
GTGAGGATCACCTACAGGGAACACTGGAGACCACGTGCCGCTGGCCAACACGCTGAAGAGGCATGGA
GTCCACCGTGGGACCGTGTGCCATCTACATGCCGCTGCCCATTTGGCTGTGGCAGCAATGCTGGCC
TGTGCCAGGATCGGAGCTGTCCACACAGTCATCTTTGCTGGCTTCAGTGCAGAGTCTTGGCTGGGAGG
ATCAATGATGCCAAGTGAAGGTGGTTATCACCTCAACCAAGGACTCCGGGGTGGGCGCGTGGTGGAG
CTGAAGAAAATAGTGGATGAGGCTGTGAAGCACTGCCCCACCGTGCAGCATGTCTGGTGGCTCACAGG
ACAGACAACAAGTCCACATGGGGGATCTGGACGTCCCGTGGAGCAGGAAATGGCCAAGGAGGACCCCT
GTTTGGCGCCCAGAGAGCATGGGCAAGTGGGACATGCTCTTCATGCTGTACACCTCAGGGAGCACCGGA
ATGCCCAAGGGCATCGTCCATACCCAGGAGGCTACCTGCTCTATGCCGCGCTGACTCACAAGCTTGTG
TTTGACCACCAGCCAGGTGACATCTTTGGCTGTGTGGCCGACATCGGTTGGATTACAGGACACAGCTAC
GTGGTGTATGGGCCTCTCTGCAATGGTGCACACGCGTCTTTTTGAGAGCACCCAGTTTATCCCAAT
GCTGGTGGTACTGGGAGACAGTAGAGAGGTTGAAGATCAATCAGTTCTATGGCGCCCCAACGGCTGTC
CGGCTGTTGCTGAAATACGGTGTGCCTGGGTGAAGAAGTATGATCGCTCCTCCCTGCGGACCCCTGGGG
TCAGTGGGAGAGCCATCAACTGTGAGGCTGGGAGTGGCTTACAGGGTGGTGGGGACAGCAGGTGC
ACGCTGGTGGACACCTGGTGGCAGACAGGTGGCATCTGCATCGCACACGCGCCCTCGGAAGAAGGGCGG
GAAATCCTCCCTGCCATGGCGATGAGGCCCTTCTTTGGCATCGTCCCGCTCCTCATGGATGAGAAGGGC
AGCGTCTGTGGAGGCGAGCAACGTCTCCGGGGCCCTGTGCATCTCCAGGCGCTGGCCGGGCATGGCCAGG
ACCATCTATGGCGACCACCAGCGATTTGTGGACGCTACTTCAAGGCCTACCCAGGCTATTACTTCACT
GGAGACGGGGCTTACCGAAGTGGGGCGGCTATTACCAGATCACAGGGCGGATGGATGATGTCATCAAC
ATCAGTGGCCACCGGCTGGGGACCCGAGAGATTGAGGACGCCATCGCCGACCACCCTGCAGTACCAGAA
AGTGTGTCATTGGCTACCCACGACATCAAAGGAGAAGCTGCCTTTGCCTTATTGTGGTGAAGAT
AGTGGGGTACTCAGATGTGGTGGTGCAGGAGCTCAAGTCCATGGTGGCCACCAAGATCGCCAAATAT
GCTGTGCCTGATGAGATCCTGGTGGTGAACGTCTTCAAAAACCAGGTCTGGGAAGGTCATGCGGCGG
CTCCTGAGGAAGATCATCACTAGTGGGCCAGGAGCTGGGAGACTACCACCTTGGAGGACCCACG
ATCATCGCAGAGATCCTGAGTGTCTACCAGAAGTGAAGGACAAGCAGGCTGCTGCTAAGTGA
  
```

Restriction Sites: Sgfl-MluI

ACCN: NM_001252675

Insert Size: 2064 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001252675.1
RefSeq Size:	3686 bp
RefSeq ORF:	2064 bp
Locus ID:	84532
UniProt ID:	Q9NUB1
Cytogenetics:	20p11.21
Protein Pathways:	Glycolysis / Gluconeogenesis, Metabolic pathways, Propanoate metabolism, Pyruvate metabolism
MW:	74.6 kDa
Gene Summary:	<p>This gene encodes a mitochondrial acetyl-CoA synthetase enzyme. A similar protein in mice plays an important role in the tricarboxylic acid cycle by catalyzing the conversion of acetate to acetyl CoA. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Nov 2011]</p> <p>Transcript Variant: This variant (2) uses an alternate splice site in the coding region, but maintains the reading frame, compared to variant 1. The encoded isoform (2) is shorter than isoform 1.</p>