

Product datasheet for **SC332222**

LGR7 (RXFP1) (NM_001253727) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LGR7 (RXFP1) (NM_001253727) Human Untagged Clone
Tag:	Tag Free
Symbol:	LGR7
Synonyms:	LGR7; RXFPR1
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

Fully Sequenced ORF: >SC332222 representing NM_001253727.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGACATCTGGTCTGTCTTCTTCTACATCTTAATTTTTGGAAAATATTTTTCTCATGGGGTGGACAG
GATGTCAAAGTGTCCCTGGCTATTTCCCTGTGGGAACATCACAAAGTGCTTGCCTCAGCTCCTGCAC
TGTAACGGTGTGGACGACTGCGGAATCAGGCCGATGAGGACAACGTGGAGACAACAATGGATGGTCT
CTGCAATTTGACAAATATTTTGCCAGTTACTACAAAATGACTTCCCAATATCCTTTTGAGGCAGAAACA
CCTGAATGTTCTTTTCATTTCTTTCTTTTATTACCCTCCTCTTCTGGTCCCTCACTGTATCATGCC
TTACCACTTCCCCTTGATTGAGTGGTGGTCTGTGCCAGTGCAATGTCTTTGCCAAGGTCTGGAGCTT
GACTGTGATGAAACCAATTTACGAGCTGTTCCATCGGTTTCTTCAAATGTGACTGCAATGCACTTCAG
TGGAACTTAATAAGAAAGCTTCCCTGATTGCTTCAAGAATTATCATGATCTTCAGAAGCTGTACCTG
CAAAACAATAAGATTACATCCATCTCCATCTATGCTTTCAGAGGACTGAATAGCCTTACTAACTGTAT
CTCAGTCATAACAGAATAACCTTCTGAAGCCGGGTGTTTTGAAGATCTTACAGACTAGAATGGCTG
ATAATTGAAGATAATCACCTCAGTCGAATTTCCCACCAACATTTTATGGACTAAATTTCTTATTCTC
TTAGTCTGATGAATAACGTCTCACCCGTTTACCTGATAAACCTCTCTGTCAACACATGCCAAGACTA
CATTGGCTGGACCTTGAAGGCAACCATATCCATAATTTAAGAAATTTGACTTTTATTTCCCTGCAGTAA
TAACTGTTTTAGTGATGAGGAAAAACAAAATTAATCACTTAAATGAAAATACTTTTGCACCTCTCCAG
AACTGGATGAATTGGATTTAGGAAGTAATAAGATTGAAAATCTTCCACCGCTTATATTCAGGACCTG
AAGGAGCTGTCACAATTGAATCTTCTATAATCCAATCCAGAAAATTCAGCAACCAATTTGATTAT
CTTGTCAAACTCAAGTCTCTCAGCCTAGAAGGGATTGAAATTTCAAATATCCAACAAAGGATGTTTGA
CCTCTTATGAATCTCTCACATATATTTAAGAAATCCAGTACTGTGGGTATGCACCACATGTTGCG
AGCTGTAACCAAACTGATGGAATTTTCACTCTAGAGAATCTTGGCAAGCATTATTCAGAGAGTA
TTTGTCTGGGTTGTATCTGCAGTTACCTGCTTTGGAAACATTTTGTCAATTTGCATGCGACCTTATATC
AGGTCTGAGAACAAGCTGTATGCCATGTCAATCATTTTCTCTGCTGTGCCGACTGCTTAATGGGAATA
TATTTATTCGTGATCGGAGGCTTTGACCTAAAGTTTTCGTGGAGAATAACAATAAGCATGCGCAGCTGTGG
ATGGAGAGTACTCATTGTCAGCTTGTAGGATCTTTGGCCATTCTGTCCACAGAAGTATCAGTTTTACTG
TTAACATTTCTGACATTGGAAAAATACATCTGCATTGTCTATCCTTTTAGATGTGTGAGACCTGGAAAA
TGCAGAACAATTACAGTTCTGATTCTCATTGGATTACTGGTTTTATAGTGGCTTTCATTCCATTGAGC
AATAAGGAATTTTCAAAACTACTATGGCACAATGGAGTATGCTTCCCTTTCATTGAGAAGATACA
GAAAGTATTGGAGCCAGATTTATTCAGTGGCAATTTTCTTGGTATTAATTTGGCCGATTTATCATC
ATAGTTTTTCTATGGAAGCATGTTTTATAGTGTTCATCAAAGTGCCATAACAGCAACTGAAATACGG
AATCAAGTAAAAAAGAGATGATCCTTGCCAAAGTTTTTTCTTTATAGTATTTACTGATGATTATGC
TGGATACCCATTTTGTAGTGAATTTCTTCACTGCTTCAAGTGAAGTACCAGGTACCATAACCTCT
TGGGTAGTGATTTTATTTCTGCCATTAACAGTCTTTGAACCAATTTCTTATACTCTGACCACAAGA
CCATTTAAGAAATGATTCATCGGTTTTGGTATAACTACAGACAAAGAAAATCTATGGACAGCAAAGGT
CAGAAAACATATGCTCCATCATTCTGCTGGTGGAAATGTGGCCACTGCAGGAGATGCCACCTGAGTTA
ATGAAGCCGACCTTTTACATACCCCTGTGAAATGCTACTGATTTCTCAATCAACGAGACTCAATTC
TATTCATGA
  
```

Restriction Sites: Sgfl-Mlul

ACCN: NM_001253727

Insert Size: 2355 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001253727.1
RefSeq Size:	3955 bp
RefSeq ORF:	2355 bp
Locus ID:	59350
Cytogenetics:	4q32.1
Protein Families:	Druggable Genome, GPCR, Transmembrane
Protein Pathways:	Neuroactive ligand-receptor interaction
MW:	90.1 kDa
Gene Summary:	<p>This gene encodes a member of the leucine-rich repeat-containing subgroup of the G protein-coupled 7-transmembrane receptor superfamily. The encoded protein plays a critical role in sperm motility, pregnancy and parturition as a receptor for the protein hormone relaxin. Decreased expression of this gene may play a role in endometriosis. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Dec 2011]</p> <p>Transcript Variant: This variant (2) uses an alternate splice site in the 5' coding region, but maintains the reading frame, compared to variant 1. The encoded isoform (2) is longer than isoform 1.</p>