

Product datasheet for **SC332234**

GPR172A (SLC52A2) (NM_001253815) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: GPR172A (SLC52A2) (NM_001253815) Human Untagged Clone
Tag: Tag Free
Symbol: GPR172A
Synonyms: BVVLS2; D15ErtD747e; GPCR41; GPR172A; hRFT3; PAR1; RFT3; RFVT2
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC332234 representing NM_001253815.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGCAGCACCCACGCCCGCCCGTCCGGTCTGACCCACCTGCTGGTGGCTCTCTTCGGCATGGGCTCC
TGGGCTGCGGTCAATGGGATCTGGGTGGAGCTACCTGTGGTGGTCAAAGAGCTTCCAGAGGGTTGGAGC
CTCCCCTTTACGTCTCTGTGCTGTGGCTCTGGGGAACCTGGGTCTGCTGGTGGTGACCCTCTGGAGG
AGGCTGGCCCCAGGAAAGGACGAGCAGGTCATCCGGTGGTGCAGGTGCTGGGCATGGTGGGCACA
GCCCTGTGGCCTCTCTGTGGCACCATGTGGCCCACTGGCAGGACAGTGCATTCTGTGGCCTTCTTA
GCACTGGCCTTTGTGCTGGCACTGGCATGCTGTGCCTCGAATGCACTTTCTGCCCTTCTGTAGCCAC
CTGCCACCTCGCTTCTTACGGTCACTTCTCCTGGGTCAAGGCCTGAGTGCCTGCTGCCCTGCGTGTG
GCCCTAGTGCAGGGTGTGGGCGCCTCGAGTGCAGCCGACCCCATCAACGGCACCCCTGGCCCCCG
CTCGACTTCTTGAGCGTTTTCCCGCCAGCACCTTCTTCTGGCACTGACTGCCCTTCTGGTTCGCTTCA
GCTGCTGCCTTCCAGGTCTTCTGCTGTGTTGCCGCCACCACCTGTACCCACAGGGGAGTTAGGA
TCAGGCCTCCAGGTGGGAGCCCCAGGAGCAGAGGAAGAGGTGGAAGAGTCCTCACCCTGCAAGAGCCA
CCAAGCCAGGCAGCAGGCACCACCCTGGTCCAGACCTAAGGCCTATCAGCTTCTATCAGCCCGCAGT
GCCTGCCTGCTGGGCCTGTTGGCCGCCACCAACGCGCTGACCAATGGCGTGTGCTGCCGTGCAGAGC
TTTTCTGCTTACCCTACGGGCGTCTGGCCTACCACCTGGTGTGGTGTGGCAGTGTGCAATCCC
CTGGCCTGCTTCCATGGGTGTGCTGTGCAGGTCTTGGCAGGGCTGGGCGCCTCTCTCTGCTG
GGCGTGTCTGTGGGGCTACCTGATGGCGCTGGCAGTCTGAGCCCTGCCCGCCCTGGTGGGCACC
TCGGCGGGGGTGGTCTCGTGGTGTGCTGTGCTGGTGTGCTTGGCGTGTCTCTACGTGAAGGTG
GCAGCCAGCTCCCTGCTGCATGGCGGGGCGCCGGCATTGCTGGCAGCCGGCGTGGCCATCCAGGTG
GGCTCTGCTCGGCGTGTGCTATGTTCCCCCGACCAGCATCTATCACGTGTCCACAGCAGAAAG
GACTGTGCAGACCCTGTGACTCTGA
  
```

Restriction Sites: SgfI-MluI
ACCN: NM_001253815
Insert Size: 1338 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001253815.1
RefSeq Size:	1806 bp
RefSeq ORF:	1338 bp
Locus ID:	79581
UniProt ID:	Q9HAB3
Cytogenetics:	8q24.3
Protein Families:	Druggable Genome, GPCR, Transmembrane
MW:	45.8 kDa
Gene Summary:	<p>This gene encodes a membrane protein which belongs to the riboflavin transporter family. In humans, riboflavin must be obtained by intestinal absorption because it cannot be synthesized by the body. The water-soluble vitamin riboflavin is processed to the coenzymes flavin mononucleotide (FMN) and flavin adenine dinucleotide (FAD) which then act as intermediaries in many cellular metabolic reactions. Paralogous members of the riboflavin transporter gene family are located on chromosomes 17 and 20. Unlike other members of this family, this gene has higher expression in brain tissue than small intestine. Alternative splicing of this gene results in multiple transcript variants encoding the same protein. Mutations in this gene have been associated with Brown-Vialetto-Van Laere syndrome 2 - an autosomal recessive progressive neurologic disorder characterized by deafness, bulbar dysfunction, and axial and limb hypotonia. [provided by RefSeq, Jul 2012]</p> <p>Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Variants 1, 2, 3, 5, 6 and 7 encode the same isoform (1).</p>