

## Product datasheet for **SC332277**

### COLEC11 (NM\_001255985) Human Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** COLEC11 (NM\_001255985) Human Untagged Clone  
**Tag:** Tag Free  
**Symbol:** COLEC11  
**Synonyms:** 3MC2; CL-11; CL-K1-I; CL-K1-II; CL-K1-IIa; CL-K1-IIb; CLK1  
**Vector:** pCMV6-Entry (PS100001)  
**Fully Sequenced ORF:** >SC332277 representing NM\_001255985.  
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGAAGAAACAAAGAGGAGTTGGTGCCTGCCTGCGCTCAGGATGAGGGGAATCTGGCCCTGGTGGGC
GTTCTAATCAGCCTGGCCTTCCTGTCACTGCTGCCATCTGGACATCCTCAGCCGGCTGGCGATGACGCC
TGCTCTGTGCAGATCCTCGTCCCTGGCCTCAAAGGGGATGCGGGAGAGAAGGGAGACAAAGGCGCCCCC
GGACGGCCTGGAAGAGTCGGCCCCACGGGAGAAAAAGGAGACATGGGGGACAAAGGACAGAAAGGCAGT
GTGGGTGCTCATGAAAAATTGGTCCCATTGGCTCTAAAGGTGAGAAAGGAGATTCCGGTGACATAGGA
CCCCCTGGTCCTAATGGAGAACCAGGCCTCCCATGTGAGTGCAGCCAGCTGCGCAAGGCCATCGGGGAG
ATGGACAACCAGGTCTCTCAGCTGACCAGCGAGCTCAAGTTCATCAAGAATGCTGTGCGCCGTGTGCGC
GAGACGGAGAGCAAGATCTACCTGCTGGTGAAGGAGGAGAAGCGCTACGCGGACGCCAGCTGTCCTGC
CAGGGCCGCGGGGGCACGCTGAGCATGCCAAGGACGAGGCTGCCAATGGCCTGATGGCCGCATACCTG
GCGCAAGCCGGCCTGGCCCGTGTCTTTCATCGGCATCAACGACCTGGAGAAGGAGGGCGCCTTCGTGTAC
TCTGACCACTCCCCATGCGGACCTTCAACAAGTGGCGCAGCGGTGAGCCCAACAATGCCTACGACGAG
GAGGACTGCGTGGAGATGGTGGCCTCGGGCGGCTGGAACGACGTGGCCTGCCACACCACCATGTACTTC
ATGTGTGAGTTTGACAAGGAGAACATGTGA
```

**Restriction Sites:** Sgfl-Mlul  
**ACCN:** NM\_001255985  
**Insert Size:** 858 bp  
**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).  
**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001255985.1](#)

**RefSeq Size:** 1595 bp

**RefSeq ORF:** 858 bp

**Locus ID:** 78989

**UniProt ID:** [Q9BWP8](#)

**Cytogenetics:** 2p25.3

**MW:** 30.2 kDa

**Gene Summary:** This gene encodes a member of the collectin family of C-type lectins that possess collagen-like sequences and carbohydrate recognition domains. Collectins are secreted proteins that play important roles in the innate immune system by binding to carbohydrate antigens on microorganisms, facilitating their recognition and removal. The encoded protein binds to multiple sugars with a preference for fucose and mannose. Mutations in this gene are a cause of 3MC syndrome-2. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Dec 2011]  
Transcript Variant: This variant (6) differs in the 5' UTR and initiates translation at an alternate start codon, compared to variant 1. The encoded isoform (f) is longer and has a distinct N-terminus, compared to isoform a.