

## Product datasheet for **SC332281**

### PSMC3IP (NM\_001256015) Human Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** PSMC3IP (NM\_001256015) Human Untagged Clone  
**Tag:** Tag Free  
**Symbol:** PSMC3IP  
**Synonyms:** GT198; HOP2; HUMGT198A; ODG3; TBPIP  
**Vector:** pCMV6-Entry (PS100001)  
**Fully Sequenced ORF:** >SC332281 representing NM\_001256015.  
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGGTGAGTGATGCTGACCTTCAAGTCCTAGATGGCAAAATCGTGGCCCTCACTGCTAAGGTGCAGAGC
TTGCAGCAGAGCTGCCGCTACATGGAGGCTGAGCTCAAGGAATTCTAGTGCCCTGACCACACCAGAG
ATGCAGAAAGAAATCCAGGAGTTAAAGAAGGAATGCGCTGGCTACAGAGAGAGATTGAAGAACATTTAAA
GCAGCTACCAATCATGTGACTCCAGAAGAGAAAGAGCAGGTGTACAGAGAGAGGCAGAAGTACTGTAAG
GAGTGGAGGAAGAGGAAGAGGATGGCTACAGAGCTGTCTGATGCAATACTTGAAGGATACCCCAAGAGC
AAGAAGCAGTTCTTTGAGGAAGTTGGGATAGAGACGGATGAAGATTACAACGTCACACTCCCAGACCCC
TGA
```

**Restriction Sites:** Sgfl-Mlul  
**ACCN:** NM\_001256015  
**Insert Size:** 417 bp  
**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).  
**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001256015.1](#)

**RefSeq Size:** 1609 bp

**RefSeq ORF:** 417 bp

**Locus ID:** 29893

**UniProt ID:** [Q9P2W1](#)

**Cytogenetics:** 17q21.2

**Protein Families:** Druggable Genome

**MW:** 16.1 kDa

**Gene Summary:** This gene encodes a protein that functions in meiotic recombination. It is a subunit of the PSMC3IP/MND1 complex, which interacts with PSMC3/TBP1 to stimulate DMC1- and RAD51-mediated strand exchange during meiosis. The protein encoded by this gene can also co-activate ligand-driven transcription mediated by estrogen, androgen, glucocorticoid, progesterone, and thyroid nuclear receptors. Mutations in this gene cause XX female gonadal dysgenesis. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Dec 2011]

Transcript Variant: This variant (4) uses an alternate splice site in an internal exon, and thus differs in its 5' UTR and uses a downstream in-frame start codon, compared to variant 2. The encoded isoform (4) is shorter at the N-terminus, compared to isoform 2. Both variants 4 and 5 encode isoform 4.