

## Product datasheet for SC332289

### ARHGAP22 (NM\_001256026) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ARHGAP22 (NM_001256026) Human Untagged Clone
Tag:	Tag Free
Symbol:	ARHGAP22
Synonyms:	RhoGAP2; RhoGap22
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC332289 representing NM_001256026. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGCCCTTCTGGCCCATCAGGTGTCTGAAGAGGAGCAGGCGTATGCCAGGGGTGGTGCCGGGGAGCGG
GAGAAGGTGCCGGCAACCCCGAGGCGCTCCTGCTCATGGCCAGCTCCCAGCGTGACATGGAGGACTGG
GTGCAGGCCATCCGCCGAGTCATCTGGGCCCGCTGGCGGAGGGATCTTTGGCAGCGCCTAGAGGAA
ACAGTCCACCACGAGCGGAAGTATGGCCCCCGCTGGCGCCCTGCTGGTGGAGCAGTGTGTGGACTTC
ATCCGGGAGCGCGGGCTCACTGAGGAGGGGCTGTTCCGCATGCCAGGCCAGGCCAACCTGGTGGAGGAC
CTGCAGGATTCTTCGACTGTGGGGAGAAGCCACTGTTTGACAGCACAAACAGACGTGCACACGGTGGCC
TCCTGCTGAAGCTGTACCTGCGGGAGCTCCCCGAGCCCGTGGTCCCCTTCGCCAGGTACGAGGACTTC
CTCAGCTGCGCCAGCTGCTCACCAAGGACGAGGGGGAGGGCACTCTGGAGTTGGCTAAACAAGTGAGC
AACCTTCTCAGGCAAATTACAACCTGCTCAGATACATCTGCAAGTTTCTGGATGAAGTTCAGGCATAC
TCAAATGTCAACAAGATGAGTGTCCAGAATCTGGCAACCGTTTTTGGACCTAACATTCTGCGGCCACAG
GTAGAGGACCCAGTAACCATCATGGAAGGCACTCCCTCGTCCAGCACCTGATGACCGTCTCATCCGC
AAACACAGCCAGCTTTACGGCACCCGGTCCCGGAAGGGCCACCTCCCCGCGGGGGCCTGCAATGC
GCAGTGGGGTGGGGCTCCGAGGAGGTACCAGGGACAGCCAAGGAGAGCCCGCGGCCCCCGGCCTGCC
GGGAGCCGGTGCAGCCCTGGGAAGAAGGTGCAGACCCTGCCAGTTGGAAGTCTCCTTCCGGCAGCCG
AGGTCCCTATCGGAAGCCGAAGGGGGCGGCTCATCCCTGGAGGTGCCCATCATCTCCTCCGGCGGG
AACTGGCTTATGAACGGGCTGTCTCCCTGCGCGGACACC CGGGCCTCGTCGGGAGACCGGCTCAAG
GACTCGGGCTCCGTGCAGAGACTCTCCACCTACGACAAATGTGCCCGCGCCGGGCTGGTCCCCGGCATA
CCCAGCGTGGCCAGTATGGCGTGGTCCGGGGCCTCGTCCAGCGAGTGTGCGGTGGGGGGCTCACTCAGC
AGCTGCACGGCCTGCCCGCCAGCGACTCGTCTGCCCGCAGTTCCCTGCACACCGACTGGGCCCTGGAG
CCCTCCCCGCTCCCCAGCAGCAGCGAGGACCCCAAGTCCCTGGACCTGGACCACAGCATGGACGAGGCG
GGCGGGGTGCCAGCAACAGCGAGCCCAGCGAGCCGACAGCCCCACCCGGAACACGCGCGCCGCTCC
GAGGCCTTACAGGGGCTGGTCACTGAGCTCAGGGCCGAGCTGTGCCGCCAGCGGACTGAGTACGAGAGG
AGTGTGAAAAGAATCGAAGAAGGGAGTGTGACCTGAGAAAACGAATGTCCCGTTAGAAGAAGAACTG
GACCAGGAAAAGAAAAATACATCATGCTGAAAATAAAGTGCAGGAACTCTGAACGGGCGCGGGAGGAT
GCGGAGAGGAGGAACAGCTGTTGCAGAGGAAAATGGAGGAGTTTTTTTCGACCCTAGGAAGCTTGACT
GTTGGGGCAAAGGTGCCAGGGCCCCAAAGTAA

```



[View online »](#)

<b>Restriction Sites:</b>	Sgfl-Mlul
<b>ACCN:</b>	NM_001256026
<b>Insert Size:</b>	1827 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001256026.1</a>
<b>RefSeq Size:</b>	2352 bp
<b>RefSeq ORF:</b>	1827 bp
<b>Locus ID:</b>	58504
<b>UniProt ID:</b>	<a href="#">Q7Z5H3</a>
<b>Cytogenetics:</b>	10q11.22-q11.23
<b>MW:</b>	66.6 kDa
<b>Gene Summary:</b>	<p>This gene encodes a member of the GTPase activating protein family which activates a GTPase belonging to the RAS superfamily of small GTP-binding proteins. The encoded protein is insulin-responsive, is dependent on the kinase Akt and requires the Akt-dependent 14-3-3 binding protein which binds sequentially to two serine residues. The result of these interactions is regulation of cell motility. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2011]</p> <p>Transcript Variant: This variant (4) differs in the 5' UTR and lacks a portion of the 5' coding region, uses an alternate in-frame splice site in a 5' exon, and initiates translation at an alternate start codon compared to variant 1. The resulting protein (isoform 4) is shorter and has a distinct N-terminus compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>