

Product datasheet for **SC332312**

ZNF85 (NM_001256171) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF85 (NM_001256171) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF85
Synonyms:	HPF4; HTF1
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC332312 representing NM_001256171. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGGACCATTGACATTTAGGGATGTGGCCATAGAATTCTCTCTGAAGGAGTGGCAATGCCTGGACACT
GCACAGCGGAATTTATATAGAAATGTGATGTTAGAGAAGTACAGAAACCTGGTCTTCTGGGTATTACT
GTTTCTAAGCCAGACCTGATCACTTGTCTGGAGCAAGGGAAGAGGCCTGGAGTATGAAGAGACATGAG
ATCATGGTGGCCAAACCCACAGGTTGTCTATCAAAAAGAATGTCGTATGAGCTATTGAGAAAGGTGTAT
ATCCTGATGTTACTGAGGAGTGTCTCTATGCCTCTGTTAGAATTATGTGTTCTCATTTTGCCTAAGAC
CTTTGGCCGGAGCAGAATATAAAGATTCTTTCCAAAAGTGACACTGAAAAGATATGGAAGATGTAGA
CATGAAAATTTACCATTAAAGAAAGGCTGTGAAAGTATGGATGAGTGTAAGATGCACAAAGAGGTTGT
AATGGACTTAACCAATGTCTCACAGCTACCCAGAGCAAAATATTTCAATGTGATAAATATGTAAAGTC
GCTCATAAATTTCAAATTCAAACAGACATGAGATAAGACATACTAAAAAGAAACCTTTCAAATGTACA
AAATGTGGCAAATCATTTGGCATGATTTTCATGCCTAACTGAACATAGCAGAATTCATACTAGAGTAAAT
TTCTACAAATGTGAAGAATGTGGAAAAGCCTTTAACTGGTCTCAACCTTACTAAACATAAGAGAATT
CATACGGGAGAGAAACCTTACAAATGTGAAGAATGTGGTAAAGCCTTTAACCAGTCTCAACCTTATT
AAACATAAGAAAATTCATACTGGAGAGAAACCTTACAAATGTGAAGAATGTGGCAAAACCTTTAACCAG
TTCTCAACTCTTACTACCCATAAGATAATTCATACTGGAGAGAAACCTTACAAATGTAAAGAATGTGGT
AAAGCTTTTAACCGATCTTCAACCTTACTACCCATAGAAAAATTCATACTGGAGAGAAACCTTACAAA
TGTGAAGAATGTGGCAAAAGCCTTTAAGCAGTCTCAACCTTACTACACATAAGATAATTCATACTGGA
GAGAAACCTTACAAATGTAAAAATGTGGAAAAGCCTTTAACCAGTCTGCACACCTTACCACACATGAG
GTAATTCATACTGGAGAGAAACCTTACAAATGTGAAAAATGTGGAAAAGCCTTTAATCATTTCTCACAC
CTTACTACACATAAGATAATTCATACTGGAGAGAAACCTTACAAATGTAAAGAATGTGGTAAAGCTTTT
AAACACTCTTCAACCTTACTAAACATAAGATAATTCATACTGGAGAGAAAGCCTTACAAATGTAAAGAA
TGTGAAAAAGCTTTTAACCAATCTCAAACTTACTGAACATAAGAAAAATTCATACTGGAGAGAAACCC
TATGAATGTGAAAAATGTGGCAAAGCTTTAACCAGTCTCAAACTTACTAGACATAAGAAAAGTCAT
ACAGAAGAGAAACCTTACAAATGTGAAGAATGTGGCAAAGGTTTTAAATGGCCCTCAACCTTACTATC
CATAAGATAATTCATACTGGAGAGAAACCTTACAAATGTGAAGAATGTGGCAAAGCTTTTAACCAATCC
TCAAACTTACCAACATAAGAAAATTCATACTGGAGAGAAACCTTACACATGTGAAGAATGTGGCAAA
GCCTTTAACCAGTCTCAACCTTACTAAACATAAGAGAATTCATACTGGAGAGAAACCTTACAAATGT
GAAGAATGTGACAAAGCTTTTAATGGTCTCAGTCTTACTAAACATAAGATAATTCATACCGGAGAA
AAATTACAAATGA
  
```



[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001256171
Insert Size:	1878 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001256171.1</u>
RefSeq Size:	2435 bp
RefSeq ORF:	1878 bp
Locus ID:	7639
Cytogenetics:	19p12
Protein Families:	Druggable Genome, Transcription Factors
MW:	72.3 kDa
Gene Summary:	May be a transcriptional repressor.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (3) represents the longest transcript and encodes the longest isoform (c).