

## Product datasheet for **SC332368**

### **E2F8 (NM\_001256371) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	E2F8 (NM_001256371) Human Untagged Clone
Tag:	Tag Free
Symbol:	E2F8
Synonyms:	E2F-8
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

**Fully Sequenced ORF:** >SC332368 representing NM\_001256371.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGAGAACGAAAAGGAAAATCTCTTTTGTGAGCCACATAAAAGGGGACTAATGAAAACACCTCTGAAA
GAATCCACCACAGCAAATATCGTGTGGCAGAGATCCAGCCTGACTTTGGCCCTTAACCCACCTACC
AAGCCCAAGGAAGGCTCTCAGGGAGACCGTGGACACCCGACCAACCTGAAAATGCTCATCAGTGCT
GTGAGCCCTGAGATCCGCAACAGAGATCAGAAAAGGGGTTTGTGTTGACAACAGAAGTGGATTACCTGAG
GCCAAAGACTGTATACACGAACACTTATCTGGAGATGAATTTGAGAAATCCCAACCAAGTCGAAAAGAG
AAAAGTTTAGGATTATTGTGTCATAAGTCTTAGCACGATATCCTAATTATCCCAACCTGCTGTGAAT
AATGACATCTGCCTTGACGAAGTGGCAGAGGAAGTAAATGTTGAACGTCGACGATTTACGATATCGTG
AACGTCCTAGAGAGTTTACATATGGTGAGCCGCTCGCCAAAAACAGGTACACTTGGCACGGGCGACAC
AATCTCAACAAAACCTTGGCACCTTGAAGAGCATCGGGGAGGAGAATAAGTACGCCGAGCAGATTATG
ATGATCAAAAAGAAAGAATATGAGCAAGAGTTTGACTTTATTAAGAGTTACAGTATAGAGGATCATATC
ATCAAATCAAACACTGGCCCAAATGGACACCCAGACATGTGTTTTGTGGAACCTCCCTGGAGTGAATTT
CGGGCAGCTTCTGTAACAGCCGCAAGACAAGTCTTTAAGGGTAAATGAGCCAGAAATTTGTGATGCTG
TTTTTGGTGTCAACGCCTCAGATAGTAAGCCTAGAAGTTGCTGCCAAGATTTAATTGGGGAGGACCAT
GTGGAAGATTTGGATAAAAGCAAGTTTAAAACAAAATTAGGAGGTTGTATGATATAGCTAATGTTCTG
AGTAGCCTGGATCTTATCAAGAAAGTTCATGTTACAGAGGAAAGAGGCCGAAAACAGCTTTCAAATGG
ACCGGCCCAGAAATCAGTCCAAATACCAGTGGCTCCAGCCAGTCATTCAATTTACTCCCTCTGATTTG
GAGGTGAGACGGTCTTCAAAGAGAAGTGTGCCAAAACCTCTTTTCCACACGTGGGAAACCAAACCTT
ACTCGACACCCATCTCTTATCAAATGGTAAAGAGTATAGAAAGTATAGAAAGTATCGGAGAAAGATAAATTTCTGCG
CCCAGTAGCCCTATCAAGACCAACAAAGCTGAGAGTTCAGAAATCTGCACCTTCCCAAGTAAAATG
GCTCAGCTCGCAGCTATTTGTAAGTGCAGTTAGAAGAGCAATCAAGTGAATCCAGACGAAAGTGAAA
GTACAGCTGGCAAGATCTGGACCTGCAAACAGTAGCCCTCTGGACCCCAAGTGAATGCTGAGATG
GAGCTGACAGCACCGTCCCTCATCCAGCCCTGGGAATGGTCCCTGATCCCCAGCCCTTGTCTATCA
GCAGTGCCCTGATCCTACCTCAGGCCCTTCCAGCCATCCTATGCCATCTACCTGCAGCCACTCAA
GCCACCAAAGTGTGACGCCACCCCAAGGCCTGAGCCCAACGGTGTGCACCACCCACTCTTCTAAAGCT
ACTGGCTCAAAGACTCCACAGATGCCACCCTGAGAAGGCAGCCAATGATACCTCAAAGGCCAGTGCC
TCTACCAGGCCTGGAAGCTTGTGCCAGCACCAGAGAGGCAAGGGGCAAAGAGCCGAACCCAGGGAGCCA
GCTGGAGAAAGAGGCTCAAAGAGGGCAAGCATGCTCGAGGACAGTGGTTCCAAAAGAAAATTTAAAGAG
GACCTAAAAGACTTGAAAATGTCTCCGCAACCTTGTCCCATCAGGATACCTAATCCCTCTCACGCAG
TGCTCATCCCTGGGGCAGAGTCCATTTGTCTGGTAAAGAAAATCAAGTGTCTTTCCCAAACCCAC
AGGATTTACAGCTCCCAATTGCAGGTGTTATTCCAGTGACATCATCTGAACTCACTGCTGTTAATTTT
CCCTCTTTTATGTAACACCGTTGAAGCTAATGGTCTCACCAACTTCCGTGGCAGCCGTACCTGTGCGG
AACAGCCCGGCTCTCGTTCAAGCCACCCTGTTCCCATCCAGAACCCAAGCTCAGCCATTGTAACCTT
ACCCTGCAGCACCTGGGACTCATCTACCCAATGTGCAGTTGTCTGCCAGCCCTGGGTCTGGAATCGTT
CCTGTGTCTCAAAGAATAGAGTCTGTTAATGTGCGACCAGAAAATGCAGGCACTCAGCAAGGAAGGGCC
ACCAACTATGACTCACCAGTCCCAGGCCAGAGCCAGCCAAATGGACAATCAGTTGCTGTGACAGGGGCA
CAACAGCTGTTCTGTGACACCCAAAGGGTCAACAATTAGTGGCCGAAAGTTTCTTCCGTACCCAGGT
GGACCCACCAAGCAACCCAGCTCATCTGCATGGATTTTGAGGGTGTAAATAAACCTCCTTAGGAACT
CTCTTTGTCCACAGCGAAAACCTGGAAGTCTCAACAGAGGATGTCCATTA
  
```

**Restriction Sites:** Sgfl-Mlul

**ACCN:** NM\_001256371

**Insert Size:** 2604 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001256371.1</u>
<b>RefSeq Size:</b>	3537 bp
<b>RefSeq ORF:</b>	2604 bp
<b>Locus ID:</b>	79733
<b>UniProt ID:</b>	<u>A0AVK6</u>
<b>Cytogenetics:</b>	11p15.1
<b>Protein Families:</b>	Druggable Genome
<b>MW:</b>	94.2 kDa
<b>Gene Summary:</b>	<p>This gene encodes a member of a family of transcription factors which regulate the expression of genes required for progression through the cell cycle. The encoded protein regulates progression from G1 to S phase by ensuring the nucleus divides at the proper time. Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Jan 2012]</p> <p>Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Variants 1, 2 and 3 encode the same protein. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>