

Product datasheet for **SC332403**

Cell adhesion molecule 2 (CADM2) (NM_001256503) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cell adhesion molecule 2 (CADM2) (NM_001256503) Human Untagged Clone
Tag:	Tag Free
Symbol:	Cell adhesion molecule 2
Synonyms:	IGSF4D; Necl-3; NECL3; SynCAM 2; synCAM2
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC332403 representing NM_001256503. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGCCTGTCAAACCTCCAAGGCATATCTCACCGTTCTGGGTGTTCTGAAAAGCCTCAGATTAGTGGA
TTCTCATCACCAGTTATGGAGGGTGACTTGATGCAGCTGACTTGCAAACATCTGGTAGTAAACCTGCA
GCTGATATAAGATGGTTCAAAAATGACAAAGAGATTAAGATGTAATAATTTAAAGAAGAGGATGCA
AATCGCAAGACATTCAGTGTGAGCAGCAGCAGTGGACTTCCGAGTGGACCGGAGTGATGATGGAGTGGCG
GTCATCTGCAGAGTAGATCACGAATCCCTCAATGCCACCCCTCAGGTAGCCATGCAGGTGCTAGAAAATA
CACTATACACCATCAGTTAAGATTATACCATCGACTCCTTTTCCACAAGAAGGACAGCCTTTAATTTTG
ACTTGTGAATCCAAAGGAAAACCACTGCCAGAACCTGTTTTGTGGACAAAGGATGCGGAGAAATACCA
GATCCTGACCGAATGGTTGTGAGTGGTAGGGAGCTAAACATTCTTTTCTGAACAAAACGGATAATGGT
ACATATCGATGTGAAGCCACAAACACCATTGGCCAAAGCAGTGCGGAATATGTTCTCATTGTGCATGAT
GTTCCCAACACTTTGCTTCCCACTACTATCATCCCTCCCTTACCCTGCAACAGTCACAACCACTGTA
GCCATAACAACAGCCCAACCACTGCAACAACCAAGCAGCAGCATCAGAGATCCTAATGCTTTGGCTGGC
CAGAATGGCCCTGACCATGCTCTCATAGGAGGAATAGTGGCTGTAGTTGTATTTGTCACGCTGTGTTCT
ATCTTTCTGCTTGGTCGATATCTGGCAAGGCATAAAGGAACGTATTTAACAATGAAGCTAAAGGAGCT
GAAGATGCACCAGATGCTGATACAGCCATTATCAATGCTGAAGGCAGCCAAGTCAATGCTGAAGAGAAA
AAAGAGTATTTTCATTAA
```

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001256503
Insert Size:	984 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001256503.1](#)

RefSeq Size: 8977 bp

RefSeq ORF: 984 bp

Locus ID: 253559

UniProt ID: [Q8N3J6](#)

Cytogenetics: 3p12.1

Protein Families: Druggable Genome, Transmembrane

MW: 35.4 kDa

Gene Summary: This gene encodes a member of the synaptic cell adhesion molecule 1 (SynCAM) family which belongs to the immunoglobulin (Ig) superfamily. The encoded protein has three Ig-like domains and a cytosolic protein 4.1 binding site near the C-terminus. Proteins belonging to the protein 4.1 family crosslink spectrin and interact with other cytoskeletal proteins. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2012]

Transcript Variant: This variant (5) differs in the 5' UTR and coding region and uses a downstream start codon compared to variant 1. The resulting protein (isoform 4) is shorter and has a shorter N-terminus compared to isoform 1. Variants 4 and 5 encode the same protein. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.