

Product datasheet for **SC332420**

Retinoid X Receptor gamma (RXRG) (NM_001256570) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Retinoid X Receptor gamma (RXRG) (NM_001256570) Human Untagged Clone
Tag:	Tag Free
Symbol:	Retinoid X Receptor gamma
Synonyms:	NR2B3; RXRC
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC332420 representing NM_001256570. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGAAC TACCCATCCACCAGCCCGGATCTCTGGTTAAACACATCTGTGCCATCTGTGGAGACAGATCC
TCAGGAAAGCACTACGGGTATACAGTTGTGAAGGCTGCAAAGGGTTCTCAAGAGGACGATAAGGAAG
GACCTCATCTACACGTGTCGGGATAATAAAGACTGCCTCATTGACAAGCGTCAGCGCAACCGTGCCAG
TACTGTCGCTATCAGAAGTGCCTTGTGTCATGGGCATGAAGAGGGAAGCTGTGCAAGAAGAAAGACAGAGG
AGCCGAGAGCGAGCTGAGAGTGAGGCAGAATGTGCTACCAGTGGTCATGAAGACATGCCTGTGGAGAGG
ATTCTAGAAGCTGAACCTGCTGTTGAACCAAAGACAGAATCCTATGGTGACATGAATATGGAGAACTCG
ACAAATGACCCTGTTACCAACATATGTCATGCTGCTGACAAGCAGCTTTTCACCCCTGTTGAATGGGCC
AAGCGTATTCCTCACTTCTCTGACCTCACCTTGGAGGACCAGGTCATTTTCTCGGGCAGGGTGGAAAT
GAATTGCTGATTGCTCTTTCTCCCACCGCTCAGTTTCCGTGTCAGGATGGCATCCTTCTGGCCACGGGT
TTACATGTCCACCGGAGCAGTGCCACAGTGTGGGGTCCGCTCCATCTTTGACAGAGTCCCTAACTGAG
CTGGTTTCCAAAATGAAAGACATGCAGATGGACAAGTCGGAACGGGATGCTGCGAGCCATTGTACTC
TTAACCAGATGCCAAGGGCTGTCCAACCCCTCTGAGGTGGAGACTCTGCGAGAGAAGGTTTATGCC
ACCCTTGAGGCTACACCAAGCAGAAGTATCCGGAACAGCCAGGAGGTTTGCCAAGCTGCTGCTGCGC
CTCCCAGCTCTGGTTCATTGGCTTCAAATGCCTGGAGCACCTTTCTTCTTCAAGCTCATCGGGGAC
ACCCCATTTGACACCTTCTCATGGAGATGTTGGAGACCCCGCTGCAGATCACCTGA
```

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001256570
Insert Size:	1023 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online >](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001256570.1](#)

RefSeq Size: 2263 bp

RefSeq ORF: 1023 bp

Locus ID: 6258

Cytogenetics: 1q23.3

Protein Families: Druggable Genome, Nuclear Hormone Receptor, Transcription Factors

Protein Pathways: Adipocytokine signaling pathway, Non-small cell lung cancer, Pathways in cancer, PPAR signaling pathway, Small cell lung cancer, Thyroid cancer

MW: 38.6 kDa

Gene Summary: This gene encodes a member of the retinoid X receptor (RXR) family of nuclear receptors which are involved in mediating the antiproliferative effects of retinoic acid (RA). This receptor forms dimers with the retinoic acid, thyroid hormone, and vitamin D receptors, increasing both DNA binding and transcriptional function on their respective response elements. This gene is expressed at significantly lower levels in non-small cell lung cancer cells. Alternatively spliced transcript variants have been described. [provided by RefSeq, Jun 2010]
Transcript Variant: This variant (3) includes an additional exon in its 5' UTR, and uses a downstream in-frame start codon, compared to variant 1. The encoded isoform (c, also known as RXRgamma2) is shorter at the N-terminus, compared to isoform a. Both variants 3 and 4 encode isoform c.