

Product datasheet for **SC332424**

TRIM46 (NM_001256599) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TRIM46 (NM_001256599) Human Untagged Clone
Tag:	Tag Free
Symbol:	TRIM46
Synonyms:	GENEY; TRIFIC
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

Fully Sequenced ORF: >SC332424 representing NM_001256599.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGAAGAACATGGAGAAGGAAGTCTGTGCCAGTGTGTCAAGAGATGTACAAGCAGCCACTGGTGCTG
CCCTGTACCCACAACGTGTGCCAGGCTGTGCCGAGAGGTCTTGGGCCAGCAGGGCTACATAGGACAT
GGTGGGGACCCAGCTCCGAGCCACCTCTCCTGCCCTCCACCCCTCCACCCGACGCCCCGCCTCTCC
CGCAGAAGTCTCCCAAGCCAGCCGCTGGACCGGCTGCTTAAGTCAGGCTTTGGGACATACCCCTGGG
AGGAAGCGAGGTGCTTTGCACCCCAAGTGATCATGTTCCCGTGCCAGCCTGCCAAGGTGATGTGGAG
CTTGGGGAGCGGGTCTGGCAGGCTTTCCGGAACCTGACCCTGGAGCGTGTGGTGGAGCGGTACCCGC
CAGAGTGTGAGTGTGGGAGGTGCCATCCTGTGCCAGTTGTGCAAGCCCCACCACTAGAGGCCACCAAG
GGCTGCACAGAGTGCCGCGCCACCTTCTGCAATGAGTGTCAAGCTCTCCACCCCTGGGGACCCAG
AAGGCCCAGCATGAGCCACCTGCCTACCCTCTCCTCCGACCCAAGGGCTTATGTGCCAGACCAC
AAGGAAGAGGTGACCCACTACTGCAAGACATGCCAACGCCTGGTATGTCAACTCTGCCGGTGGCGGC
ACCCACAGCGGGACAAGATCACACCAGTGTCTAGTGCCTACCAGGCCCTCAAGGACAAGTGACAAAG
AGCCTGACATACATCCTGGGAAACCAGGACACGGTACAGACCCAGATCTGTGAGCTGGAGGAGCCGTG
AGGCACACCGAGGTGAGTGGTACAGCAGGCAAGGAGGAGGTGTCCAGCTGGTGGGGGGCTGGGGGCT
GTGCTGGAGGAGAAGCGGGCATCACTGCTTCAGGCCATTGAAGAATGCCAGCAGGAGCGGCTGGCCCGT
CTCAGCGCCAGATCCAGGAGCACCGGAGCCTGTGGATGGCTCAGGCTGGTGGGCTATGCCCAGGAA
GTACTTAAGGAAACAGACCAGCCTTGTCTTGTGCAAGCCGCAAGCAGCTGCACAACAGGATTGCCCGA
GCCACTGAAGCCCTCCAGACATTCGGCCAGCTGCCAGCTCCTCCTCCGCCATTGCCAGCTCGACGTG
GGACGTGAGATGAAGCTGCTGACAGAGCTTAACCTCCTGCGAGTGCCTGAGGCCCCCGTCATTGACACC
CAGCGCACCTTTGCCTATGATCAGATCTCCTGTGCTGGCGCTGCCCCCAATTCACCACCTGCCTGG
CACTATACCGTTGAGTTCGGCGCACGGATGTGCCTGCTCAGCCAGGCCCAAGCCGCTGGCAGCGGGCGG
GAGGAGGTGAGGGGACCCAGTGCCTGCTTGAAGACCCGACACGGGCTCTGTGTATGTGCTGCGTGTG
CGCGGCTGCAACAAGCCGGCTACGCGCAATACAGTGAAGATGTGCACCTGCACACGCCCCGGCACCT
GTCCTGCCTTCTCCTCGATAGCCGCTGGGGCGCAAGCCGAGAGCGGCTGGCTATCAGCAAGGACCAG
CGAGCAGTACGGAGTGTCCAGGGCTGCCCTGCTGCTGGCTGCTGACCGGCTGCTGACCGGCTGCCAC
CTGAGTGTGGATGTGGTCTGGGCGACGTGGCTGTGACCCAGGGCCGACGCTACTGGGCTGCGCCGTA
GACCCAGCCTCCTACTTGGTCAAGGTGGGCGTGGGCTGGAGAGCAAGCTTCAAGAAAGTTCCAGGGT
GCCCCGATGTGATCAGCCCAGGTACGACCCGGACAGCGGGCACGACAGCGGTGCCGAGGATGCCACA
GTGGAGGGCTGCGCACCTTCGCTTTCCTAACCATTGGCATGGGCAAGATCCTGTGGGGTGGGGGCA
AGCTCAAACGCAGGGCTGACAGGGAGGGATGGCCACAGCCGGCTGCACAGTGCCTGCCACCCCGC
CTGGGCATCTGCCTGGACTATGAGCGGGGCCGGGTTTCCTTCTGGATGTGTTTCTTCCGTGGGCTC
TTGGAGTGCCCTGGACTGCTCAGGGCCTGTGTGCCCTGCCTTTTGTTCATCGGGGTGGCGCAGTA
CAGCTCCAGGAGCCAGTGGGCACTAAGCCTGAGAGGAAAGTACCATTGGGGGCTTCGCCAAGCTGGAC
TGA
  
```

Restriction Sites: Sgfl-Mlul

ACCN: NM_001256599

Insert Size: 2211 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001256599.1
RefSeq Size:	3251 bp
RefSeq ORF:	2211 bp
Locus ID:	80128
UniProt ID:	Q7Z4K8
Cytogenetics:	1q22
Protein Families:	Druggable Genome
MW:	80.9 kDa
Gene Summary:	<p>Microtubule-associated protein that is involved in the formation of parallel microtubule bundles linked by cross-bridges in the proximal axon. Required for the uniform orientation and maintenance of the parallel microtubule fascicles, which are important for efficient cargo delivery and trafficking in axons. Thereby also required for proper axon specification, the establishment of neuronal polarity and proper neuronal migration.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (2) differs in the 5' UTR and initiates translation at a downstream, in-frame start codon, compared to variant 1. The encoded isoform (2) has a shorter N-terminus, compared to isoform 1.</p>