

Product datasheet for **SC332437**

ZNF43 (NM_001256650) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF43 (NM_001256650) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF43
Synonyms:	HTF6; KOX27; ZNF39L1
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

Fully Sequenced ORF: >SC332437 representing NM_001256650.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGATGTGGCCATAGAATTCTGTCTGGAGGAGTGGCAATGCCTGGACATTGCACAGCAGAATTTATAT
AGGAATGTGATGTTAGAGAACTACAGAACTGGTCTTCTGGGTATTGCTGTCTCTAAGCCAGACCTG
ATCACCTGTCTGGAGCAAGAAAAAGACCTTGGGAGCCTATGAGGAGACATGAAATGGTAGCCAAACCC
CCAGTTATGTGTTCTCATTTTACCCAAGACTTTTGGCCAGAGCAGCATATAAAAGATCCTTTCCAAAAA
GGGACACTGAGAAGATATAAAACTGTGAACATAAAAATGTACATTTAAAAAAGACCATAAAAAGTGTG
GATGAGTGAAGGTGCACAGAGGAGTTATAATGGATTTAACCAATGTTTGCCAGCTACCCAGAGCAAA
ATATTTCTATTTGATAAATGTGTGAAAGCCTTTCATAAATTTTCAAATCAAACAGACATAAGATAAGC
CATACTGAAAAAAACTTTTCAAATGCAAAGAATGTGGCAATCATTTTGCATGCTTCCACATCTAGCT
CAACATAAAATAATTCAACCAGAGTGAATTTCTGCAAAATGTGAAAAATGTGGAAAAGCTTTTAACTGC
CCTTCAATCATCACTAACATAAGAGAATTAATACTGGAGAGAAACCTACACATGTGAAGAATGTGGC
AAAGTCTTTAATTGGTCCTCACGCCTTACTACACATAAAAAAATATACTAGATACAAACTCTACAAA
TGTGAAGAATGTGGCAAAGCTTTTAAACAAGTCTCAATCCTTACTACCCATAAGATAATTCGCACTGGA
GAGAAATCTACAAATGTAAAGAATGTGCCAAAGCTTTTAAACCAATCCTCAAACCTTACTGAACATAAG
AAAATTCATCCTGGAGAGAAACCTTACAAATGTGAAGAATGTGGCAAAGCCTTTAACTGGCCCTCAACT
CTTACTAAACATAAGAGAATTCATACTGGAGAGAAACCTACACATGTGAAGAATGTGGCAAAGCCTTT
AACCAGTTCTCAAACCTTACTACACATAAGAGAATCCATACTGCAGAGAAATTCATAAAATGTACAGAA
TGTGGTGAAGCTTTTAGCCGGTCTCAAACCTTACTAAACATAAGAAAATTCATACTGAAAAGAAACCC
TACAAATGTGAAGAATGTGGCAAAGCTTTTAAAGTGGTCTCAAAGCTTACTGAACATAAGTTAACTCAT
ACTGGAGAGAAACCTACAAATGTGAAGAATGTGGCAAAGCCTTTAACTGGCCCTCAACCCCTTACTAAA
CATAACAGAATTCATACTGGAGAGAAACCTACAAATGTGAAGTATGTGGCAAAGCCTTTAACCAGTTC
TCAAACCTTACTACACATAAGAGAATTCATACTGCAGAAAAACCTACAAATGTGAAGAATGTGGCAAA
GCTTTTAGCCGGTCTCAAACCTTACTAAACATAAGAAAATTCACATTGAAAAGAAACCTACAAATGT
GAAGAATGTGGCAAAGCTTTTAAAGTGGTCTCAAAGCTTACTGAACATAAGATAACTCATACTGGAGAG
AAACCTACAAATGTGAAGAATGTGGCAAAGCTTTTAAACATTTCTCAATCCTTACCAAACATAAGAGG
ATTCATACTGGAGAGAAACCTACAAAGTGTGAAGAATGTGGCAAAGCTTTTACCCAATCTCAAACCTT
ACTACACATAAGAAAATTCATACTGGAGAGAAATTCATAAAATGTGAAGAATGTGGCAAAGCTTTTACC
CAATCTTCAAACCTTACTACACATAAAAAAATTCATACTGGAGGAAAACCTACAAATGTGAAGAATGT
GGCAAAGCTTTTAAACAGTTCTCAACTTACTAAACATAAGATAATTCACACTGAGGAGAAACCTAC
AAATGTGAAGAATGTGGCAAAGCTTTTAAAGTGGTCTCAAACCTTACTAAACATAAGATAATTCATACT
GGAGAGAAACCTACAAATGTGAAGAATGTGGCAAAGCTTTTAAACTGTCCTCAAACCTTTCTACACAT
AAGATTATTCATACTGGAGAGAAACCTACAAATGTGAAAAATGTGGCAAAGCTTTTAAACCGATCCTCA
AACCTTATTGAACATAAGAAAATTCATACTGGAGAGCAACCTACAAATGTGAAGAATGTGGCAAAGCA
TTTAACTATTCTCACACCTTAATACACATAAGAGAATTCATACTAAAGAGCAACCTACAAATGTAAA
GAATGTGGCAAAGCTTTCAACCAATATTCAAACCTTACTACACATAACAAAATTCATACTGGAGAGAAA
CTCTACAAACCTGAAGATGTGACAGTGATTTTGACAACACCTCAAACCTTTTCAAACATAAAATAA
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001256650

Insert Size: 2412 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001256650.1](#)

RefSeq Size: 5508 bp

RefSeq ORF: 2412 bp

Locus ID: 7594

UniProt ID: [P17038](#)

Cytogenetics: 19p12

Protein Families: Transcription Factors

MW: 93.5 kDa

Gene Summary: This gene belongs to the C2H2-type zinc finger gene family. The zinc finger proteins are involved in gene regulation and development, and are quite conserved throughout evolution. Like this gene product, a third of the zinc finger proteins containing C2H2 fingers also contain the KRAB domain, which has been found to be involved in protein-protein interactions. [provided by RefSeq, Jul 2008]

Transcript Variant: This variant (3) differs in the 5' UTR compared to variant 1. Variants 1, 2 and 3 encode the same isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.