

## Product datasheet for **SC332454**

### SSBP2 (NM\_001256732) Human Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** SSBP2 (NM\_001256732) Human Untagged Clone  
**Tag:** Tag Free  
**Symbol:** SSBP2  
**Synonyms:** HSPC116; SOSS-B2  
**Vector:** pCMV6-Entry (PS100001)  
**Fully Sequenced ORF:** >SC332454 representing NM\_001256732.  
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGTACGGCAAAGGCAAGAGTAACAGCAGCGCCGTCCCGTCCGACAGCCAGGCCCGGGAGAAGTTAGCA
CTCTACGTATATGAATATCTGCTCCATGTAGGAGCTCAGAAATCAGCTCAAACATTTTTATCAGAGATA
AGATGGGAAAAAACATCACATTGGGGAAACCACCAGGATTCTTACATTCTTGGTGGTGTGATTTTTGG
GATCTCTACTGTGCAGCTCCAGAGAGACGTGAAACATGTGAACACTCAAGTGAAGCAAAAGCCTCCAT
GATTACAGTGTGCAGCAGCTCCAGTCCAGTGTAGGAAACATCCCCCAGGAGATGCCATGCCAGTA
GGTCTGTACCACCAGGGTTCTTTCAGCCTTTTATGTCACCTCGTACCCTGGAGGTCCAAGGCCCCCA
TTGAGGATACCTAATCAGGCACTTGGAGGTGTCCAGGAAGTCAGCCATTACTCCCCAGTGGAAATGGAT
CCAACCTGACAACAAGGACATCCAAATATGGGTGGGCAATGCAGAGAATGACTCCTCCAAGAGGAATG
GTGCCCTTAGGACCACAGTCTGACCCTTGGTTATCATTACAGAATATGGAGGTGCAATGAGACCCCA
CTGAATGCTTTAGGTGGCCCTGGAATGCCTGGAATGAACATGGGTCCAGGTGGTGGTAGACCTTGGCCA
AACCCAACAAATGCCAATTCATACCATACTCCTCAGCATCTCCTGGGAATTATGTAGGTCTCCAGGA
GGTGGAGGGCCACCAGGAACACCCATCATGCCTAGTCCAGCAGATTCAACCAACTCTGGTGATAACATG
TATACTTTAATGAATGCAGTACCTCCTGGACCTAACAGACCTAATTTTCCAATGGGTCTGGGTGAGAT
GGTCCCATGGGTGGATTAGGAGGAATGGAGTCACATCACATGAATGGCTCTTTAGGCTCAGGAGATATG
GACAGTATTTCCAAGAATTCTCCAATAATATGAGCCTGAGTAATCAACCGGGCACTCCAAGGGATGAT
GGCGAAATGGGGGAAATTTCTTAAATCCTTTTCAGAGTGAGAGTACTCCCCTAGCATGACAATGAGC
GTGTGA
```

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001256732  
**Insert Size:** 1110 bp  
**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001256732.1</a></u>
<b>RefSeq Size:</b>	4465 bp
<b>RefSeq ORF:</b>	1110 bp
<b>Locus ID:</b>	23635
<b>UniProt ID:</b>	<u><a href="#">P81877</a></u>
<b>Cytogenetics:</b>	5q14.1
<b>Protein Families:</b>	Transcription Factors
<b>MW:</b>	38.8 kDa
<b>Gene Summary:</b>	<p>This gene encodes a subunit of a protein complex that interacts with single-stranded DNA and is involved in the DNA damage response and maintenance of genome stability. The encoded protein may also play a role in telomere repair. A variant of this gene may be associated with survival in human glioblastoma patients. [provided by RefSeq, Sep 2016]</p> <p>Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>