

Product datasheet for **SC332459**

CNOT10 (NM_001256741) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CNOT10 (NM_001256741) Human Untagged Clone
Tag:	Tag Free
Symbol:	CNOT10
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

Fully Sequenced ORF: >SC332459 representing NM_001256741.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGGCTGCAGACAAGCCTGCAGATCAGGGAGCAGAGAAACATGAAGGCACAGGTCACTCTGGGATC
ACTGATCAAGAGAAGGAGTTATCCACCAATGCTTTCCAAGCTTTCACATCTGGAAATTATGATGCCTGT
CTACAACACCTTGCCTGTCTACAAGATATAAACAAAGATGATTAAAAATAATTTTGAATACAGCAGTA
GCTGAGTTTTTTAAAAGTAACCAAAACAACAACAGATAATTTGAGACAAACACTTAACCAGCTGAAGAA
CAGGTCCACTCAGCTGTTGAAGAAATGGATGGATTAGATGATGTTGAAAACAGCATGTTGTACTATAAT
CAAGCAGTCATTCTTTATCATCTGCGGCAGTATACAGAAGCCATATCAGTTGGTAAAAACTTTATCAG
TTCATAGAGCCTTTTGAAGAAAAATTTGCCAAGCAGTGTGTTTTTGTCTGTAGACCTGTATATATTA
ACCTACCAAGCTGAGAAAGCTTTGCATCTTCTTGTCTCCTAGAAAAATGATTTACAGGGTAACAAT
AACAAAAATGGAAGAATGAGACTGGTAATAACAACAACAAGATGGATCTAATCATAAAGCTGAAAGT
GGAGCTCTAATAGAAGCTGCAAAATCAAAGATACATCAGTACAAAGTACGAGCTTATATCCAAATGAAG
TCTCTGAAAGCATGTAAGGGAAATCAAGTCAGTCATGAATACAGCTGGAATTCGACCCCTCTCTC
TTTCTTAAAAGCAATTTTGAAGTACTTAAGAGGTAATTATCGAAAAGCCGGAAGCTATTAATAGTTCA
AACATTGCTGAGCATCCAGGATTCATGAAAACAGGTGAATGCTTGAGATGCATGTTCTGGAATAACCTT
GTTTGCATCCATTTTGCATGAGCAAGCACAATTTGGGAATATTCTACTTTAAAAAGGCTCTGCAAGAG
AATGACAAATGTCTGTGCACAGCTCAGTGCAGGTAGCACTGATCCAGGTAAAAATTTTTCAGGAAGACCC
ATGTGTACGTTACTAACCAATAAGAGATATGAGTTGCTGTATAACTGTGGAATTCAGCTTCTTACATT
GGAAGGCTCTTGTCTGCCTTGAATGTCTGATTGAAGCTGTTCAAGTTTATCATGCAAAATCCTCGCCTC
TGGCTACGGCTGGCTGAATGTCTGATTGCTGCAATAAGGGGACTTCTGAACAAGAACTAAAGGCTT
CCCAGCAAAAAAGGAATTGTACAGTCTATTGTTGGTCAAGGCTATCATCGTAAAAATAGTTTTGGCATCA
CAGTCTATACAGAATACTGTTTATAATGATGGGCAGTCTTCGGCCATTCTGTAGCCAGTATGGAGTTT
GCAGCCATATGTCTCAGAAATGCCTTGTGCTGCTACCTGAAGAACAGCAAGATCCAAAGCAGGAAAAT
GGGGCTAAAAATAGTAATCAATTAGGTGGGAACACAGAGAGCAGCGAAAAGCAGTGAAACTTGCAAGTGC
TCCATACTTGCTTGCAGTGCCTACGTGGCTCTGGCTTTGGGTGATAACCTCATGGCTTTGAATCATGCA
GATAAACTTCTTACAGCAGCCCAAGCTGTGAGGATCTCTTAAGTTTTTGGGACATTTATATGCTGCAGAA
GCCCTCATCTCTCGACAGAATATCTGATGCCATTACTCACTTGAACCCGGAGAATGTCAGTGTGTC
TCCTTAGGGATCTCTTCAAATGAGCAGGACCAAGGATCAGACAAAGGTGAAAATGAAGCAATGGAATCC
TCTGGTAAGCGGGCCCTCAGTGTACCCAGTCCGTCAACTCTGCCAGGACTGTGATGCTGTTCAAC
CTTGGCAGCGCTTACTGCCTGAGGAGCGAATATGACAAAGCCGAAAGTGTCTCCACCAGCGGCTTCA
ATGATCCATCCTAAAGAGGTGCCCCCTGAGGCCATCTTGCTGGCAGTCTACCTTGAAGTGCAGAATGGT
AATACTCAGCTGGCCTTACAGATCATCAAAGGAATCAGCTGCTCCCTGCAGTGAAAAACACTCTGAA
GTGAGAAAGAAGCCAGTGTTCAGCCTGTCCACCCGATCCAGCCCATCAAATGCCGGCTTTCACCACT
GTGCAGAGAAAGTGA
```

Restriction Sites: Sgfl-Mlul

ACCN: NM_001256741

Insert Size: 2154 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001256741.1](#)

RefSeq Size: 2748 bp

RefSeq ORF: 2154 bp

Locus ID: 25904

UniProt ID: [Q9H9A5](#)

Cytogenetics: 3p22.3

Protein Pathways: RNA degradation

MW: 79.4 kDa

Gene Summary: Component of the CCR4-NOT complex which is one of the major cellular mRNA deadenylases and is linked to various cellular processes including bulk mRNA degradation, miRNA-mediated repression, translational repression during translational initiation and general transcription regulation. Additional complex functions may be a consequence of its influence on mRNA expression. Is not required for association of CNOT7 to the CCR4-NOT complex.

[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (2) lacks an exon in the coding region, but maintains the reading frame, compared to variant 1. The encoded isoform (2) is shorter than isoform 1.