

Product datasheet for **SC332462**

PEX5 related protein (PEX5L) (NM_001256751) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PEX5 related protein (PEX5L) (NM_001256751) Human Untagged Clone
Tag:	Tag Free
Symbol:	PEX5 related protein
Synonyms:	PEX5R; PEX5RP; PXR2; PXR2B; TRIP8b
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC332462 representing NM_001256751. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGTACCAGGGACACATGCAGGGAAAAGGCTCTAGGGCGGCAGATAAGGCTGTTGCCATGGTGATGAAG
GAGATACCAGGGAGGAGTCTGTGAAGAAAAGCCCTCCTTACTATGACATCACAGCTGGTGAATGAG
CAACAAGAAAGCAGACCCCTCCTGAGTCCCTCCATCGATGACTTCTCTGTGAAACCAAACCGGAAGCA
ATAGCAAGGCCAGTAACATCCAATACAGCTGTATTGACCACTGGCTTAGATCTCCTCGACCTGAGTGAA
CCAGTCTCTCAAACCCAAACCAAAGCCAAGAAGTCAGAGCCCTCATCAAAAACCTCATCCCTCAAGAAA
AAGGCCGATGGATCTGACCTCATCAGCAGGATGCTGAGCAGAGAGGCCAGCCTCTCAGAGTCCCGGAG
ACTTCATCCTTAGATCTAGACATTCAAACACAACCTGGAATAATGGGACGATGTTAAGTTTCATGGAGAT
CGAAATACCAAGGGACATCCAATGGCAGAGAGAAAATCATCCTCATCTAGAAGTGGATCAAAAGAGCTC
TTATGGTCCTCAGAACACAGATCTCAACCAGAAGTGGTGGAAAAAGCGCCCTCAACTCTGAGTCG
GCTTCAGAATTGGAATTAGTGGCTCCGACTCAGGCTCGACTGACCAAGAACATCGCTGGGGAAGCGCA
TTACTTTCTAGAAACCACTCCTTAGAAGAAGAGTTTGAAGGGCAAAAAGCAGCAGTGGAGTCAGATACA
GAGTTTTGGGATAAGATGCAAGCAGAATGGGAAGAAATGGCTCGGAGGAAGTGGATATCTGAGAACCAA
GAAGCCCAGAACCAAGTAACCATCTCGGCTAGTGAGAAGGGATATTACTTTCACACTGAAAACCCCTTC
AAGGACTGGCCTGGAGCATTGGAAGAAGGCTTAAAAGGCTGAAGGAAGGGGATCTGCCAGTCACCATC
CTGTTTCATGGAAGCAGCAATTCTCAGGACCCTGGAGATGCAGAGGCATGGCAGTTCCTCGGGATAACC
CAGGCGGAGAATGAAAATGAACAAGCAGCTATTGTGCGCCCTCCAGAGGTGCTTAGAATTACAGCCCAAC
AACTTAAAAGCTTTGATGGCCTTGGCTGTGAGTTATACTAACACTGGCCATCAGCAGGATGCCTGTGAC
GCTCTGAAGAATTGGATTAAGCAAAATCCAAAGTACAAATACCTTGTGAAAAGCAAGAAGGGATCTCCA
GGCCTCACCCGGCGGATGTCTAAGTCCCAGTTGATAGCTCTGTTCTGGAAGGGGTGAAGGAATTATAT
CTGGAAGTGGCCACCAAAATGGAGATATGATCGACCCAGACCTGCAGACAGGTCTAGGGGTTCTGTTT
CACCTGAGTGGAGAATTTAATAGAGCAATAGATGCATTTAACGCTGCCTTAACTGTTTCGGCCAGAGGAC
TATTCATATGGAACCGCTCGGGGCGACCTTGGCGAACGGAGACCGCAGCGAGGAAGCCGTGGAGGCC
TATACGGGAGCACTGGAGATTCAGCCAGGATTCATCCGGTCCAGATACAACCTAGGAATAAGCTGCATC
AACCTGGGCGCCTACAGAGAAGCGGTGACGAATTTTCTACTGCCCTCAGTTTCAAAGAAAAGAGCAGG
AATCAGCAGCAAGTTCCTCATCTGCAATCTCTGGGAATATCTGGGCTGCCCTCAGAATTGCGCTCTCT
CTGATGGACCAACCAGAACTCTCCAGGCGGCTAATCTTGGTGACCTGGATGTCTCTTAAGAGCTTTC
AACTTGGATCCTTGA

```



[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001256751
Insert Size:	1809 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001256751.1
RefSeq Size:	9133 bp
RefSeq ORF:	1809 bp
Locus ID:	51555
UniProt ID:	Q8IYB4
Cytogenetics:	3q26.33
Protein Families:	Druggable Genome
MW:	67 kDa
Gene Summary:	<p>Accessory subunit of hyperpolarization-activated cyclic nucleotide-gated (HCN) channels, regulating their cell-surface expression and cyclic nucleotide dependence.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (3) encodes isoform 3. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>