

Product datasheet for **SC332463**

PEX5 related protein (PEX5L) (NM_001256752) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PEX5 related protein (PEX5L) (NM_001256752) Human Untagged Clone
Tag:	Tag Free
Symbol:	PEX5 related protein
Synonyms:	PEX5R; PEX5RP; PXR2; PXR2B; TRIP8b
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC332463 representing NM_001256752. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGTACCAGGGACACATGCAGAAAAGTAAAGAACAAGGATATGGAAAACCTAAGCAGTGATGAAGACCTC
GAAATAATTGTTGATCAAAAGCAGCTGGTGAATGAGCAACAAGAAAGCAGACCCCTCCTGAGTCCCTCC
ATCGATGACTTTCTGTGAAACCAATCGGAAGCAATAGCAAGGCCAGTAAACATCCAATACAGCTGTA
TTGACCACTGGCTTAGATCTCCTCGACCTGAGTGAACCACTCTCAAACCCAAACCAAGCCAAGAAG
TCAGAGCCCTCATCAAAAACCTCATCCCTCAAGAAAAAGGCCGATGGATCTGACCTCATCAGCACGGAT
GCTGAGCAGAGAGGCCAGCCTCTCAGAGTCCCGGAGACTTCATCCTTAGATCTAGACATTCAAACACAA
CTGGAAAAATGGGACGATGTTAAGTTTCATGGAGATCGAAATACCAAGGGACATCCAATGGCAGAGAGA
AAATCATCCTCATCTAGAACTGGATCAAAAAGAGCTTTATGGTCCTCAGAACACAGATCTCAACCAGAA
CTGAGTGGTGGAAAAAGCGCCCTCAACTCTGAGTCGGCTTCAAGATTGGAATTAGTGGCTCCGACTCAG
GCTCGACTGACCAAGAACATCGCTGGGGAAGCGCATTACTTTCTAGAAACCACTCCTTAGAAGAAGAG
TTTGAAGGGCAAAAGCAGCAGTGGAGTACAGATACAGAGTTTGGGATAAGATGCAAGCAGAATGGGAA
GAAATGGCTCGGAGGAAGTGGATATCTGAGAACCAAGAAGCCAGAACCAAGTAACCATCTCGGCTAGT
GAGAAGGGATATTACTTTCACACTGAAAACCCCTTCAAGGACTGGCCTGGAGCATTGGAAGAAGCCTTA
AAAAGGCTGAAGGAAGGGGATCTGCCAGTCACCATCCTGTTTATGGAAGCAGCAATCTTCAGGACCCCT
GGAGATGCAGAGGCATGGCAGTTCTCGGGATAACCCAGGCGGAGAAATGAAAATGAACAAGCAGCTATT
GTCGCCCTCCAGAGGTGCTTAGAATTACAGCCCAACAACCTAAAAGCTTTGATGGCCTTGCTGTGAGT
TATACTAACACTGGCCATCAGCAGGATGCCTGTGACGCTCTGAAGAATTGGATTAAGCAAAATCCAAAG
TACAAATACCTTGTGAAAAGCAAGAAGGGATCTCCAGGCCTCACCCGGCGGATGTCTAAGTCCCCAGTT
GATAGCTCTGTTCTGGAAGGGTGAAGGAATTATCTGGAAGCTGCCACCAAAATGGAGATATGATC
GACCCAGACCTGCAGACAGGTCTAGGGTTCTGTTCCACCTGAGTGGAGAATTTAATAGACCAATAGAT
GCATTTAACGCTGCCTTAACTGTTCCGCCAGAGGACTATTCATATGGAACCGCCTCGGGGCGACCTTG
GCGAACGGAGACCGCAGCGAGGAAGCCGTGGAGGCTATACGCGAGCACTGGAGATTCAGCCAGGATTC
ATCCGGTCCAGATACAACCTAGGAATAAGCTGCATCAACCTGGGCGCCTACAGAGAAGCGGTGAGCAAT
TTTCTCACTGCCCTCAGTTTCAAAGAAAAGAGCAGGAATCAGCAGCAAGTTCTCATCCTGCAATCTCT
GGGAATATCTGGGCTGCCTCAGAAATTGCGCTCTCTCTGATGGACCAACCAGAACTTCCAGGCGGCT
AATCTTGGTGACCTGGATGTCCTCTTAAGAGCTTCAACTGGATCCTTGA

```

Restriction Sites: SgfI-MluI



[View online »](#)

ACCN:	NM_001256752
Insert Size:	1776 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001256752.1
RefSeq Size:	9100 bp
RefSeq ORF:	1776 bp
Locus ID:	51555
UniProt ID:	Q8IYB4
Cytogenetics:	3q26.33
Protein Families:	Druggable Genome
MW:	65.9 kDa
Gene Summary:	<p>Accessory subunit of hyperpolarization-activated cyclic nucleotide-gated (HCN) channels, regulating their cell-surface expression and cyclic nucleotide dependence.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (4) encodes isoform 4. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>