

Product datasheet for SC332465

PEX5 related protein (PEX5L) (NM_001256754) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PEX5 related protein (PEX5L) (NM_001256754) Human Untagged Clone
Tag:	Tag Free
Symbol:	PEX5 related protein
Synonyms:	PEX5R; PEX5RP; PXR2; PXR2B; TRIP8b
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC332465 representing NM_001256754. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGTGATGAAGGAGATACCGAGGGAGGAGTCTGCTGAAGAAAAGCCCTCCTTACTATGACATCACAG
CTGGTGAATGAGCAACAAGAAAGCAGACCCCTCCTGAGTCCATCGATGACTTTCTGTGAAACC
AAATCGGAAGCAATAGCAAGGCCAGTAACATCCAATACAGCTGTATTGACCACTGGCTTAGATCTCCTC
GACCTGAGTGAACAGTCTCTCAAACCCAAACCAAGCCAAGAAGTCAGAGCCCTCATCAAAAACCTCA
TCCCTCAAGAAAAGGCCGATGGATCTGACCTCATCAGCACGGATGCTGAGCAGAGAGGCCAGCCTCTC
AGAGTCCCGGAGACTTTCCTTAGATCTAGACATTCAAACACAAGTGGAAAAATGGGACGATGTTAAG
TTTCATGGAGATCGAAATACCAAGGGACATCCAATGGCAGAGAGAAAATCATCCTCATCTAGAAGTGA
TCAAAAAGAGCTCTTATGGTCTCAGAACACAGATCTCAACCAGAAGTGGTGGAAAAAGCGCCCTC
AACTCTGAGTCGGCTTACAGAATTGGAATTAGTGGCTCCGACTCAGGCTCGACTGACCAAAGAATCGC
TGGGGAAGCGCATTACTTTCTAGAAACCACTCCTTAGAAGAAGAGTTTGAAGGGGCAAAGCAGCAGTG
GAGTCAGATACAGAGTTTTGGGATAAGATGCAAGCAGAATGGGAAGAAATGGCTCGGAGGAAGTGGATA
TCTGAGAACCAAGAAGCCGAGAACCAAGTAACCATCTCGGCTAGTGAGAAGGGATATTACTTTACACT
GAAAACCCCTTCAAGGACTGGCCTGGAGCATTGGAAGAAGGCTTAAAAAGGCTGAAGGAAGGGGATCTG
CCAGTACCATCCTGTTTCATGGAAGCAGCAATTCTTCAGGACCCGGAGATGCAGAGGCATGGCAGTTT
CTCGGGATAACCCAGGCGGAGAATGAAAATGAACAAGCAGCTATTGTCCGCCCTCCAGAGGTGCTTAGAA
TTACAGCCCAACAACTTAAAAAGCTTTGATGGCCTTGGCTGTGAGTTATACTAACACTGGCCATCAGCAG
GATGCCTGTGACGCTCTGAAGAATTGGATTAAGCAAAAATCCAAAGTACAAATACCTTGTGAAAAGCAAG
AAGGGATCTCCAGGCCTCACCGCGGATGTCTAAGTCCCGAGTTGATAGCTCTGTTCTGGAAGGGGTG
AAGGAATTATATCTGGAAGCTGCCACCAAAAATGGAGATATGATCGACCCAGACCTGCAGACAGGTCTA
GGGTTTCTGTTCCACCTGAGTGGAGAATTTAATAGAGCAATAGATGCATTTAACGCTGCCTTAACTGTT
CGGCCAGAGGACTATTTACTATGGAACCGCCTCGGGGCGACCTTGCGCAACGGAGACCGCAGCGAGGAA
GCCGTGGAGGCCTATACGCGAGCACTGGAGATTAGCCAGGATTCATCCGGTCCAGATACAACCTAGGA
ATAAGCTGCATCAACCTGGGCGCTACAGAGAAGCGGTGACGAATTTTCTCACTGCCTCAGTTTGCAA
AGAAAAGAGCAGGAATCAGCAGCAAGTTCCTCATCTGCAATCTCTGGGAATATCTGGGCTGCCCTCAGA
ATTGCGCTCTCTGATGGACCAACCAGAAGTCTTCCAGGCGGCTAATCTTGGTACCTGGATGTCTCT
TTAAGAGCTTTCAACTTGGATCTTGA

```

Restriction Sites: Sgfl-Mlul



[View online »](#)

ACCN:	NM_001256754
Insert Size:	1752 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001256754.1</u>
RefSeq Size:	8935 bp
RefSeq ORF:	1752 bp
Locus ID:	51555
UniProt ID:	<u>Q8IYB4</u>
Cytogenetics:	3q26.33
Protein Families:	Druggable Genome
MW:	65 kDa
Gene Summary:	<p>Accessory subunit of hyperpolarization-activated cyclic nucleotide-gated (HCN) channels, regulating their cell-surface expression and cyclic nucleotide dependence.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (6), as well as variant 18, encodes isoform 6. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>