

Product datasheet for **SC332483**

ZNF195 (NM_001256824) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: ZNF195 (NM_001256824) Human Untagged Clone
Tag: Tag Free
Symbol: ZNF195
Synonyms: HRF1; ZNFP104
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC332483 representing NM_001256824.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGCTGGTGCTCAGACTCTGTTGACGTT CAGGGATGTGGCCATAGAATTCTCCCTGGAGGAGTGGAAA
TGCCTGGACCTCGCTCAGCAGAAATTTGTACAGGGATGTGATGTTGGAGA ACTACAGAACTTGTCTCC
GTTGGTCTCACTGTCTGTAAGCCAGGCCTGATCACCTGCCTGGAGCAACGAAAAGAGCCCTGGAATGTG
AAGAGACAGGAGGCAGCAGACGGACATCCAGAGATGGGATTTACCATGCTACTCAGGCTTGTCTTGAA
CTCCTGGGCTCAAGCGATCTGCCTGCCTCAGCCTCCCAAAGTGTGGGATTACAGCTATGTCTTCTCAT
TTTACCCAAGACCTTCTGCCAGAGCAGGGCATAACAAGATGCATCCCAAAGAATACTGAGAGGATAT
GGAAATTTGTGGCCTTGATAATTTATATTTAAGGAAAGACTGGGAAAGTTTAGATGAGTGAAGTTGCAA
AAAGATTATAATGGACTTAACCAATGTTTCATCAACTACCCATAGCAAAATCTTTCAATATAATAAATAT
GTTAAATCTTTGATAACTTTTCAAATTTACATAGACGTAATATAAGTAATACTGGAGAGAAACCTTTC
AAATGTCAAGAAATGTGGCAAATCCTTTCAAATGCTCTCATTCTAACTGAACATCAGAAAATTCACACT
GGAAAAAATTCAAAAATGTGGAGAATGTGGCAAACCTTTATCCAGTGCTCACACTTTACTGAACCT
GAGAACATTGACTGGAGAGAAACCTTACAAGTGTCAAGAATGTAACAACGTCATTAACCTTGTCTCA
GTCCTTACTAAAAATAGAATTTATGCCGAGGGGAACATTACAGATGTGAAGAATTTGGCAAAGTATTT
AACCAGTGCTCCACCTTACTGAACATGAGCATGGTACTGAGGAAAAACCTGCAAAATGAAGAGTGC
AGCAGTGCTTTATCTCTTGTCTCAAGCCTTTCTAATCAACAGATGATTCTTGCTGGAGAGAAGCTCTCC
AAATGTGAAACATGGTACAAAGTTTTAAACCACAGCCAAATCCTTCCAACACCCAGAGAAATGAGATT
GGAGGGAAACCTTTCAAATGTGAGGAATGTGACAGCATCTTCAAGTGGTTCTCAGACCTTACTAAACAT
AAGAGAATTCACACTGGTGAGAAACCATACAAATGTGACGAATGTGGAAAGCCTATACACAGTCTCA
CACCTCAGTGAACACAGGAGGATTCACACCGGAGAGAAACCTACCAATGTGAAGAATGTGGGAAGGTC
TTCAGAACTTGTCAAGCCTTTCTAACCATAAGAGA ACTCATTCTGAAGAAAAACCTACACGTGTGAA
GAATGTGGCAACATCTTTAAGCAGTTATCAGACCTCACTAAGCATAAGAAAAACCCATACTGGAGAGAAG
CCCTACAATGTGACGAATGTGGAAAAA ACTTTACCCAGTCTCCAACCTTATTGTACATAAGAGAATT
CATACTGGAGAGAAACCTACAAGTGTGAAGAATGTGGCAGAGTCTTCATGTGGTTCTCAGACATTACC
AAACATAAGAAAACCCATACTGGAGAGAAACCTACAATGTGACGAATGTGGAAAAA ACTTTACCCAG
TCCTCAAACCTTATTGTACATAAGAGAATT CATACTGGAGAGAAACCTACAAGTGTGAAAAGTGTGGC
AAAGCCTTCAACCCAGTCTCACACCTGACTGTACATGAAAGCATT CATACTGA
  
```

Restriction Sites: SgfI-MluI



[View online »](#)

ACCN:	NM_001256824
Insert Size:	1779 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001256824.1
RefSeq Size:	3239 bp
RefSeq ORF:	1779 bp
Locus ID:	7748
UniProt ID:	O14628
Cytogenetics:	11p15.4
Protein Families:	Transcription Factors
MW:	68.1 kDa
Gene Summary:	<p>This gene encodes a protein belonging to the Krueppel C2H2-type zinc-finger protein family. These family members are transcription factors that are implicated in a variety of cellular processes. This gene is located near the centromeric border of chromosome 11p15.5, next to an imprinted domain that is associated with maternal-specific loss of heterozygosity in Wilms' tumors. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2012]</p> <p>Transcript Variant: This variant (11) differs in its 5' UTR, uses an alternate start codon, uses an alternate in-frame splice site and lacks an alternate in-frame exon in the 3' coding region, compared to variant 1. The encoded isoform (8) has a distinct N-terminus and is shorter than isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>