

Product datasheet for **SC332484**

ZNF195 (NM_001256825) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: ZNF195 (NM_001256825) Human Untagged Clone
Tag: Tag Free
Symbol: ZNF195
Synonyms: HRF1; ZNFP104
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC332484 representing NM_001256825.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGCTGGTGCTCAGACTCTGTTGACGTT CAGGGATGTGGCCATAGAATTCTCCCTGGAGGAGTGGAAA
TGCCTGGACCTCGCTCAGCAGAAATTTGTACAGGGATGTGATGTTGGAGA ACTACAGAACTTGTCTCC
GTTGGTCTCACTGTCTGTAAGCCAGGCCTGATCACCTGCCTGGAGCAACGAAAAGAGCCCTGGAATGTG
AAGAGACAGGAGGCAGCAGACGGACATCCAGCTATGTCTTCTCATTTTACCCAAGACCTTCTGCCAGAG
CAGGGCATAACAAGATGCATTCCAAAAAGAATACTGAGAGGATATGGAAATGTGGCCTTGATAATTTA
TATTTAAGGAAAGACTGGGAAAGTTTAGATGAGTGAAGTTGCAAAAAGATTATAATGGACTTAACCAA
TGTTTCATCAACTACCCATAGCAAAATCTTTCAATAATAATAATATGTTAAAATCTTTGATAACTTTTCA
AATTTACATAGACGTAATAATAAGTAATACTGGAGAGAAACCTTTCAAATGTCAAGAATGTGGCAAATCC
TTTCAAATGCTCTCATTCCCTAACTGAACATCAGAAAAATTCACACTGGAAAAAAATTCAAAAATGTGGA
GAATGTGGCAAACCTTTATCCAGTGCTCACACTTTACTGAACTGAGAACATTGACACTGGAGAGAAA
CCTTACAAGTGTCAAGAATGTAACAACGTCATTA AAACTTGCTCAGTCCTTACTAAAAATAGAATTTAT
GCCGGAGGGGAACATTACAGATGTGAAGAATTTGGCAAAGTATTTAACCAAGTGCTCCACCTTACTGAA
CATGAGCATGGTACTGAGGAAAAACCCTGCAAAATGAAGAGTGCAGCAGTGTCTTTATCTCTTGCTCA
AGCCTTTCTAATCAACAGATGATTTCTTGCTGGAGAGAAGCTCTCAAATGTGAAACATGGTACAAAGGT
TTTAACACAGCCAAAATCCTTCCAAACACCAGAGAAATGAGATTGGAGGGAAAACCTTTCAAATGTGAG
GAATGTGACAGCATCTTCAAGTGGTTCTCAGACCTTACTAAACATAAGAGAATTCACACTGGTGAGAAA
CCATACAAATGTGACGAATGTGGGAAAGCCTATACACAGTCCTCACACCTCAGTGAACACAGGAGGATT
CACACCGGAGAGAAAACCTACCAATGTGAAGAATGTGGGAAGGTCTTCAGAACTTGCTCAAGCCTTTCT
AACCATAAAGAGAACTCATTCTGAAGAAAAACCTACACGTGTGAAGAATGTGGCAACATCTTTAAGCAG
TTATCAGACCTCACTAAGCATAAGAAAACCCATACTGGAGAGAAGCCCTACAAATGTGACGAATGTGGA
AAAACTTTACCCAGTCCTCAACCTTATTGTACATAAGAGAATTCATACTGGAGAGAAAACCTACAAG
TGTGAAGAATGTGGCAGAGTCTTCATGTGGTTCTCAGACATTACCAAACATAAGAAAACCCATACTGGA
GAGAAAACCTACAAATGTGACGAATGTGGAAAAA ACTTTACCCAGTCCTCAAACCTTATTGTACATAAG
AGAATTCATACTGGAGAGAAAACCTACAAGTGTGAAAAGTGTGGCAAAGCCTTACCCAGTTCTCACAC
CTGACTGTACATGAAAGCATTCACTTGA
  
```

Restriction Sites: Sgfl-Mlul
ACCN: NM_001256825



[View online »](#)

Insert Size:	1686 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001256825.1</u>
RefSeq Size:	3165 bp
RefSeq ORF:	1686 bp
Locus ID:	7748
UniProt ID:	<u>O14628</u>
Cytogenetics:	11p15.4
Protein Families:	Transcription Factors
MW:	65 kDa
Gene Summary:	<p>This gene encodes a protein belonging to the Krueppel C2H2-type zinc-finger protein family. These family members are transcription factors that are implicated in a variety of cellular processes. This gene is located near the centromeric border of chromosome 11p15.5, next to an imprinted domain that is associated with maternal-specific loss of heterozygosity in Wilms' tumors. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2012]</p> <p>Transcript Variant: This variant (12) differs in its 5' UTR, uses an alternate start codon, and lacks two alternate in-frame exons in the 3' coding region, compared to variant 1. The encoded isoform (6) has a distinct N-terminus and is shorter than isoform 1. Both variants 6 and 12 encode isoform 6. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>