

Product datasheet for **SC332491**

DNAJC6 (NM_001256864) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DNAJC6 (NM_001256864) Human Untagged Clone
Tag:	Tag Free
Symbol:	DNAJC6
Synonyms:	DJC6; PARK19
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

Fully Sequenced ORF: >SC332491 representing NM_001256864.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGAGCCTCTCGGGAGCTACCGGAAAAAGACCAGCAACGATGGTTATGAATCTTTCAGCTGGTGGAC
AGTAACGGGGACTTAAGTGCGGGAAGCGCGGGGTTGGCGCAAGCAGAGAGTGAACGCCGGGGCAGCG
GCGGGAGTCCC GCCCGACAGCCTCCGGACCGCCAGCACCATGGACAGCTCAGGTGCCTCATCTCCA
GACATGGAGCCCAGCTATGGGGAGGTCTCTTTGACATGGTAAAAGGAGGTGCAGGGAGGCTCTTTAGT
AACCTAAAGGACAACCTTGAAGACACCCTCAAAGACACATCTTCTAGAGTGATACAATCTGTGACCAGC
TACACAAAAGGAGATTTAGACTTCACTTATGTTACCTCCAGAATTATTGTGATGTCCTTTCTCTGGAC
AATGTTGACATAGGATTCAGGAATCAGGTTGATGACATTCGAAGCTTTTTGGATTCCAGACATCTTGAC
CACTACACAGTATACAATCTGTCACCTAAGTCTTATCGAACTGCCAAGTTTCACAGCCGGGTCTCAGAA
TGCAGTTGGCCATTAGGCAGGCTCCAGTCTGCACAACCTTTTTGCTGTGTGTCGGAATATGTATAAC
TGGCTACTGCAGAATCCAAAAATGTCTGTGTTGCCACTGCTTGGATGGACGGGCGGCATCATCAATT
CTGTTGGTGCTATGTTCATTTTCTGTAATCTCTACTCTACTCTGCCCAGCCATTGATTGCTATAT
GCAAAGCGACCAGGAATTGGACTTTCACCATCCATAGGAGATACCTGGGCTATATGTGTGACCTACTG
GCAGACAAGCCCTACCGCCTCACTTCAAGCCTCTCACAATTAAGTCGATCACTGTCAAGTCAATACCC
TTTTTCAACAAAACAGAGGAATGGATGTCGCCCTTACTGTGATGACTCATTGGAGAAACAAAATATAT
TCGACTTGCACAGATTTTGAACGAATGAAAGAATATCGTGTCCAAGATGGAAAAATCTTCATTCCCTTG
AACACTACTGTGCAAGGAGACGTGGTGTTCATGATCACTTGAGGTCAACCATTGGGAGCCGGCTA
CAGGCTAAGGTGACCAACACACAGATATCCAGCTTCACTTTCACACTGGATTCACTACTGGACACA
ACAGTTTTAAAGTTCACCAAGCCTGAGTTAGATGCATGTGATGTACCAGAAAAATATCCTCAGCTATTT
CAGGTGACACTGGATGTAGAACTACAGCCCATGACAAAAGTAAGACTTAACTCCACCATGGGAACAT
TACTGCACAAAAGATGTCAATCCAGCATCTCTCTCTCTCACCAGGAACATCAAGATACGCTGGCC
TTAGGAGGACAGGCTCCAATAGATATCCCTCAGACAACCCAGGCATTACGGACAAAAGTGGTTTTCTTT
GCCTCTCTCTGTTGGCAAGATCAGAAATCGGAGAAGTCATTCTGTGAGGAGGACCACGCTGCCCTAGTG
AATCAGGAAAGTGAGCAATCAGATGATGAACCTTCTGACACTTTCAGTCCGCATGGCAATGCCAATGGT
GACAAGCCTCATGGAGTCAAGAAGCCAGCAAAAAGCAGCAGGAGCCAGCAGCCCTCCACCCCTGAG
GATGTGGACCTTTTGGGCTGGAAGGTCTGCAATGAGTAACAGCTTCTCTCCGCCAGCGGCTCTCCC
ACCAATTCTGAACACTGAGTGACCTGTTGGGGTGGAGGTGCAGCTGGTCCCACCCAGGCTGGACAG
TCAGGAGTGGAAGATGTGTTTCATCCTAGTGGACCTGCGTCTACCCAGTCAACACCACGCCGCTCTGCC
ACCTCCACCTCTGCGTCTCCAACCCTAAGAGTGGGAGAAGTGCCACCTTGACCATTTGGAGCACCT
TCTAAACCATCAGGTGAGGATTTGCTGGTCTTTTCTGAACACATCCAGTGCTTCCAGTGACCCCTTT
CTCCAGCCAACAAGAAGTCTTCGCCACAGTACATGCTTCTAGTACGCCTGCTGTGAACATTCAGCCA
GATGTTTCTGGAGTTGGGACTGGCATGCTAAACCAGGAGGCTTTGGAATGGGAAGCAAGTCAGCTGCC
ACCAGCCAACCGGATCCTCGCATGGTACTCCCACCCATCAAAGCAAACCCAGACTCTGGATCCTTTT
GCCGACCTTGGGACACTAGGTAGTCTTCTCTTGGCAGCAAACCCACCACTGGATTGGGTGGA
GGATCCCCTCTCAGCTCGCCACAGAAGGCGTCTCCCAGCCTATGGGTGGCGGGTGGCAGCAGGGA
GGTGCCTACAACGTGAGCTTCTCAGCCATGCCTGGGGCCAGAACGAACGTGGGAAAGGATCAAGT
AATTTGGAAGGGAAACAAAAGCAGCTGATTTTGAAGACCTACTCTCTGGTCAAGTTTTCAATGCTCAC
AAAGACAAAAGGGCCCTCGGACAATAGCTGAGATGAGAAAGGAGGAAATGGCCAAGGAAATGGATCCT
GAGAAATTAAGATTCTGGAATGGATTGAAGGCAAAAGAAAGAAATATCAGAGCCCTTCTTCCACGATG
CATACCGTACTATGGGCTGGGAGACCAAGTGGAAACCAGTTGGCATGGCAGACCTGGTAACACCAGAG
CAGGTGAAGAAGGTGTACAGGAAGGCTGTCTGGTGGTGCACCCAGATAAAGCTACTGGCAACCCAT
GAACAATACGCAAAGATGATTTTCATGGAGCTCAATGATGCTGGTCTGAATTTGAAAACCAAGGCCAA
AAGCCCTATATAA
  
```

Restriction Sites: Sgfl-Mlul
ACCN: NM_001256864
Insert Size: 2913 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001256864.1
RefSeq Size:	5968 bp
RefSeq ORF:	2913 bp
Locus ID:	9829
UniProt ID:	O75061
Cytogenetics:	1p31.3
Protein Families:	Druggable Genome
Protein Pathways:	Endocytosis
MW:	105.7 kDa
Gene Summary:	<p>DNAJC6 belongs to the evolutionarily conserved DNAJ/HSP40 family of proteins, which regulate molecular chaperone activity by stimulating ATPase activity. DNAJ proteins may have up to 3 distinct domains: a conserved 70-amino acid J domain, usually at the N terminus, a glycine/phenylalanine (G/F)-rich region, and a cysteine-rich domain containing 4 motifs resembling a zinc finger domain (Ohtsuka and Hata, 2000 [PubMed 11147971]).[supplied by OMIM, Mar 2008]</p> <p>Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>