

Product datasheet for **SC332506**

DDX11 (NM_001257144) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DDX11 (NM_001257144) Human Untagged Clone
Tag:	Tag Free
Symbol:	DDX11
Synonyms:	CHL1; CHLR1; KRG2; WABS
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

Fully Sequenced ORF: >SC332506 representing NM_001257144.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGCTAATGAAACACAGAAGGTTGGTGCCATCCATTTTCCTTTCCCTTCACACCCTATTCCATCCAG
GAAGACTTCATGGCAGAGCTGTACCGGTTTTGGAGGCTGGCAAGATTGGGATATTTGAGAGTCCAAT
GGCACTGGGAAGTCCCTTAAGTCTATTTGTGGGGCCCTCTTTGGCTCCGTGACTTTGAACAGAAGAAG
CGTGAAGAAGAGGCACGACTCCTTGAACCTGGAACGGCCCTTACATGATGAGAAGATGAATCCCTG
TGTCTGTCTTCTTCTGCGAAGGGGCTGCAGGCACCCGAGGCCTGCTGGAGAACCAGGCCTGGGTTACT
CAGTTTGTGCAGAAGAAAGAAGAGAGGGACCTGGTGGACCGACTAAAGCGGAGCAGGCCAGGAGGAAG
CAGCGAGAAGAACGCCTGCAGCAGCTGCAGCACAGGGTGCAGCTCAAGTATGCAGCCAAGCGCTGAGG
CAGGAAGAAGAAGAAAGAGAGAATCTCTCCGCTCAGCAGGGAGATGCTAGAGACAGGCCCGGAGGCT
GAGCGGCTGGAGCAGCTGGAGTCTGGGGAGGAGGAGTGGTCTCGCGGAATACGAGAGTATGAGGAG
AAAAAGTGGCGAGCAGAGTGGATGAGGATGAGGATGACCTGGAGGAAGAACACATAACTAAGATTTAT
TACTGTAGTCGGACACTCCCAGCTGGCCAGTTTGTGCATGAGGTGAAGAAGAGCCCTTTGGCAAG
GATGTTCCGGCTGGTCTCCCTTGGCTCCCGCAGAACCTTTGTGTAATGAAGACGTGAAAAGCCTAGGT
TCTGTGCAGCTTATCAACGACCGCTGTGTGGACATGCAGAGAAGCAGGCACGAGAAGAAGAAAGGAGCT
GAGGAGGAGAAGCCAAAGAGGAGGAGGCAGGAGAAGCAGGCAGCCTGCCCTTCTACAACCAGCAGCAG
ATGGGCCTTCTCCGGATGAGGCCCTGGCAGAGGTGAAGGACATGGAGCAGCTGCTGGCCCTTGGGAAG
GAGGCCCGGGCTGTCCCTATTACGGGAGCCGCTTGGCATCCCTGCAGCCAGCTGGTGGTGTGCC
TATCAGATGCTGCTGCATGCGGCCACTCGGCAGCCGCGGGCATCCGGCTGCAGGACCAGGTGGTATC
ATCGACGAGGCGCACAACTGATCGACACCATCACGGGCATGCACAGCGTGGAGGTCAGCGGCTCCCAG
CTCTGCCAGGCCATCCCAGCTGCTGCAGTACGTGGAGCGATACGGGAAGCGTTTGAAGCCAAAGAAC
CTGATGTACCTGAAGCAGATCCTGTATTTGTGGAGAAATTCGTGGCTGTCTAGGGGGAAACATTAAG
CAAAATCCAATACACAGAGTCTGTACAGACAGGGACGGAGCTGAAGACCATCAACGACTTCTCTTC
CAGAGCCAGATCGACAACATCAACCTGTTCAAGGTGCAGCGATACTGTGAGAAGAGCATGATCAGCAGA
AAGCTCTTTGGATTCACTGAACGGTACGGAGCAGTGTCTCATCCGGGAGCAGCCAAACTGGCTGGG
TTTCAGCAATTCCTGCAGAGCTGCAGCCAGGACGACTGAAGCTTTCAGCCCTGCAGACGAGAGT
CAGGCCAGCACCTGCGACCAGTCTCCACTGATGCACATCCAAGGCTTCTGGCAGCTCTCACTACG
GCCAACCCAGGACGGCAGGGTCACTCTGAGCCGCAAGGCAGCCTCAGTCAGAGCACCTGAAGTTTTTG
CTCCTGAATCCAGCTGTGCACTTTGCCAAGTGGTGAAGGAATGCCGGGAGTGGTATTGCGGGGGT
ACCATGCAGCCGGTGTCTGACTTCGGCAGCAGCTGCTGGCCTGTGCCGGGTGGAAGCTGAGCGCGTG
GTGGAGTTTTCTGTGGTCACGTGATCCCTCCAGACAACATCCTGCCCTCGTCACTGCAGCGGATC
TCCAACCAGCCGCTGGAATTCACGTTCCAGAAAAGAGAGCTGCCTCAGATGATGGACGAGGTGGTGC
ATTCTCTGTAACCTGTGCGGTGTGGTCTCTGGAGGGTGGTCTGTTTTCTCCCTCTACAGTACCTG
CGCCAGGTCCATGCCACTGGGAGAAGGGTGGCCTGCTGGGCCGTCTGGCTGCCAGGAAGAAGATATTC
CAGGAACCTAAGAGCGCACACCAGGTGGAGCAGGTGCTGCTGGCATATCCAGGTGCATCCAGGCCCTGT
GGCCAGGAGAGAGGCCAGGTGACAGGGGCCCTGCTCCTCTGTGGTTGGAGGAAGATGAGTGAAGGG
ATCAACTTCTCTGACAACCTAGGCCGGTGTGTGGTGTGGTGGCATGCCCTTCCCAACATCAGGTCT
GCAGAGCTGCAGGAGAAGATGGCCTACTTGGATCAAACCCTCAGCCCCGCCAGGCACCCCAAGGAA
GGCTCTGGTGGAGAACCCTGTGCATGAAGCCGTCACCAGTCCATAGGCAGGGCCATCAGGCACCCAGAA
GGATTTTGCCAGCGTAGTCTCCTGGACCAGCATATGCCGGCCCCCTGTCTGGCCAAGCTGCCGGC
CTGGATCCGAGCCGTGTGGAGGTCAAAGCTACCTTTGGCCCCGCAATTGCTGCTGTGCAGAAGTTTCA
CCGGGAGAAGTCGGCCTTCTCTGATGGGCAACCACCACTGCCTGGCGCCGTGCCCTTCTTTGTCC
TGCCCCGTGGAGACAGTGTTCGTGGCGTGGTCTGCGGGGATCCTGTTACAAAGGTGAAACCCAGG
AGGAGAGTGTGGAGTCCAGAGTGTGCCAGGACCCAGGCACAGGCCTTAGCTCCCGTAGGAGAAAATGG
GGGAATCCTGAATGA
  
```

Restriction Sites: SgfI-MluI
ACCN: NM_001257144
Insert Size: 2913 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001257144.1
RefSeq Size:	3823 bp
RefSeq ORF:	2913 bp
Locus ID:	1663
UniProt ID:	Q96FC9
Cytogenetics:	12p11.21
Protein Families:	Stem cell - Pluripotency
MW:	108.3 kDa
Gene Summary:	<p>DEAD box proteins, characterized by the conserved motif Asp-Glu-Ala-Asp (DEAD), are putative RNA helicases. They are implicated in a number of cellular processes involving alteration of RNA secondary structure such as translation initiation, nuclear and mitochondrial splicing, and ribosome and spliceosome assembly. Based on their distribution patterns, some members of this family are believed to be involved in embryogenesis, spermatogenesis, and cellular growth and division. This gene encodes a DEAD box protein, which is an enzyme that possesses both ATPase and DNA helicase activities. This gene is a homolog of the yeast CHL1 gene, and may function to maintain chromosome transmission fidelity and genome stability. Alternative splicing results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (4) differs in the 5' UTR compared to variant 3. Variants 3 and 4 encode the same protein. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>