

Product datasheet for **SC332542**

IL4R (NM_001257406) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	IL4R (NM_001257406) Human Untagged Clone
Tag:	Tag Free
Symbol:	IL4R
Synonyms:	CD124; IL-4RA; IL4RA
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

Fully Sequenced ORF: >SC332542 representing NM_001257406.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGGGTGGCTTTGCTCTGGGCTCCTGTTCCCTGTGAGCTGCCTGGTCTGCTGCAGGTGGCAAGCTCT
GGGAACATGAAGGTCTTGCAGGAGCCCACCTGCGTCTCCGACTACATGAGCATCTCTACTTGGCAGTGG
AAGATGAATGGTCCCACCAATTGCAGCACCGAGCTCCGCTGTTGTACCAGCTGGTTTTCTGCTCTCC
GAAGCCCACACGTGTATCCCTGAGAACAACGGAGGCGGGGTGCGTGTGCCACCTGCTCATGGATGAC
GTGGTCAGTGCGGATAACTATACACTGGACCTGTGGGCTGGGCAGCAGCTGCTGTGGAAGGGCTCCTTC
AAGCCCAGCGAGCATGTGAAACCCAGGGCCCCAGGAAACCTGACAGTTACACCAATGTCTCCGACACT
CTGCTGCTGACCTGGAGCAACCCGTATCCCCCTGACAATTACCTGTATAATCATCTCACCTATGCAGTC
AACATTTGGAGTGAAAACGACCCGGCAGATTTCAAGATCTATAACGTGACCTACCTAGAACCCTCCCTC
CGCATCGCAGCCAGCACCTGAAGTCTGGGATTTCTACAGGGCACGGGTGAGGGCTGGGCTCAGTGC
TATAACACCACCTGGAGTGAGTGGAGCCCCAGCACCAAGTGGCACAACCTCTACAGGGAGCCCTCGAG
CAGCACCTCCTGTGGGCTCAGCGTTTCTGCATTGTCATCCTGGCCGTCTGCCTGTTGTGCTATGTC
AGCATCACAAGATTAAGAAAGAATGGTGGGATCAGATTCACAACCCAGCCCGCAGCCGCTCGTGGCT
ATAATAATCCAGGATGCTCAGGGGTACAGTGGGAGAAGCGGTCCCGAGGCCAGGAACCAGCCAAGTGC
CCACACTGGAAGAATTGTCTTACCAAGCTCTTGCCTGTTTTCTGGAGCACAACATGAAAAGGGATGAA
GATCCTCACAAGGCTGCCAAAGAGATGCCTTCCAGGGCTCTGGAAAATCAGCATGGTGGCCAGTGGAG
ATCAGCAAAGACAGTCTCTGGCCAGAGAGCATCAGCGTGGTGCGATGTGTGGAGTTGTTTGGGCCCCG
GTGGAGTGTGAGGAGGAGGAGGAGGTAGAGGAAGAAAAGGGAGCTTCTGTGCATCGCCTGAGAGCAGC
AGGGATGACTTCCAGGAGGAAGGGAGGGCATTGTGGCCCGGCTAACAGAGAGCCTGTTCTGGACCTG
CTCGGAGAGGAGAATGGGGCTTTTCCAGCAGGACATGGGGAGTTCATGCCTTCTTCCACCTTCGGGA
AGTACGAGTGTACATGCCCTGGGATGAGTCCCAAGTGCAGGGCCCAAGGAGGCACCTCCCTGGGGC
AAGGAGCAGCCTCTCCACCTGGAGCCAAGTCTCCTGCCAGCCGACCCAGAGTCCAGACAACCTGACT
TGCACAGAGACGCCCTCGTTCATCGCAGGCAACCTGCTTACCGCAGCTTACGAACTCCCTGAGCCAG
TCACCGTGTCCCAGAGAGCTGGGTCCAGACCCACTGCTGGCCAGACACCTGGAGGAAGTAGAACCCGAG
ATGCCCTGTGTCCCCAGCTCTCTGAGCCAACCACTGTGCCCAACCTGAGCCAGAAACCTGGGAGCAG
ATCCTCCGCCGAAATGTCTCCAGCATGGGGCAGCTGCAGCCCCGCTCTCGGCCCCACCAGTGGCTAT
CAGGAGTTGTACATGCGGTGGAGCAGGTGGCACCCAGGCCAGTGGGTGGTGGGCTGGGTCCCCCA
GGAGAGGTGGTTACAAGCCTTCTCAAGCCTGCTTGCCAGCAGTGTGTGTCCCAGAGAAATGTGGG
TTTGGGGCTAGCAGTGGGAAGAGGGGTATAAGCCTTTCAAGACCTCATTCTGGCTGCCCTGGGGAC
CCTGCCCCAGTCCCTGTCCCCTTGTTCACCTTTGGACTGGACAGGGAGCCACCTCGCAGTCCGCAGAGC
TCACATCTCCCAAGCAGCTCCCCAGAGCACCTGGGTCTGGAGCCGGGGGAAAAGGTAGAGGACATGCCA
AAGCCCCCACTTCCCCAGGAGCAGGCCACAGACCCCTTGTGGACAGCCTGGGCAGTGGCATTGTCTAC
TCAGCCCTTACCTGCCACCTGTGGCGCCACCTGAAACAGTGTGATGGCCAGGAGGATGGTGGCCAGACC
CCTGTGATGGCCAGTCTTGTGTGGTGTGCTGTGGAGACAGGTCTCGCCCCCTACAACCCCTG
AGGGCCCCAGACCCCTCTCCAGTGGGGTTCCTACTGGAGGCCAGTCTGTGTCCGGCCTCCCTGGCACCC
TCGGGCATCTCAGAGAAGAGTAAATCCTCATCATCCTTCCATCCTGCCCTGGCAATGCTCAGAGCTCA
AGCCAGACCCCAAAATCGTGAACCTTGTCTCCGTGGGACCCACATACATGAGGGTCTCTAG
  
```

Restriction Sites: Sgfl-Mlul

ACCN: NM_001257406

Insert Size: 2478 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001257406.1
RefSeq Size:	3577 bp
RefSeq ORF:	2478 bp
Locus ID:	3566
UniProt ID:	P24394
Cytogenetics:	16p12.1
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	Cytokine-cytokine receptor interaction, Hematopoietic cell lineage, Jak-STAT signaling pathway
MW:	89.7 kDa
Gene Summary:	<p>This gene encodes the alpha chain of the interleukin-4 receptor, a type I transmembrane protein that can bind interleukin 4 and interleukin 13 to regulate IgE production. The encoded protein also can bind interleukin 4 to promote differentiation of Th2 cells. A soluble form of the encoded protein can be produced by proteolysis of the membrane-bound protein, and this soluble form can inhibit IL4-mediated cell proliferation and IL5 upregulation by T-cells. Allelic variations in this gene have been associated with atopy, a condition that can manifest itself as allergic rhinitis, sinusitis, asthma, or eczema. Polymorphisms in this gene are also associated with resistance to human immunodeficiency virus type-1 infection. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Apr 2012]</p> <p>Transcript Variant: This variant (3) differs in the 5' UTR compared to variant 1. Both variants 1 and 3 encode the same isoform (a).</p>