

Product datasheet for **SC332543**

IL4R (NM_001257407) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	IL4R (NM_001257407) Human Untagged Clone
Tag:	Tag Free
Symbol:	IL4R
Synonyms:	CD124; IL-4RA; IL4RA
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

Fully Sequenced ORF: >SC332543 representing NM_001257407.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGCAGAAGGATGCCAGAAGAGAAGGGAACATGAAGGTCTTGCAGGAGCCCACCTGCGTCTCCGACTAC
ATGAGCATCTCTACTTGCAGTGAAGATGAATGGTCCCACCAATTGCAGCACCGAGCTCCGCCTGTTG
TACCAGCTGGTTTTCTGCTCTCCGAAGCCACACGTGTATCCCTGAGAACAACGGAGGCGCGGGTGC
GTGTGCCACCTGCTCATGGATGACGTGGTCACTGCGGATAACTATAACACTGGACCTGTGGGCTGGGCAG
CAGCTGCTGTGGAAGGGCTCCTTCAAGCCAGCGAGCATGTGAAACCCAGGGCCCCAGGAAACCTGACA
GTTACACCAATGTCTCCGACACTCTGCTGCTGACCTGGAGCAACCCGTATCCCCCTGACAAATTACCTG
TATAATCATCTCACCTATGCAGTCAACATTTGGAGTGAAAACGACCCGGCAGATTCAGAATCTATAAC
GTGACCTACCTAGAACCCTCCCTCCGCATCGCAGCCAGCACCCCTGAAGTCTGGGATTTCTACAGGGCA
CGGGTGAAGGGCTGGGCTCAGTGTATAACACCACCTGGAGTGAGTGGAGCCCCAGCACCAAGTGGCAC
AACTCCTACAGGGAGCCCTTCGAGCAGCACCTCCTGCTGGGCGTCAGCGTTTCTGCATTGTCATCTG
GCCGTCTGCCTGTTGTGCTATGTCAGCATACCAAGATTAAGAAAGAATGGTGGGATCAGATTCACAAC
CCAGCCCCGAGCCGCTCGTGGCTATAATAATCCAGGATGCTCAGGGTACAGTGGGAGAAGCGGTCC
CGAGGCCAGGAACAGCCAAGTCCCACACTGGAAGAATTGTCTTACCAAGCTCTTGCCCTGTTTTCTG
GAGCACAACATGAAAAGGGATGAAGATCCTCACAAGGCTGCCAAAGAGATGCCTTTCCAGGGCTCTGGA
AAATCAGCATGGTGCCAGTGGAGATCAGCAAGACAGTCTCTGGCCAGAGAGCATCAGCGTGGTGCGA
TGTGTGGAGTTGTTTGAAGCCCCGGTGGAGTGTGAGGAGGAGGAGGAGTGGAGGAGGAGGAGGAGGAGG
TTCTGTGCATCGCTGAGAGCAGCAGGGATGACTTCCAGGAGGGAAGGGAGGGCATTGTGGCCCCGGTA
ACAGAGAGCCTGTTCTGGACCTGCTCGGAGAGGAGAAATGGGGCTTTTGCCAGCAGGACATGGGGGAG
TCATGCCTTTCCACCTTCGGGAAGTACGAGTGTCCATGCCCTGGGATGAGTTCCTCAAGTGCAGGG
CCCAAGGAGGCACCTCCCTGGGCAAGGAGCAGCCTCTCCACCTGGAGCCAAGTCTCCTGCCAGCCCCG
ACCCAGAGTCCAGACAACCTGACTTGCACAGAGACGCCCCCTCGTATCGCAGGCAACCTGCTTACCCG
AGCTTCAGCAACTCCCTGAGCCAGTCCCGTGTCCAGAGAGCTGGGTCCAGACCCACTGCTGGCCAGA
CACCTGGAGGAAGTAGAACCCGAGATGCCTGTGTCCCCAGCTCTCTGAGCCAACCACTGTGCCCAA
CCTGAGCCAGAAACCTGGGAGCAGATCCTCCGCCGAAATGTCCTCCAGCATGGGGCAGCTGCAGCCCC
GTCTCGGCCCCCACCAGTGGCTATCAGGAGTTGTACATGCGGTGGAGCAGGGTGGCACCCAGGCCAGT
GCGGTGGTGGGCTGGGTCCCCAGGAGAGGCTGGTTACAAGGCCTTCTCAAGCCTGCTTCCAGCAGT
GCTGTGTCCCCAGAGAAATGTGGGTTGGGGCTAGCAGTGGGGAAGAGGGGTATAAGCCTTTCCAAGAC
CTATTCTGGTGCCTGGGGACCCTGCCCCAGTCCCTGTCCCCTTGTTCACCTTTGGACTGGACAGG
GAGCCACCTCGCAGTCCGAGAGCTCACATCTCCAAGCAGCTCCCAGAGCACCTGGGTCTGGAGCCG
GGGAAAAGGTAGAGGACATGCCAAAGCCCCACTTCCCAGGAGCAGGCCACAGACCCCTTGTGGAC
AGCCTGGGCGAGTGGCATTGTCTACTCAGCCCTTACCTGCCACCTGTGCGGCCACCTGAAACAGTGTGAT
GGCCAGGAGGATGGTGGCCAGACCCCTGTGATGGCCAGTCTTGTGTGGCTGCTGCTGTGGAGACAGG
TCCTCGCCCCCTACAACCCCTGAGGGCCCCAGACCCCTCTCCAGGTGGGGTCCACTGGAGGCCAGT
CTGTGTCCGGCCTCCCTGGCACCTCGGGCATCTCAGAGAAGAGTAAATCCTCATCATCTTCCATCCT
GCCCTGGCAATGCTCAGAGCTCAAGCCAGACCCCAAAATCGTGAACCTTGTCTCCGTGGGACCCACA
TACATGAGGGTCTCTAG
  
```

Restriction Sites: Sgfl-Mlul

ACCN: NM_001257407

Insert Size: 2433 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001257407.1
RefSeq Size:	3821 bp
RefSeq ORF:	2433 bp
Locus ID:	3566
UniProt ID:	P24394
Cytogenetics:	16p12.1
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	Cytokine-cytokine receptor interaction, Hematopoietic cell lineage, Jak-STAT signaling pathway
MW:	88.3 kDa
Gene Summary:	<p>This gene encodes the alpha chain of the interleukin-4 receptor, a type I transmembrane protein that can bind interleukin 4 and interleukin 13 to regulate IgE production. The encoded protein also can bind interleukin 4 to promote differentiation of Th2 cells. A soluble form of the encoded protein can be produced by proteolysis of the membrane-bound protein, and this soluble form can inhibit IL4-mediated cell proliferation and IL5 upregulation by T-cells. Allelic variations in this gene have been associated with atopy, a condition that can manifest itself as allergic rhinitis, sinusitis, asthma, or eczema. Polymorphisms in this gene are also associated with resistance to human immunodeficiency virus type-1 infection. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Apr 2012]</p> <p>Transcript Variant: This variant (4) differs in the 5' UTR and 5' coding region, and initiates translation at an alternate start codon, compared to variant 1. The encoded isoform (c) has a distinct N-terminus and is shorter than isoform 1.</p>