

Product datasheet for **SC332572**

IL4I1 (NM_001258017) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: IL4I1 (NM_001258017) Human Untagged Clone
Tag: Tag Free
Symbol: IL4I1
Synonyms: FIG1; hIL4I1; LAAO; LAO
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC332572 representing NM_001258017.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGCCCAACGATGACTTCTGTCCTGGGCTAACCATAAAGGCCATGGGTGCTGAGAGAGCCCCCAGAGG
CAGCCATGCACCCTGCACCTCCTCGTCCTCGTCCCATCCTCCTCAGCCTGGTGGCTCCCAGGACTGG
AAGGCTGAACGCAGCCAAGACCCCTTCGAGAAATGCATGCAGGATCCTGACTATGAGCAGCTGCTCAAG
GTGGTGACCTGGGGGCTCAATCGGACCTGAAGCCCCAGAGGGTATTGTGGTTGGCGTGGTGTGGCC
GGGCTGGTGGCCCAAGGTGCTCAGCGATGCTGGACACAAGGTACCATCCTGGAGGCAGATAACAGG
ATCGGGGGCCGCATCTTACCTACCGGGACCAGAACACGGGCTGGATTGGGGAGCTGGGAGCCATGCGC
ATGCCCAGCTCTACAGGATCCTCCACAAGCTCTGCCAGGGCCTGGGGCTCAACCTGACCAAGTTACC
CAGTACGACAAGAACACGTGGACGGAGGTGCACGAAGTGAAGCTGCGCAACTATGTGGTGGAGAAGGTG
CCCAGAAAGTGGGCTACGCCTTGCCTCCCAGGAAAAGGGCCACTCGCCGAAGACATCTACCAGATG
GCTCTCAACCAGGCCCTCAAAGACCTCAAGGCACTGGGCTGCAGAAAGGCGATGAAGAAGTTTAAAGG
CACACGCTCTTGAATATCTTCTCGGGGAGGGAACTGAGCCGGCCGGCCGTGCAGCTTCTGGGAGAC
GTGATGTCGAGGATGGCTTCTCTATCTCAGCTTCGCCGAGGCCCTCCGGGCCACAGCTGCCTCAGC
GACAGACTCCAGTACAGCCGATCGTGGTGGCTGGGACCTGCTGCCGCGCGCTGCTGAGCTCGCTG
TCCGGGTTGTGCTGTTGAACGCGCCCGTGGTGGCGATGACCCAGGGACCGCACGATGTGCACGTGCAG
ATCGAGACCTCTCCCCGCGCGGAATCTGAAGGTGCTGAAGGCCGACGTGGTGTGCTGACGCGGAGC
GGACCGCGGTGAAGCGCATCACCTTCTCGCCGCGCTGCCCGCCACATGCAGGAGGCGTGCAGGAGG
CTGCACTACGTGCCGGCCACCAAGGTGTTCTTAAGCTTCCGACGGCCCTTCTGGCGGAGGACACATT
GAAGGCGGCCACTCAAACACCGATCGCCCGTCCGCGATGATTTTCTACCCGCGCCGCGGAGGGCGCG
CTGCTGCTGGCCTCGTACACGTGGTCCGACGCGGCGGACGCTTCCGCGGCTTGAAGCGGGAAGAGGCG
TTGCGCTTGGCGCTCGACGACGTGGCGGATTGCACGGGCTGCTGCGCCAGCTCTGGGACGGCACC
GGCGTCTCAAGCGTTGGGCGGAGGACCAGCACAGCCAGGGTGGCTTTGTGGTACAGCCCGCGGCTC
TGGCAAACCGAAAAGGATGACTGGACGGTCCCTTATGGCCGATCTACTTGGCCGCGAGCACACCGCC
TACCCGACCGGCTGGGTGGAGACGGCGGTAAGTCCGGCGTGCAGCCGCGCATCAAGATCAACAGCCGG
AAGGGGCTGCATCGGACACGGCCAGCCCCGAGGGGACGCATCTGACATGGAGGGGAGGGGATGTG
CATGGGGTGGCCAGCAGCCCTCGCATGACCTGGCAAAGGAAGAAGGCAGCCACCCTCCAGTCCAAGGC
CAGTTATCTCTCCAAAACACGACCCACAGGACCTCGCATTAA
  
```

Restriction Sites: SgfI-MluI



[View online »](#)

ACCN:	NM_001258017
Insert Size:	1770 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001258017.1
RefSeq Size:	2335 bp
RefSeq ORF:	1770 bp
Locus ID:	259307
UniProt ID:	Q96RQ9
Cytogenetics:	19q13.33
Protein Families:	Druggable Genome
Protein Pathways:	Alanine, aspartate and glutamate metabolism, Cysteine and methionine metabolism, Metabolic pathways, Phenylalanine, tyrosine and tryptophan biosynthesis, Phenylalanine metabolism, Tryptophan metabolism, Tyrosine metabolism, Valine, leucine and isoleucine degradation
MW:	65.3 kDa
Gene Summary:	<p>This gene encodes a secreted L-amino acid oxidase protein which primarily catabolizes L-phenylalanine and, to a lesser extent, L-arginine. The expression of this gene is induced by the cytokine interleukin 4 in B cells. This gene is also expressed in macrophages and dendritic cells. This protein may play a role immune system escape as it is expressed in tumor-associated macrophages and suppresses T-cell responses. This protein also contains domains thought to be involved in the binding of flavin adenine dinucleotide (FAD) cofactor. Multiple transcript variants encoding different isoforms have been found for this gene. Some transcripts of this gene share a promoter and exons of the 5' UTR with the overlapping NUP62 gene. [provided by RefSeq, Jul 2020]</p> <p>Transcript Variant: This variant (3) differs in the 5' UTR compared to variant 2. Variants 2, 3 and 4 encode the same protein.</p>