

## Product datasheet for **SC332631**

### TMEM200A (NM\_001258278) Human Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** TMEM200A (NM\_001258278) Human Untagged Clone  
**Tag:** Tag Free  
**Symbol:** TMEM200A  
**Synonyms:** KIAA1913; TTMA; TTMC  
**Vector:** pCMV6-Entry (PS100001)  
**Fully Sequenced ORF:** >SC332631 representing NM\_001258278.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGATAGCAACTGGTGGAGTGATAACTGGCCTGGCCGCTTGAAAAGGCAAGACTCTGCCAGATCACAG
CAGCATGTCAACCTCAGCCCGTCTCCTGCTACCCAAGAGAAGAAGCCCATCAGGCGCCGGCCCCGGGCA
GATGTTGTGGTTGTCGTGGCAAATCCGGCTTTATTCCCATCTGGTTTTTTCTATTTTAGGAGTG
CTCATCTCCATTATAGGAATTGCTATGGCCGTTCTTGGATATTGGCCCCAAAAGAACATTTTATTGAT
GCTGAAACAACACTGTCAACAAATGAAACTCAGGTCATTGCGAATGAAGGCGGTGTGGTGGTTCGCTTC
TTTGAGCAGCATTGCAATCTGATAAGATGAAAATGCTTGGCCCATTCACCATGGGGATTGGCATTTC
ATTTTCATTTGTGCTAATGCCATTCTTCATGAAAACCGTGACAAAAGAGACCAAAATCATACATAGAGG
GATATCTATTCCACAGTCATTGACATTCACACGCTAAGAATCAAGGAGCAAAAGGCAAATGAACGGCATG
TACACTGGTTTGTATGGGAGAAACAGAAGTAAAACAGAATGGGAGCTCCTGTGCCTCGAGATTGGCAGCA
AATACGATCGCCTCTTTCTCGGGTTTTCGGAGCAGTTTTCGAATGGACAGCTCCGTGGAGGAGGATGAA
CTTATGTTAAATGAAGGTAAGAGTTCTGGGCATCTTATGCCCCCTTTGCTCTCTGACAGCTCTGTGTCT
GTCTTTGGCCTCTATCCACCTCCTTCCAAGACAACACTGATGATAAGACCAGCGGCTCTAAGAAATGTGAA
ACCAAGTCAATTGTGTCATCGTCCATCAGTGCTTTTACATTGCCTGTGATCAAACCTAATAACTGTGTT
ATTGATGAGCCCAGTATAGATAACATCACTGAAGATGCTGACAACCTCAAAGTAGTCAAGGAATTTG
TCAATGGATTCCCTTGTGGTTCCTTTGCCAACACCAGTGAATCCTTCCAGCCCGTCAGCACAGTGCTA
CCAAGGAATAATTCCATTGGGGAGTCGTTGTCGAGTCAGTACAAGTCATCTATGGCTCTCGGACCTGGG
GCTGGACAGCTCTGTCTCCTGGGGCTGCCAGAAGACAGTTTGGGTCCAATACATCCTTGCATTTGCTC
TCGTCACACTCAAAGTCTTGGACTTAGACCGGGTCCCTCCACTCTAACTGTTCAAGGAGAACACCGG
AAACATCCAAGTTGGCCTAGGTTGGATCGGAACAACAGCAAGGGATATATGAAACTAGAGAACAAGAA
GACCCGATGGATAGGTTGCTTGTCCCCAAGTTGCCATCAAAAAGGACTTACCAATAAGGAGAAGCTT
CTTATGATTTCAAGATCTACAATAATTTGAGTTTTGAACATGATGAGTTTTGAGTAACAACCTAAAG
AGGGGAACCTCTGAAACAAGTTTAA
  
```

**Restriction Sites:** Sgfl-Mlul  
**ACCN:** NM\_001258278  
**Insert Size:** 1476 bp



[View online »](#)

|                               |   |
|-------------------------------|---|
| <b>OTI Disclaimer:</b>        | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).  |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).  |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol> |
| <b>RefSeq:</b>                | <u><a href="#">NM_001258278.1</a></u>   |
| <b>RefSeq Size:</b>           | 2943 bp   |
| <b>RefSeq ORF:</b>            | 1476 bp   |
| <b>Locus ID:</b>              | 114801  |
| <b>UniProt ID:</b>            | <u><a href="#">Q86VY9</a></u>   |
| <b>Cytogenetics:</b>          | 6q23.1  |
| <b>Protein Families:</b>      | Transmembrane   |
| <b>MW:</b>                    | 54.4 kDa  |