

## Product datasheet for **SC332633**

### LINGO2 (NM\_001258282) Human Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** LINGO2 (NM\_001258282) Human Untagged Clone  
**Tag:** Tag Free  
**Symbol:** LINGO2  
**Synonyms:** LERN3; LRRN6C  
**Vector:** pCMV6-Entry (PS100001)  
**Fully Sequenced ORF:** >SC332633 representing NM\_001258282.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGCTTCACACGGCCATATCATGCTGGCAGCCATTCTGGGTCTGGCTGTGGTGTAAATCTTCATGGGA
TCCACCATTGGCTGCCCGCTCGCTGTGAGTGCTCTGCCAGAACAATCTGTTAGCTGCACAGAAGG
CGATTGATCGCCATCCCAGAGGGCATTCCCATCGAAACAAAATCTTGGACCTCAGTAAAAACAGGCTA
AAAAGCGTCAACCCTGAAGAATTCATATCATATCCTCTGCTGGAAGAGATAGACTTGAGTGACAACATC
ATTGCCAATGTGGAACCAGGAGCATTCAACAATCTCTTTAACCTGCGTTCCTCCGCTAAAAGGCAAT
CGTCTAAAGCTGGTCCCTTTGGGAGTATTCACGGGGCTGTCCAATCTCACTAAGCTTGACATTAGTGAG
AATAAGATTGTCATTTACTAGACTACATGTTCCAAGATCTACATAACCTGAAGTCTCTAGAAGTGGGG
GACAATGATTTGGTTTATATATCACACAGGGCATTTCAGTGGGCTTCTTAGCTTGGAGCAGCTCACCTG
GAGAAATGCAACTTAACAGCAGTACCAACAGAAGCCCTCTCCACCTCCGCAGCCTCATCAGCCTGCAT
CTGAAGCATCTCAATATCAACAATATGCCTGTGTATGCCTTTAAAAGATTGTTCCACCTGAAACACCTA
GAGATTGACTATTGGCCTTACTGGATATGATGCCTGCCAATAGCCTCTACGGTCTCAACCTCACATCC
CTTTCAGTCACCAACACCAATCTGTCTACTGTACCCTTCCTTGCCTTTAAACACCTGGTATACCTGACT
CACCTAACCTCTCCTACAATCCCATCAGCACTATTGAAGCAGGCATGTTCTCTGACCTGATCCGCCTT
CAGGAGCTTCATATAGTGGGGGCCAGCTTCGCACCATTTGAGCCTCACTCCTTCCAAGGGCTCCGCTTC
CTACGCGTGCTCAATGTGTCTCAGAACCTGCTGAAACTTTGGAAGAGAATGTCTTCTCCTCCCCTAGG
GCTCTGGAGGTCTTGAGCATTAAACAACAACCCTCTGGCCTGTGACTGCCGCTTCTCTGGATCTTGAG
CGACAGCCACCCTGCAGTTTGGTGGCCAGCAACCTATGTGTGCTGGCCAGACACCATCCGTGAGAGG
TCTTTCAAGGATTTCCATAGCACTGCCCTTTCTTTTACTTTACCTGCAAAAAACCCAAAATCCGTGAA
AAGAAGTTGCAGCATCTGCTAGTAGATGAAGGGCAGACAGTCCAGCTAGAATGCAGTGCAGATGGAGAC
CCGCAGCCTGTGATTTCTGGGTGACACCCGAAGGCGTTTCATCACCACCAAGTCCAATGGAAGAGCC
ACCGTGTGGGTGATGGCACCTTGAAAATCCGCTTTGCCAGGATCAAGACAGCGGGATGTATGTTTGC
ATCGCTAGCAATGCTGCTGGGAATGATACCTTCACAGCCTCCTTAACTGTGAAAGGATTCGTTTCAGAT
CGTTTTCTTTATGCGAACAGGACCCCTATGTACATGACCGACTCCAATGACACCATTTCCAATGGCACC
AATGCCAATACTTTTCCCTGGACCTTAAAACAATACTGGTGTCTACAGCTATGGGCTGTTTCACATTC
CTGGGAGTGGTTTTATTTGTTTTCTTCTCTTTTTGTGTGGAGCCGAGGGAAAGGCAAGCACAAAAAC
AGCATTGACCTTGAGTATGTGCCAGAAAAACAATGGTGTCTGTTGTGGAAGGGGAGGTAGCTGGACCC
AGGAGGTTCAACATGAAAATGATTGA
  
```



[View online »](#)

<b>Restriction Sites:</b>	Sgfl-RsrII
<b>ACCN:</b>	NM_001258282
<b>Insert Size:</b>	1821 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_001258282.1</a></u>
<b>RefSeq Size:</b>	3156 bp
<b>RefSeq ORF:</b>	1821 bp
<b>Locus ID:</b>	158038
<b>UniProt ID:</b>	<u><a href="#">Q7L985</a></u>
<b>Cytogenetics:</b>	9p21.2-p21.1
<b>Protein Families:</b>	Druggable Genome, Transmembrane
<b>MW:</b>	68.1 kDa