

Product datasheet for **SC332659**

ECT2 (NM_001258316) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ECT2 (NM_001258316) Human Untagged Clone
Tag:	Tag Free
Symbol:	ECT2
Synonyms:	ARHGEF31
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

Fully Sequenced ORF: >SC332659 representing NM_001258316.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGCTGAAAATAGTGTATTAACATCCACTACTGGGAGGACTAGCTTGGCAGACTCTTCCATTTTTGAT
TCTAAAGTTACTGAGATTTCCAAGGAAAACCTACTTATTGGATCTACTTCATATGTAGAAGAAGAGATG
CCTCAGATTTAAACAAGAGTGATATTGGTTCAAGAAGCTGGAAAAACAAGAAGAACCTATAAAAAGCCTTA
AAGGACATTAAGTGGGCTTTGTAAAGATGGAGTCAGTGGAAAGAATTTGAAGTTTGGATCTCCGGAA
TTTGAAAATGATTTGTAGTCACGGACTTTCAGGATTCTGTCTTTAATGACCTCTACAAGGCTGATTGT
AGAGTTATTGGACCACCAGTTGTATTAATTGTTCAAAAAAGGAGAGCCTTTGCCATTTTCATGTCGC
CCGTTGTATTGTACAAGTATGATGAATCTAGTACTATGCTTTACTGGATTTAGGAAAAAAGAAGAACTA
GTCAGGTTGGTGACATTGGTCCATCACATGGGTGGAGTTATTCGAAAAGACTTTAATTCAAAAGTTACA
CATTTGGTGGCAAATGTACACAAGGAGAAAAATTCAGGGTTGCTGTGAGTCTAGTACTCCAATTATG
AAGCCAGAATGGATTTATAAAGCTTGGGAAAGGCGGAATGAACAGGATTTCTATGCAGCAGTTGATGAC
TTTAGAAATGAATTTAAAGTTCCCTCATTTCAAGATTGATTTTAAAGTTTCTGGGATTTTCAGATGAA
GAGAAAACCAATATGGAAGAAATGACTGAAATGCAAGGAGGTAATATTTACCGCTTGGAGATGAAAGA
TGCACTCACCTTTGATTTGAAGAGAATATAGTAAAAGATCTTCCCTTTGAACCTTCAAAGAACTTTAT
GTTGTCAAGCAAGAGTGGTTCTGGGGAAGCATTCAAATGGATGCCCGAGCTGGAGAAACTATGATTTTA
TATGAAAAGGCAAATACTCCTGAGCTCAAGAAATCAGTGTCAATGCTTTCTCAAATACCCCTAACAGC
AATCGCAAACGACGTCGTTTAAAAGAAACACTTGTCTCAGCTTTCAAGAGAGACAGACGTCACCATTT
CCACCCCGTAAGCGCCCATCAGCTGAGCATTCCCTTTCCATAGGGTCACTCCTAGATATCTCAACACA
CCAGAGCTAGCATTAACTATGGAGACACCCCAAAGTCTTGTACTAAGTCTTCTAAAAGTCCACTCCA
GTTCTTCAAAGCAGTCAGCAAGGTGGCAAGTTGCAAAAAGAGCTTTATCAAAGTAAAAGTAAATATGTT
AATATATTGGCAACAATTATTCAGTTATTTCAAGTACCATTGGAAGAGGAAGGACAACGTTGGTGGACCT
ATCCTTGCACCAGAGGATTAAGACTATTTTTGGTAGCATCCCAGATATCTTTGATGTACACACTAAG
ATAAAGGATGATCTTGAAGACCTTATAGTTAATTGGGATGAGAGCAAAAGCATTGGTGACATTTTTCTG
AAATATTCAAAAGATTTGGTAAAAACCTACCCTCCCTTTGTAACCTTCTTTGAAATGAGCAAGGAAACA
ATTATTAATGTGAAAAACAGAAACCAAGATTTTCATGCTTTTCTCAAGATAAACCAAGCAAAACCAGAA
TGTGGACGGCAGAGCCTTGTGAACCTTATCCGACCAGTACAGAGGTTACCCAGTGTTCATTACTT
TTAAATGATCTTAAGAAGCATACAGCTGATGAAAATCCAGACAAAAGCACTTTAGAAAAAGCTATTGGA
TCACTGAAGGAAGTAATGACGCATATTAATGAGGATAAGAGAAAAACAGAAGCTCAAAGCAAATTTTT
GATGTTGTTTATGAAGTAGATGGATGCCAGCTAATCTTTATCTTCTCACCGAAGCTTAGTACAGCGG
GTTGAAACAATTTCTCTAGGTGAGCACCCCTGTGACAGAGGAGAACAAGTAACTCTTCTCCTTCAAT
GATTGCCTAGAGATAGCAAGAAAACGGCACAAGGTTATTGGCACTTTTAGGAGTCCCTATGGCCAAACC
CGACCCCCAGCTTCTCTTAAGCATATTCACCTAATGCCTCTTCTCAGATTAAGAAGGATTGGACATA
AGAGAGACAGAAGATTGCCATAATGCTTTTGCCTTGTGTGAGGCCACCAACAGAGCAGGCAAATGTG
CTACTCAGTTCCAGATGACATCAGATGAACCTTCCAAAAGAAAAGTGGCTAAAGATGCTGTGTCGACAT
GTAGCTAACACCATTTGTAAAGCAGATGCTGAGAATCTTATTTACTGCTGATCCAGAATCCTTTGAA
GTTACAAGAGCATTCTCTTTCTCCAAAAGTCCAAAAGAGCTTTCGAAGGGCTTTATGACATCCCAC
GGCTCAGTGGAGGAAGAAGTCCTTCCAGCAATGATAAGCATGTAATGAGTCGTCTTTCTAGCACATCA
TCATTAGCAGGTATCCCTTCTCCCTCCCTTGTGACGCTTCTTCTTCTTTGAAAGGAGAAGTCATACG
TTAAGTAGATCTACAACCTATTTGATA TGA
  
```

Restriction Sites: SgfI-MluI
ACCN: NM_001258316
Insert Size: 2652 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001258316.1</u>
RefSeq Size:	4422 bp
RefSeq ORF:	2652 bp
Locus ID:	1894
UniProt ID:	<u>Q9H8V3</u>
Cytogenetics:	3q26.31
Protein Families:	Druggable Genome
MW:	100.1 kDa
Gene Summary:	<p>The protein encoded by this gene is a guanine nucleotide exchange factor and transforming protein that is related to Rho-specific exchange factors and yeast cell cycle regulators. The expression of this gene is elevated with the onset of DNA synthesis and remains elevated during G2 and M phases. In situ hybridization analysis showed that expression is at a high level in cells undergoing mitosis in regenerating liver. Thus, this protein is expressed in a cell cycle-dependent manner during liver regeneration, and is thought to have an important role in the regulation of cytokinesis. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Mar 2017]</p> <p>Transcript Variant: This variant (2), as well as variant 3, encodes isoform b.</p>