

## Product datasheet for **SC332681**

### EPS15R (EPS15L1) (NM\_001258376) Human Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** EPS15R (EPS15L1) (NM\_001258376) Human Untagged Clone  
**Tag:** Tag Free  
**Symbol:** EPS15R  
**Synonyms:** EPS15R  
**Vector:** pCMV6-Entry (PS100001)  
**Fully Sequenced ORF:** >SC332681 representing NM\_001258376.  
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGGCGGCGCCGCTCATCCCCTCTCCAGCAGATTCCCACTGGAAATTCGTTGTATGAATCTTATTAC
AAGCAGGTCGATCCGGCATACACAGGGAGGGTGGGGGCGAGTGAAGCTGCGCTTTTTCTAAAGAAGTCT
GGCCTCTCGGACATTATCCTTGGGAAGATATGGGACTTGGCCGATCCAGAAGGTAAGGGTCTTGGAC
AAACAGGGTTTCTATGTTGCACTGAGACTGGTGGCCTGTGCACAGAGTGGCCATGAAGTTACCTTGAGC
AATCTGAATTTGAGCATGCCACCGCCTAAATTTACAGACACCAGCAGCCCTCTGATGGTCACACCCGCC
TCTGCAGAGGCCACTGGGCTGTGAGGGTGAAGAAAAGGCCAAATTTGATGGGATTTTTGAAAGCCTC
TTGCCCATCAATGGTTTGTCTCTGAGAGACAAAGTCAAGCCAGTCTCATGAACCTAAAGCTGCCTCTT
GATGTCCTGGGCAGGGTCTGGGACCTCAGTGACATTGACAAGGATGGGCCTTGGATCGAGATGAGTTC
GCTGTGGCCATGCACTTGGTGTACCGAGCCCTGGAGAAGGAGCCCGTGCCTCCGCCCTGCCCCCGTCC
CTCATCCCACCCTCCAAGAGAAAGAAGACTGTGTTCCCTGGCGCCGTCCCGTCTGCCTGCCAGCCCC
CCACAAAAGACAGCCTCCGCTCCACGCCGTCCACGGCAGCGTCAGCAGCCTCAACAGCACAGGGAGC
CTGTCCCCAAGCACAGCCTCAAGCAAACACAGCAACAGTGAAGTGGTGGTGGCCGTGGCAGACAAG
ATGCGATTTGATGAGATATTCCTGAAGACCGACCTGGACCTGGATGGCTACGTGAGTGGCAGGAGGTG
AAGGAGATCTTCATGCACTCGGGCCTCACCCAGAACCTTCTAGCACACATATGGGCCCTGGCCGATACG
AGGCAAACGGGGAAGTTAAGCAAAGACCAATTTCGCTTAGCTATGTATTTTCATTACAGAGAAGGTCAGT
AAAGGCATCGACCCTCTCAAGTCTCTCGCCGACATGGTCCCGCTTCGGAGAGAGGCACGCCCGGC
CCGGACAGTTCAGGCTCTCTCGGCTCCGGGGAGTTACTGGCGTGAAGGAGCTTGTGACATCAGTCAA
GAGATTGCCAGTTACAAGAGAGAAAATTTCACTGGAACAAGACATTCGAGAAAAGGAAGAGGCAATC
AGACAGAAAACCAGCGAGGTGCAGGAATTACAAAATGACCTAGACCGGAAAACAAGCAGTTTGCAGGAG
CTCGAGGCTCAGAAACAGGATGCTCAAGACCGCTGGACGAGATGGACCAGCAGAAGGCCAAGCTCCGA
GACATGCTGAGCGACGTCCGGCAGAAGTCCAGGATGAGACTCAGATGATCTCATCACTGAAAACGCAA
ATCCAATCTCAGGAATCTGACTTAAAGTCCCAGGAAGACGATCTGAACCGAGCCAAGTCGGAGCTGAAC
CGATTGCAGCAGGAGAAACCCAGCTGGAGCAGAGCATTAGGCTGGGCGAGTCCAGCTGGAACCATC
ATCAAGTCCCTGAAGTCAACGCAAGACGAAATCAACCAGGCAAGGAGCAAACCTTTCCAGCTGCATGAA
AGCCGCCAGGAGGCCACAGGAGCCTGGAGCAGTATGACCAGGTGCTCGATGGAGCCCATGGTGCCAGC
CTGACCGACCTGGCCAACCTGAGCGAAGGGCTCTCCCTGGCAGAGAGGGGCAGTTTTGGAGCCATGGTA
AAGGTTGAATAA
```



[View online »](#)

<b>Restriction Sites:</b>	Sgfl-Mlul
<b>ACCN:</b>	NM_001258376
<b>Insert Size:</b>	1806 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001258376.1</a>
<b>RefSeq Size:</b>	3566 bp
<b>RefSeq ORF:</b>	1806 bp
<b>Locus ID:</b>	58513
<b>UniProt ID:</b>	<a href="#">Q9UBC2</a>
<b>Cytogenetics:</b>	19p13.11
<b>Protein Families:</b>	Druggable Genome
<b>MW:</b>	66.5 kDa
<b>Gene Summary:</b>	<p>Seems to be a constitutive component of clathrin-coated pits that is required for receptor-mediated endocytosis. Involved in endocytosis of integrin beta-1 (ITGB1) and transferrin receptor (TFR); internalization of ITGB1 as DAB2-dependent cargo but not TFR seems to require association with DAB2.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (4) differs in the 3' UTR, compared to variant 1. The encoded isoform (4) is shorter and has a distinct C-terminus, compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>