

Product datasheet for **SC332727**

WDR4 (NM_001260474) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: WDR4 (NM_001260474) Human Untagged Clone
Tag: Tag Free
Symbol: WDR4
Synonyms: GAMOS6; hWH; MIGSB; TRM82; TRMT82; Wuho
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC332727 representing NM_001260474.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGGCGGGCTCTGTGGGACTGGCGTTGTGCGGGCAGACGTTGGTGGTGC GGCGGGCGGCAGCCGATTCTCTG
GCCACCTCCATAGCAAGCAGTGATGATGACAGCCTCTTCATCTATGACTGCAGTGCTGCAGAAAAGAAG
TCACAAGAAAATAAAGGGGAGGACGCGCCCTTGGACCAGGGGAGCGGTGCGATTCTGGCGTCCACCTTC
TCCAAGTCTGGCAGCTATTTTGCTTAACCGATGACAGTAAGCGTCTGATTCTTTTCCGTACAAAACCA
TGGCAATGTCTGAGTGTGAGGACCGTGGCAAGGAGGTGTACAGCCCTGACTTTCATAGCCTCGGAGGAG
AAGGTCTTGGTGGCCGACAAGTCTGGAGACGTCTACTCCTTTTCGGTGTGGAGCCACACGGGTGTGGC
CGTCTAGAGCTGGGGCACCTGTCTATGCTGTAGATGTGGCTGTGAGTCCCTGATGACCGCTTCATCCTC
ACTGCCGACCGGGACGAGAAGATCCGAGTCAGCTGGGCCGCGGCCCCCATAGCATCGAGTCCCTCTGC
TTGGGGCACACAGAGTTTGTGAGCCGTATCTCCGTGGTGCCAACTCAGCCCGGGCTGCTTCTGCTCTCC
TCTGGGGACGGCACCCCTGAGGCTCTGGGAGTACAGGAGCGGCCCGCAGCTGCACTGCTGTACCTGGCC
AGTCTGCAGGAGCTGGTGGACCCCAAGCCCGCCAGTTTGCCCGCTCCAGGATTGCATTCTGGTGCCAG
GAGAACTGCGTGGCGCTCCTGTGCGACGGCACTCCTGTGGTCTACATCTCCAGCTGGACGCCCGCAGA
CAGCAGTTGGTGTACAGGCAGCAGCTGGCGTCCAGCACCAAGTGTGGGACGTGGCTTTCGAGGAGACC
CAGGGGCTGTGGGTGCTCCAGGACTGCCAGGAAGCCCCCTGGTGTCTACAGCCCTGTGGCGACCCAG
TGGCAGTCTGTTCTGAGAGCACCGTGTAAAGAAAGTCTCTGGTGTCTTCGTGGGAACTGGGCCATG
CTGGAAGGCTCTGCCGGCGCAGACGCCAGCTTCAGCAGTCTCTACAAGGCCACGTTTCGACAACGTGACC
TCTACCTGAAGAAGAAAGAGGAGAGACTGCAGCAGCAGCTAGAGAAGAAGCAGCGGCCGCGGAGTCCC
CCGCTGGGCCGACGGGCATGCCAAGAAGATGAGACCGGGGAGGCGACGCTAAGTTGCTGA
```

Restriction Sites: SgfI-MluI
ACCN: NM_001260474
Insert Size: 1236 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001260474.1</u>
RefSeq Size:	2165 bp
RefSeq ORF:	1236 bp
Locus ID:	10785
UniProt ID:	<u>P57081</u>
Cytogenetics:	21q22.3
MW:	45.4 kDa
Gene Summary:	<p>This gene encodes a member of the WD repeat protein family. WD repeats are minimally conserved regions of approximately 40 amino acids typically bracketed by gly-his and trp-asp (GH-WD), which may facilitate formation of heterotrimeric or multiprotein complexes. Members of this family are involved in a variety of cellular processes, including cell cycle progression, signal transduction, apoptosis, and gene regulation. This gene is excluded as a candidate for a form of nonsyndromic deafness (DFNB10), but is still a candidate for other disorders mapped to 21q22.3 as well as for the development of Down syndrome phenotypes. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, May 2012]</p> <p>Transcript Variant: This variant (3) uses an alternate in-frame acceptor splice site at an internal coding exon compared to variant 1, resulting in a shorter isoform (2) missing a single aa compared to isoform 1.</p>