

Product datasheet for **SC332732**

ZNF155 (NM_001260487) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: ZNF155 (NM_001260487) Human Untagged Clone
Tag: Tag Free
Symbol: ZNF155
Synonyms: pHZ-96
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC332732 representing NM_001260487.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGACCACATTCAAGGAGGCAGTGACCTTCAAGGATGTGGCTGTGGTCTTCACTGAGGAGGAGCTGGGG
CTGCTGGACCCTGCCAGAGGAAGCTGTACCGAGATGTGATGCTGGAGAAGCTTCAGGAACCTGCTCTCA
GTGGGGCATCAACCGTTCCACCAAGATACTTGCCACTTCTAAGGGAAGAAAAGTTTTGGATGATGGGG
ACAGCAACCCAAAGAGAAGGGAATTCAGGAGGCAAGATCCAACTGAGTTGGAGTCTGTTCCAGAAGCA
GGAGCACATGAAGAGTGGTCTGCCAGCAAATCTGGGAACAAATGCAAAAGACTTAACCCAGGTCTCAG
GACTCTATCATAAATAACTCTCAGTTCTTTGAAAATGGTATGTCCCCTCCAGGTTGAAGCAGGACTA
CCCACAATTCATACAGGACAGAAACCTTCCCAGGGTGGGAAGTGTAAACAGTCCCTCAGTGATGTTCCC
ATCTTTGATCTTCTCAGCAGTTATACTCAGAAGAGAAGTCTTATACATGTGATGAGTGTGAAAAAGC
ATCTGTTACATCTCAGCTCTTTCATGTTTCATCAGAGAGTCCACGTGGGAGAGAAAACCTTTTATGTGTGAT
GTGTGTGGCAAGGAATTTAGTCAAAGCTCACATCTGCAAACCTCATCAGAGAGTCCACACTGGAGAGAAA
CCATTCAAATGTGAGCAATGTGGGAAAGGTTTCAGTCGTAGATCAGCACTTAATGTTTCATCGTAAATTA
CACACAGGAGAGAAAACCTTACATTTGTGAGGCATGTGGGAAGGCCTTCATTCATGATCCCAGCTTAAG
GAACATAAGAGAATCCATACTGGGAGAAAACCTTCAAATGTGATATATGTGGTAAGACCTTCTATTTT
AGGTCAAAGACTTAAGAGCCATTCATGGTTCACACAGGAGAAAACCTTTAGGTGTGATACATGTGAT
AAGAGCTTTTCATCAGAGATCAGCACTTAATAGGCATTGCATGGTCCACACAGGAGAGAAAACCGTACAGA
TGTGAGCAGTGTGAAAAAGGCTTTATTGGTAGGCTAGATTTTTATAAGCATCAGGTGGTCCACACAGGA
GAGAAACCATATAAATTGTAAGAATGTGGGAAGAGCTTCAGATGGTCTCATGCCTTTTGAACCATCAG
CGAGTCCACAGTGGAGAAAAAGCTTCAAATGTGAAGAATGTGGAAAGGGATTTTATACAAATTCACAA
CTGTCTTCCCATCAGAGATCCCACAGTGGTAAAAAGCCATATAAATGTGAGGAGTGTGGGAAAGGCTAT
GTTACTAAGTTTAACTTCTGACTTGCACCAGAGGGTCCACACGGGAGAGAGACCTTATAATTGTAAGGAA
TGTGGGAAGAATTTAGCCGGGCTCAAGTATTTTGAATCATAAGAGACTCCACTGCCAGAAAAACCA
TTCAAATGTGAGGACTGTGAAAGAGGCTTGTACACAGGACATACCGTAAAGACCAGCCGAGAGACTAT
AGTGGGAAAACCATCCAAATGTGAGGATTGTGGGAGACGCTACAAGAGGCGCTGAATCTGGATATA
CTTTTATCATTATTTCTAAATGACACATAA
```

Restriction Sites: SgfI-MluI
ACCN: NM_001260487



[View online »](#)

Insert Size:	1617 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001260487.1</u>
RefSeq Size:	2615 bp
RefSeq ORF:	1617 bp
Locus ID:	7711
UniProt ID:	<u>Q12901</u>
Cytogenetics:	19q13.31
Protein Families:	Transcription Factors
MW:	62.1 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (4) differs in the 5' UTR compared to variant 1. Variants 1, 2, 3 and 4 encode the same isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.