

## Product datasheet for **SC332733**

### ZNF155 (NM\_001260488) Human Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** ZNF155 (NM\_001260488) Human Untagged Clone  
**Tag:** Tag Free  
**Symbol:** ZNF155  
**Synonyms:** pHZ-96  
**Vector:** pCMV6-Entry (PS100001)  
**Fully Sequenced ORF:** >SC332733 representing NM\_001260488.  
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGATAAAGAGAAGCAAAGTACTTGTATTCTGTTCAAAGATTATGGAGGAGGCAGTGACCTTCAAGGAT
GTGGCTGTGGTCTTCACTGAGGAGGAGCTGGGGCTGCTGGACCCTGCCAGAGGAAGCTGTACCGAGAT
GTGATGCTGGAGAACTTCAGGAACCTGCTCTCAGTGGGGCATCAACCGTTCACCAAGATACTTGCCAC
TTCTAAGGGAAGAAAAGTTTTGGATGATGGGGACAGCAACCCAAAGAGAAGGGAATTCCAGGAGCAAG
ATCCAAACTGAGTTGGAGTCTGTTCCAGAAGCAGGAGCACATGAAGAGTGGTCTGCCAGCAAATCTGG
GAACAAAATGCAAAGACTTAACCAGGTCTCAGGACTCTATCATAAATAACTCTCAGTTCTTTGAAAAT
GGTGATGTCCCCTCCAGGTTGAAGCAGGACTACCCACAATTCATACAGGACAGAAAACCTCCCAGGGT
GGGAAGTGTAACAGTCTTTCAGTGATGTTCCCATCTTTGATCTTCTCAGCAGTTATACTCAGAAGAG
AAGTCTTATACATGTGATGAGTGTGAAAAAGCATCTGTTACATCTCAGCTCTTCATGTTTCATCAGAGA
GTCCACGTGGGAGAGAACTCTTTATGTGTGATGTGTGTGGCAAGGAATTTAGTCAAAGCTCACATCTG
CAAACATCAGAGAGTCCACACTGGAGAGAAACCATTCAAATGTGAGCAATGTGGGAAAGGTTTCAGT
CGTAGATCAGCACTTAATGTTTCATCGTAAATTACACACAGGAGAGAAACCTTACATTTGTGAGGCATGT
GGGAAGGCCTTCATTCATGATCCAGCTTAAGGAACATAAGAGAATCCATACTGGGAGAAACCATTC
AAATGTGATATATGGTGAAGACCTTCTATTTTAGTCAAGCTTAAGAGCCATTCATGGTTTCACACA
GGAGAAAAACCATTTAGGTGTGATACATGTGATAAGAGCTTTCATCAGAGATCAGCACTTAATAGGCAT
TGATGGTCCACACAGGAGAGAAACCGTACAGATGTGAGCAGTGTGAAAAGGCTTTATTGGTAGGCTA
GATTTTTATAAGCATCAGGTGGTCCACACAGGAGAGAAACCATATAATTGTAAGAATGTGGGAAGAGC
TTCAGATGGTCTCATGCCTTTTGAACCATCAGCGAGTCCACAGTGGAGAAAAAGCTTCAAATGTGAA
GAATGTGGAAAGGGATTTTATACAAATTCACAACTGTCTTCCCATCAGAGATCCACAGTGGTGAAGG
CCATATAAATGTGAGGAGTGTGGGAAGGGCTATGTTACTAAGTTAATCTTGACTGCACCAGAGGGTC
CACACGGGAGAGAGACCTTATAATTGAAGGAATGTGGGAAGAACTTTAGCCGGGCTCAAGTATTTTG
AATCATAAGAGACTCCACTGCCAGAAAAACCATCAAATGTGAGGACTGTGAAAAGAGGCTTGTACAC
AGGACATACCGTAAAGACCAGCCGAGAGACTATAGTGGGAAAAACCCATCAAATGTGAGGATTGTGGG
AGACGCTACAAGAGGCGCTTGAATCTGGATATACTTTTATCATTATTTCTAAATGACACATAA
```

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001260488



[View online »](#)

<b>Insert Size:</b>	1650 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_001260488.1</a></u>
<b>RefSeq Size:</b>	2718 bp
<b>RefSeq ORF:</b>	1650 bp
<b>Locus ID:</b>	7711
<b>UniProt ID:</b>	<u><a href="#">Q12901</a></u>
<b>Cytogenetics:</b>	19q13.31
<b>Protein Families:</b>	Transcription Factors
<b>MW:</b>	63.4 kDa
<b>Gene Summary:</b>	<p>May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (5) differs in the 5' UTR and initiates translation at an alternate start codon, compared to variant 1. The encoded isoform (2) is longer and has a distinct N-terminus, compared to variant 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>