

## Product datasheet for **SC332738**

### Radixin (RDX) (NM\_001260493) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Radixin (RDX) (NM_001260493) Human Untagged Clone
Tag:	Tag Free
Symbol:	Radixin
Synonyms:	DFNB24
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC332738 representing NM_001260493. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGCCGAAACCAATCAACGTAAGAGTAACATACTGATGCTGAGCTGGAATTTGCCATTAGCCCAAT
ACAACGGCAAACAACCTTTTGGACCAGGTGGTGAACAGTTGGTTTGCCTGAGGTCTGGTTTTTGGG
CTGCAGTATGTAGACAGCAAAGTTATTCTACATGGCTTAACTAAATAAAAAGGTAACACAGCAGGAT
GTTAAAAAAGAGAATCCTTTACAGTTCAAGTTTAGAGCTAAATCTTTCTGAAGATGTTTCTGAGGAA
TTAATTCAGAAATAACCCAGAGACTCTTCTTCTGCAAGTTAAGAAGCCATCTAAATGATGAGATA
TATTGCCCGCCAGAACTGCAGTCTTTTGGCTTCTATGCTGTCCAAGCCAAGTATGGAGATTACAAT
AAAGAGATTCATAAGCCAGGCTACTGGCTAATGATAGACTCCTACCCAGCGTGTATTGGAACAACAC
AAACTAACAAAAGAACAGTGGGAAGAAAGAATACAGAAGTGGCATGAAGAACATAGAGGAATGTTAAGG
GAGGATTCATGATGGAATACCTGAAGATTGCACAAGTCTAGAAATGTATGGAGTCAACTATTTTGA
ATAAAAAATAAAAAGGAACTGAATTGTGGCTAGGTGTTGATGCTTTGGGTCTGAATATTTATGAGCAT
GACGACAAGTTAACACCTAAAATTTGGTTTTCCCTGGAGTGAATCAGAAATATTTTCAATTAATGACAAA
AAATTTGTTATAAAGCCAATCGACAAAAAGGCACCTGATTTTGTGTTTTATGCACCTCGTCTGAGAATC
AATAAGCGGATTTTGGCCTTATGTATGGGAAACCATGAACTATACATGCGAAGAAGGAAGCCTGATACT
ATTGAAGTACAACAGATGAAGGCTCAGGCTAGGGAGGAGAAACATCAGAAGCAGTTGAAAGGGCACAA
TTAGAGAATGAAAAGAAGAAAAGAGAAAATAGCAGAAAAGGAAAAGAAAAGAAATAGAAGCTGAAAAGGAA
GAGCTAATGGAACGTCTAAAACAATTGAAGAGCAGACAATTAAGCTCAGAAAAGAACTAGAAGAACAG
ACTCGAAAAGCTCTAGAAGTGGATCAAGAACGAAAACGAGCAAAAAGAAAGCAGAACGACTTGAAGG
GAGCGTCGAGCTGCTGAAGAGGCAAGTCTGCCATAGCAAAAACAGCTGCCGACCAGATGAAGAATCAG
GAGCAGCTAGCAGCAGAACTTGTGAATCACTGCCAAGATTGCACCTCTAGAGGAAGCCAAGAAGAAA
AAGGAAGAGGAAGCTACTGAGTGGCAACACAAAGCTTTTGCAGCCAGGAAGACTTGGAAAAGACCAAA
GAAGAGTTAAAACCTGTGATGTCTGCCCCCTCCACCTCCACCACCAGTCATTCTCCAACAGAA
AACGAACATGATGAACACGATGAGAATAATGCTGAAGCTAGTGTGAATTATCAAATGAAGGGGTAAATG
AACCATAGAAGCGAGGAAGAAGCTGTAACCGAAACACAGAAAAATGAGCGTGTTAAGAAGCAACTCAG
GCATTAAGTTCAGAATTAGCCCAAGCCAGAGATGAAACCAAGAAAACACAAAATGATGTTCTTCATGCT
GAGAATGTTAAAGCAGGCCGTGATAAGTACAAGACTCTGCGACAGATTCGACAAGGCAATACAAAGCAG
CGTATCGATGAGTTTGAAGCAATGTGGGGACCCAACTGTATGCATTGTTCCAGATGCGCTCTTGCCAA
AGCAGCATAAAGCAAATGTAA

```



[View online »](#)

<b>Restriction Sites:</b>	Sgfl-Mlul
<b>ACCN:</b>	NM_001260493
<b>Insert Size:</b>	1815 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001260493.1</a>
<b>RefSeq Size:</b>	2511 bp
<b>RefSeq ORF:</b>	1815 bp
<b>Locus ID:</b>	5962
<b>UniProt ID:</b>	<a href="#">P35241</a>
<b>Cytogenetics:</b>	11q22.3
<b>Protein Families:</b>	Druggable Genome
<b>Protein Pathways:</b>	Regulation of actin cytoskeleton
<b>MW:</b>	71 kDa
<b>Gene Summary:</b>	<p>Radixin is a cytoskeletal protein that may be important in linking actin to the plasma membrane. It is highly similar in sequence to both ezrin and moesin. The radixin gene has been localized by fluorescence in situ hybridization to 11q23. A truncated version representing a pseudogene (RDXP2) was assigned to Xp21.3. Another pseudogene that seemed to lack introns (RDXP1) was mapped to 11p by Southern and PCR analyses. Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2012]</p> <p>Transcript Variant: This variant (2) differs in the 3' UTR exon, compared to variant 1. Both variants 1 and 2 encode the same isoform 1.</p>