

Product datasheet for **SC332739**

Radixin (RDX) (NM_001260494) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Radixin (RDX) (NM_001260494) Human Untagged Clone
Tag: Tag Free
Symbol: Radixin
Synonyms: DFNB24
Vector: pCMV6-Entry (PS10001)
Fully Sequenced ORF: >SC332739 representing NM_001260494.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGCTGAGCTGGAATTTGCCATTGACCCCAATACAACCTGGCAAACAACCTTTTACCAGTGTATTGGAA
CAACACAACTAACAAAAGAACAGTGGGAAGAAAGAATACAGAAGTGGCATGAAGAACATAGAGGAATG
TTAAGGGAGGATTCTATGATGGAATACCTGAAGATTGCACAAGATCTAGAAATGTATGGAGTCAACTAT
TTTGAAATAAAAAATAAAAAAGGAAGTGAATTGTGGCTAGGTGTTGATGCTTTGGTCTGAATATTTAT
GAGCATGACGACAAGTTAACACCTAAAAATTGGTTTTCCCTGGAGTGAAATCAGAAATATTTCAATTAAT
GACAAAAATTTGTTATAAAGCCAATCGACAAAAAGGCACCTGATTTTGTGTTTTATGCACCTCGTCTG
AGAATCAATAAGCGGATTTTGGCCTTATGTATGGGAAACCATGAACTATACATGCGAAGAAGGAAGCCT
GATACTATTGAAGTACAACAGATGAAGGCTCAGGCTAGGGAGGAGAAACATCAGAAGCAGTTGAAAGG
GCACAATTAGAGAATGAAAAGAAGAAAAGAGAAAATAGCAGAAAAGGAAAAGGAAAGAAATAGAACGTGAA
AAGGAAGAGCTAATGGAACGTCTAAAACAAATTGAAGAGCAGACAATTAAGCTCAGAAAGAACTAGAA
GAACAGACTCGAAAAGCTCTAGAAGTGGATCAAGAACGAAAACGAGCAAAAAGAAGCAGAACGACTT
GAAAAGGAGCGTCGAGCTGCTGAAGAGGCAAAGTCTGCCATAGCAAAAACAAGCTGCCGACCAGATGAAG
AATCAGGAGCAGCTAGCAGCAGAAGTCTGCTGAATTCAGTCCCAAGATTGCACTTCTAGAGGAAGCCAAG
AAGAAAAAGGAAGAGGAAGCTACTGAGTGGCAACACAAAGCTTTTGCAGCCAGGAAGACTTGGAAAAG
ACCAAAGAAGAGTTAAAAACTGTGATGTCTGCCCCCTCCACCTCCACCACCAGTCATTCCTCCA
ACAGAAAACGAACATGATGAACAGATGAGAATAATGCTGAAGCTAGTGCTGAATTATCAATGAAGGG
GTAATGAACCATAGAAGCGAGGAAGAAGCTGTAACCGAAACACAGAAAAATGAGCGTGTAAAGAAGCAA
CTTCAGGCATTAAGTTCAGAATTAGCCCAAGCCAGAGATGAAACCAAGAAAAACAAAAATGATGTTCTT
CATGCTGAGAATGTTAAGCAGGCCGTGATAAGTACAAGACTCTGCGACAGATTCGACAAGGCAATACA
AAGCAGCGTATCGATGAGTTTGAAGCAATGTGA
  
```

Restriction Sites: SgfI-MluI
ACCN: NM_001260494
Insert Size: 1344 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001260494.1
RefSeq Size:	4123 bp
RefSeq ORF:	1344 bp
Locus ID:	5962
UniProt ID:	P35241
Cytogenetics:	11q22.3
Protein Families:	Druggable Genome
Protein Pathways:	Regulation of actin cytoskeleton
MW:	52.7 kDa
Gene Summary:	<p>Radixin is a cytoskeletal protein that may be important in linking actin to the plasma membrane. It is highly similar in sequence to both ezrin and moesin. The radixin gene has been localized by fluorescence in situ hybridization to 11q23. A truncated version representing a pseudogene (RDXP2) was assigned to Xp21.3. Another pseudogene that seemed to lack introns (RDXP1) was mapped to 11p by Southern and PCR analyses. Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2012]</p> <p>Transcript Variant: This variant (4) lacks two internal coding exons and differs in the 3' region, compared to variant 1. The resulting isoform (3) lacks an internal segment and has a shorter C-terminus, compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>