

Product datasheet for **SC332763**

BCAR3 (NM_001261410) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	BCAR3 (NM_001261410) Human Untagged Clone
Tag:	Tag Free
Symbol:	BCAR3
Synonyms:	AND-34; MIG7; NSP2; SH2D3B
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

Fully Sequenced ORF: >SC332763 representing NM_001261410.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGCCTAAGGAATGCAGTGCTTTCCACGCCCTCTCAGCAGCTCTGTGCTGCTTCTATCACCGCAAGTCT
TTCATTGGAGTCAAGTTCTCCAAGGAGAGGCACATCATGGACAGGACCCCGAGAACTGAAGAAGGAG
CTGGAGGAGGAGCTGCTCCTGAGCAGCGAGGACCTGCGCAGCCATGCCTGGTACCACGGCCGCATCCCC
CGACAGGTGTCTGAAAACCTTGTGCAGCGAGATGGTGACTTCTAGTTCGTGACTCTCTGTCCAGCCCT
GGGAACTTTGTCTGACCTGTCACTGGAAGAACCCTGCTCAGCACTTCAAATCAACCGGACAGTTCTG
CGACTCAGCGAGGCTACAGCCGCTGCACTACAGTTGAGATGGAGAGCTTCGACTCCATCCCCGGC
CTGGTGCCTGCTACGTGGGCAACCGCCGGCCATCTCCAGCAGAGTGGCGCCATCATCTTCCAGCCC
ATCAACAGGACGGTGCCTCTGCGGTGCCTGGAGGAGCATTATGGCACCTCCCCAGGCCAGGCCGGGAG
GGCAGCCTACCAAGGGAAGGCCGGATGTGGCCAAGAGGCTGAGCCTCACCATGGGTGGCGTCCAGGCC
CGAGAGCAGAATTTGCCAGGGGAAACCTCTCAGAAACAAGAAAAGAGTGGTAGCCAGCCCGCTGC
CTGGATCATATGCAGGACAGAAGAGCCTTGTCCCTCAAAGCCCACCAGTCAAGAGACTACCTGCCGATT
GGCTGCAAGCTGCCACCTCAGTCTCGGTGTGGACACAAGCCCTGCCAAACTCACCTGTGTTACAGG
ACGGGAAGCGAGCTGCCCTGAGCCAGAGTGGTTTCGGAGGGTCTCTCAGACGCCAGGGCTGGGGAG
GCGCTGAGGGGATCAGACAGTCAACTGTGCCCTAAGCCCCGCCTAAGCCCTGCAAGGTGCCGTTCTC
AAGGTTCCCTCGTCTCCCTCTGCCTGGCTCAACTCAGAGGCCAATACTGTGAAGTGAACCCAGCGTTT
GCCACAGGCTGCGGCAGGGGAGCAAAGTACCCTCATGTGCCAGGGAAGCCACACAGAAGTGTCCACA
GCCAAGCAGAATGAGGCGCCAGGTCCCCGGAACCTGGCGTCAACTACTTGATCCTTGATGATGATGAC
AGGGAAGACCTTGGGAACCTGCGGCAGCTCAGATGGAGAAGGGGAGTGGGACAAGGGCGAGTTTGTG
ACGCCCCCTTGGAGACTGTCTCTCTCCTCAGGCCAACGAGTTTGAGTCAAAGTCTCTCCCCCTGAG
AATAAGCCCTGGAAACAGCAATGTTGAAACGTGCAAAAGAAGTGTTCACCAACAACGCCCAAGGTC
ATCGCCACGACGTAAGTACTGAGCATGGACTGCAGGGTTGCTAGGATACTGGAGTCTCTGAAGAGATGAGG
AGGAACATGGGGGTGAGCTCAGGCTGGAACCTATTACCTTGCCTCACGGACACCAGTGCCTGAGC
ATAATTGAAAGACACAACAATGGCCATCGGCATTGCACTGGACATTCTGGGATGCACGGGCACTTTG
GAGGACCGAGCGGCCACTCTGAGTAAGATCATCCAGGTGGCGGTGGAAGTGAAGGATTCCATGGGGGAC
CTCTATTCTTCTCAGCTCTCATGAAAGCCCTGGAAATGCCACAGATCACAAGGTTAGAAAAGACGTGG
ACTGCTCTGCGGCACCAGTACACCCAAACTGCCATTCTCTATGAGAAACAGCTGAAGCCCTTCCAGCAA
CTCCTGCATGAAGGCAGAGATCCACATGTGTTCCCCAAACAATGTATCAGTCCCAGTCTGATGCCG
CTTGAGCTTAATGGAGCGCCAGGCTGTGACTTTGAAAGAACCGACATGTGGGAAAAAACGACCAG
AGCTGTGAAATCATGCTGAACATTTGGCAACAGCGGATTCATGGCCGAGGCTGCAGACAGCTACCGG
ATGAATGCTGAGAGGATCCTGGCAGTTTTCAACCAGATGAAGAAATGAATGAAATCTGCAAGACTGAA
TTTCAAATGCGATTGCTATGGGGCAGCAAAGGTGCACAAGTCAATCAGACAGAGAGATATGAGAAATTC
AACCAGATTTAACTGCCCTCTCGCGTAAATTGGAACCTCTCTGTAAGCAGGCGAGAGCTTGA
```

Restriction Sites: SgfI-MluI

ACCN: NM_001261410

Insert Size: 2205 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001261410.1
RefSeq Size:	2893 bp
RefSeq ORF:	2205 bp
Locus ID:	8412
UniProt ID:	O75815
Cytogenetics:	1p22.1
Protein Families:	Druggable Genome
MW:	82.4 kDa
Gene Summary:	<p>Breast tumors are initially dependent on estrogens for growth and progression and can be inhibited by anti-estrogens such as tamoxifen. However, breast cancers progress to become anti-estrogen resistant. Breast cancer anti-estrogen resistance gene 3 was identified in the search for genes involved in the development of estrogen resistance. The gene encodes a component of intracellular signal transduction that causes estrogen-independent proliferation in human breast cancer cells. The protein contains a putative src homology 2 (SH2) domain, a hall mark of cellular tyrosine kinase signaling molecules, and is partly homologous to the cell division cycle protein CDC48. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2012]</p> <p>Transcript Variant: This variant (4) differs in the 5' UTR, lacks a portion of the 5' coding region and initiates translation at an alternate start codon compared to variant 1. The resulting protein (isoform 2) is shorter and has a distinct N-terminus compared to isoform 1.</p>