

Product datasheet for **SC332769**

OS9 (NM_001261420) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: OS9 (NM_001261420) Human Untagged Clone
Tag: Tag Free
Symbol: OS9
Synonyms: ERLEC2; OS-9
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC332769 representing NM_001261420.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGCGGCGAAACGCTGCTGTCCAGTTTGTAGGACTGCTGCTTCTGGGACTCCTGTTACCCGCAAGT
CTGACCGCGGTGTGGGAGCCTGAACCTGGAGGAGCTGAGTGAGATGCGTTATGGGATCGAGATCCTG
CCGTTGCCTGTGATGGGAGGGCAGAGCCAATCTCGGACGTGGTATTGTCTCCTCTAAGTACAAACAG
CGCTATGAGTGTCGCTGCCAGCTGGAGCTATTCACCTCCAGCGTAAAGGGAGGAGAAACACCTGCT
TACCAAGGGCCTGGGATCCCTGAGTTGTTGAGCCCAATGAGAGATGCTCCCTGCTTGTGAAGACAAAG
GACTGGTGGACATATGAATTCTGTTATGGACGCCACATCCAGCAATACCACATGGAAGATTGAGAGATC
AAAGGTGAAGTCTCTATCTCGGCTACTACCAATCAGCCTTCGACTGGGATGATGAAACAGCCAAGGCC
TCCAAGCAGCATCGTCTTAAACGCTACCACAGCCAGACCTATGGCAATGGGTCCAAGTGGCAGCTTAAT
GGGAGGCCCGGGAGGCCGAGGTTGCGTTCCTCTGTGACGAGGGTGCAGGTATCTCTGGGACTACATC
GATCGCGTGGACGAGCCCTTGTCTGCTTATGTGCTGACCATTCGCACTCCTCGGCTCTGCCCCAC
CCTCTCTCCGGCCCCACCCAGTGTGCACCCGAGGCCATCCTCTGTCAACCTTCCCTACAGCCTGAG
GAGTACATGGCCTACGTTAGAGGCAAGCCGTAGACTCAAAGCAGTATGGAGATAAAATCATAGAGGAG
CTGCAAGATCTAGGCCCAAGTGTGGAGTGAGACCAAGTCTGGGTGGCACCCCAAAAGATGCCAGGT
GCGAGCCCGACCAAGGATGACAGTAAGGACTCAGATTTCTGGAAGATGCTTAATGAGCCAGAGACCAG
GCCCCAGGAGGGGAGGAGGTGCCGCTGAGGAGCAGGACCAAGCCTGAGGCAGCAGATTGAGCTTCT
GGTCTCCAATGATTTTCAGAACAACGTGCAGGTCAAAGTCATTGCAAGCCCTGCGGATTTGATTGCA
TTCATAGAGGAGCTGAAAGGTGGAACAAAAAGGGGAAGCCAAATATAGGCCAAGAGCAGCCTGTGGAT
GATGCTGCAGAAGTCCCTCAGAGGGAACAGAGAAGGAAAGGGGTGATCCAGAACGCAGAGAGAGATG
GAAGAAGAGGAGGATGAGGATGAGGATGAGGATGAAGATGAGGATGAACGGCAGTTACTGGGAGAATTT
GAGAAGGAACTGGAAGGGATCCTGCTTCCGTGAGCCGAGACCGGCTCCGTTGAGGATGAAGGCTGGC
ATGGAGCGGAACTGGAACATCATCCAGGAGACAGAAAAGAGCTGGACCCAGATGGGCTGAAGAAG
GAGTCAGAGCGGGATCGGCAATGCTGGCTCTCACATCCACTCTCAACAACTCATAAAAGACTGGAG
GAAAAACAGAGTCCAGAGCTGGTGAAGAAGCACAAAGAAAAGAGGGTGTCCCAAAAAGCTCCCCCA
TCACCCCAACCTACAGGGAAAATTGAGATCAAAATTGTCCGCCATGGGCTGAAGGGACTGAAGAGGGT
GCACGTTGGCTGACTGATGAGGACACGAGAAAACCTCAAGGAGATCTTCTCAATATCTTGGTGCCGGGA
GCTGAAGAGGCCAGAAAGGACGCCAGCGGCAGAAAAGAGCTGGAGAGCAATTACCGCCGGGTGTGGGGC
TCTCCAGTGGGGAGGGCACAGGGGACCTGGACGAATTTGACTTCTGA
  
```



[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001261420
Insert Size:	1842 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001261420.1
RefSeq Size:	2733 bp
RefSeq ORF:	1842 bp
Locus ID:	10956
UniProt ID:	Q13438
Cytogenetics:	12q13.3-q14.1
Protein Families:	Transmembrane
MW:	69.3 kDa
Gene Summary:	<p>This gene encodes a protein that is highly expressed in osteosarcomas. This protein binds to the hypoxia-inducible factor 1 (HIF-1), a key regulator of the hypoxic response and angiogenesis, and promotes the degradation of one of its subunits. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (5) has an alternate splice site in the central coding region and lacks an exon in the 3' coding region, compared to variant 1. The resulting isoform (5) has an additional nt in the central region and lacks a segment in the C-terminal region, compared to isoform 1.</p>